



Item Number SO182
Item Sodium Metabisulfite, Granular, NF, EP, BP,
CAS Number 7681-57-4
Molecular Formula $\text{Na}_2\text{S}_2\text{O}_5$
Molecular Weight 190.11

Assay (SO_2)..... 65.0-67.4%
Assay ($\text{Na}_2\text{S}_2\text{O}_5$)..... 95.0-100.5%
pH (5% Solution)..... 3.5-5.0

MAXIMUM LIMITS

Chloride (Cl)..... 0.05%
Iron (Fe)..... 20 ppm
Thiosulfate (S_2O_3)..... To pass test
Heavy Metals (as Pb)..... 20 ppm
Appearance of Solution..... To pass test
Arsenic..... 4 ppm
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843