


VWR[®] Pierceable Silicone Microplate Sealing Mats

01. FOR USE WITH 96-WELL PLATES

02. MINIMIZE SAMPLE EVAPORATION
& CONTAMINATION




VWR® PIERCEABLE SILICONE MICROPLATE SEALING MATS

- Airtight sealing of storage assay, and separation plates
- Easy penetration
- Self-healing even after multiple piercings
- Eliminates cross-contamination between wells
- Resists coring and tearing
- High chemical resistance


Silicone microplate sealing mats are suitable for high-throughput applications and high-temperature storage. All mats are reusable; autoclave and follow with 10% bleach wash and ethanol rinse.

Certifications: Certified Rnase-, Dnase-, DNA- and Pyrogen-free.


For Use With	Cat. No.	Unit
96-Well Deep Well Plates with 7 mm Round Holes	76311-636	Case of 50 Pack of 10
96-Well Deep Well Plates with 8.5 mm Round Holes	76311-634	Case of 50 Pack of 10
96-Well Deep Well Plates with Square Holes	76311-638	Case of 50 Pack of 10
96-Well PCR Plates	76311-640	Case of 50 Pack of 10



76311-636



76311-634



76311-640