



 **avantor**[™] | **VWR.COM**
delivered by **VWR**[™]

AURELIA[®] **GLOVES**
Exam Glove Catalog





We hereby certify that our Aurelia Gloves range complies and conforms with the following:

- The general safety and performance requirements of FDA Medical Device Regulation for Class 1 medical devices (21 CFR 880.6250) and complies with all general controls (section 513(a)(1)(A) of the Federal Food, Drug, and Cosmetic Act (Act); 21 U.S.C. 360c(a)(1)(A)), which include: Establishment Registration Device Listing Premarket Notification 510(k) Labeling Quality System and Medical Device Reporting (MDR)*
- The general safety and performance requirements of Annex I, Medical Device Regulation (EU) 2017/745 for Class I medical devices. Classification: Class I based on Rule 5 transient use, Annex VIII of the Medical Device Regulation (EU) 2017/745 - With the national standard transposing harmonized standard EN455-1, EN455-2, EN455-3 and EN455-4 and is self-certified as a Class I non-sterile medical device*
- The standard specification requirements for ASTM D6319- Freedom From Holes, Physical Dimensions Test, Physical Requirement Test and Packaging and ASTM F1671-07 – Phi-x174 Bacteriophage Penetration Test*
- EEC regulations concerning the conformity of materials and products that are allowed to come into contact with food. In accordance with Regulation EEC 1935/2004, Regulation EC 10/2011 & Regulation (EC) No 2023/2006. The gloves are suitable for contact with dry, fatty, alcoholic and aqueous food for short term contact.*
- The gloves are certified manufactured according to ISO 9001:2015 and ISO 13485:2016 Quality Management Systems*

AURELIA® GLOVES

NITRILE GLOVES – THE RIGHT FIT FOR YOU

A show of hands—which attributes are most important to you in an examination glove? Stretch & Comfort? Improved gripping potential? Greater tactile sensitivity? Whatever you're looking for, Aurelia Gloves offers products to meet your individual preferences. With 15 different styles to choose from, you're sure to find exactly what makes your hands happy.



Amazing®



VWR Item No.	Size
75997-412	Extra Small
75997-414	Small
75997-416	Medium
75997-418	Large
75997-420	Extra Large

- Advanced thin nitrile formulation for enhanced tactility, flexibility and all-day comfort
- High performance lightweight min 3.2mil finger/2.0mil palm thickness
- Eco-friendly packed 300 pieces per box
- Tested & Approved for Chemotherapy drugs resistance (ASTM D6978-05)



Sonic®



VWR Item No.	Size
75999-306	Extra Small
75999-308	Small
75999-290	Medium
75999-292	Large
75999-294	Extra Large

- Proprietary ultra-comfort nitrile formulation for high tactility
- 3.5mil finger/2.2mil palm thickness for enhanced lightweight durability
- Eco-friendly packed 300 pieces per box
- Tested and Approved for Chemotherapy drugs resistance (ASTM D6978-05)

Transform®



VWR Item No.	Size
75997-496	Extra Small
75997-498	Small
75997-500	Medium
75997-502	Large
75997-504	Extra Large

- Soft Premium Medical Nitrile
- 4mil finger/3.2mil palm thickness
- Enhanced puncture and abrasion resistance nitrile formulation
- 200 pieces per box
- Tested and Approved for Chemotherapy drugs resistance (ASTM D6978-05)



Robust®



VWR Item No.	Size
75997-432	Small
75997-434	Medium
75997-436	Large
75997-438	Extra Large

Bold®



VWR Item No.	Size
75997-460	Extra Small
75997-462	Small
75997-464	Medium
75997-466	Large
75997-468	Extra Large

- Heavy duty nitrile formulation for enhanced strength, durability and protection
- Premium min 5.5mil finger/4.5mil palm thickness with full texture for maximum grip
- Tested and Approved for Chemotherapy drugs resistance (ASTM D6978-05)

- Industrial strength medical grade nitrile formulation for enhanced puncture and abrasion resistance
- Min 6mil finger/5mil palm thickness construction
- High-contrast black finish
- Tested and Approved for Chemotherapy drugs resistance (ASTM D6978-05)



Sonic 100®



VWR Item No.	Size
76492-700	Extra Small
76492-702	Small
76492-704	Medium
76492-706	Large
76492-708	Extra Large

- Proprietary ultra-comfort soft nitrile formulation for high tactility
- Min 3.5mil finger/2.2mil palm thickness for enhanced lightweight durability
- Cobalt Blue
- Tested and Approved for Chemotherapy drugs resistance (ASTM D6978-05)

EMPOWERING PEOPLE TO PERFORM AT THEIR BEST ACROSS MANY INDUSTRIES.

For over 15 years, you've been the driving force behind the success of Aurelia Gloves, making us the preferred provider of examination gloves in all of North America. You told us what's important to you—quality hand protection that provides the right comfort, fit, and sensitivity needed to deliver superior patient care. So we've made it our mission to manufacture a premium quality glove line that provides the protection you need and the comfort you deserve, without charging premium cost. We look forward to continuing to serve as a valuable resource to the healthcare community and a trusted partner to everyone who touches our products.



NO MATTER THE INDUSTRY, WE HAVE YOU COVERED.



DENTAL



MEDICAL



EMS



LABORATORY



INDUSTRIAL



JANITORIAL



AUTOMOTIVE



FOOD SERVICE

QUALITY BY DESIGN

COMFORT

Examination gloves need to feel good over an extended period of time. Our gloves do just that, ensuring a proper fit that flexes and stretches with the movements of your hand and reduces muscle fatigue. With sizes that range from extra small to extra extra large, and a fit that conforms to the subtle nuances of fingers and palms, our gloves are built to feel and function like the skin you're in.

SAFETY

Aurelia gloves are manufactured with the attention to quality and detail that has earned the respect of healthcare professionals everywhere. Our gloves conform to the highest standards set by the American Society of Testing and Materials (ASTM) and are subject to our own rigorous quality assurance testing to ensure dependable barrier protection.

SENSITIVITY

Tactile sensitivity allows you to manipulate instruments and materials effectively and to experience realistic feeling during the most intricate procedures. Our patented Honeycomb® Textured surfacing provides the highest tactile sensitivity and ensures exceptional gripping action in both wet and dry conditions. We offer products constructed from the newest ultra-thin, yet durable, textured material, providing a difference that you can feel.



ISO 9001 Certification
ISO 13485:2003 ISO 13485:2012
Directive 89/686/EEC
ASTM Certification

Current Good Manufacturing Practice (GMP)
Standard Malaysian Gloves (SMG)
Health Canada
FDA Certified (CFG)

GLOBAL STORY. LOCAL IMPACT.

Aurelia Gloves is a division of Supermax, one of the largest glove manufacturers in the world. Healthcare professionals in more than 150 countries across the globe depend on our products for their safety, comfort, quality, and value.



AURELIA CARES

Our corporate social responsibility program ensures that a percentage of our profits is channeled back into the community. From donating free examination gloves to aid dental health in Guatemala to supporting the Oral Cancer Foundation, we are proud to be touching lives and improving health around the world.



MARKETS SERVED

At Aurelia Gloves, we are committed to providing gloves for all markets, including Dental, Medical, EMS, Laboratory, Industrial, Janitorial, Automotive & Food Service Industries



STREAMLINED SUPPLY CHAIN

The Supermax Healthcare difference is our capability in offering alternatives in offering alternatives for product delivery that we customize to our customer's supply chain model, providing a streamlined process that ensures the most cost effective distribution.



GREEN INITIATIVE

Whether it's utilizing solar energy to power our leed gold certified green North American headquarters or reducing the amount of gas used in our manufacturing process, Aurelia is taking important steps to minimize our environmental footprint.



AURELIA® GLOVES



VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2021 Avantor, Inc. All rights reserved.