



# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Revision Date 12/17/2015

Version 1.8

## SECTION 1. Identification

### Product identifier

Product number	800849
Product name	Ethanolamine for synthesis
CAS-No.	141-43-5

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Chemical for synthesis
-----------------	------------------------

### Details of the supplier of the safety data sheet

Company	EMD Millipore Corporation   290 Concord Road, Billerica, MA 01821, United States of America   General Inquiries: +1-978-715-4321   Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5)
---------	--

Emergency telephone	800-424-9300 CHEMTREC (USA) +1-703-527-3887 CHEMTREC (International) 24 Hours/day; 7 Days/week
---------------------	--

## SECTION 2. Hazards identification

### GHS Classification

Flammable liquid, Category 4, H227  
Acute toxicity, Category 4, Oral, H302  
Acute toxicity, Category 4, Inhalation, H332  
Acute toxicity, Category 4, Dermal, H312  
Skin corrosion, Category 1B, H314  
Serious eye damage, Category 1, H318

For the full text of the H-Statements mentioned in this Section, see Section 16.

### GHS-Labeling

#### Hazard pictograms



#### Signal Word

Danger

#### Hazard Statements

H227 Combustible liquid.

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 800849  
Product name Ethanolamine for synthesis

Version 1.8

H302 + H312 + H332 Harmful if swallowed, in contact with skin or if inhaled.  
H314 Causes severe skin burns and eye damage.

### *Precautionary Statements*

P210 Keep away from heat/sparks/open flames/hot surfaces. No smoking.  
P261 Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray.  
P264 Wash skin thoroughly after handling.  
P270 Do not eat, drink or smoke when using this product.  
P271 Use only outdoors or in a well-ventilated area.  
P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.  
P301 + P312 IF SWALLOWED: Call a POISON CENTER or doctor/ physician if you feel unwell.  
P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.  
P303 + P361 + P353 IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/ shower.  
P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.  
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P310 Immediately call a POISON CENTER or doctor/ physician.  
P322 Specific measures (see supplemental first aid instructions on this label).  
P363 Wash contaminated clothing before reuse.  
P370 + P378 In case of fire: Use dry sand, dry chemical or alcohol-resistant foam for extinction.  
P403 + P235 Store in a well-ventilated place. Keep cool.  
P405 Store locked up.  
P501 Dispose of contents/ container to an approved waste disposal plant.

### **Other hazards**

None known.

---

## **SECTION 3. Composition/information on ingredients**

Formula	$\text{NH}_2\text{CH}_2\text{CH}_2\text{OH}$	$\text{C}_2\text{H}_7\text{NO}$ (Hill)
Molar mass	61.08 g/mol	

### **Hazardous ingredients**

*Chemical Name (Concentration)*

CAS-No.

*ethanolamine (>= 90 % - <= 100 %)*

141-43-5

Exact percentages are being withheld as a trade secret.

---

## **SECTION 4. First aid measures**

### **Description of first-aid measures**

*General advice*

First aider needs to protect himself.

*Inhalation*

After inhalation: fresh air. Call in physician.

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

800849

Version 1.8

Product name

Ethanolamine for synthesis

---

### *Skin contact*

In case of skin contact: Take off immediately all contaminated clothing. Rinse skin with water/shower. Call a physician immediately.

### *Eye contact*

After eye contact: rinse out with plenty of water. Immediately call in ophthalmologist.

### *Ingestion*

After swallowing: make victim drink water (two glasses at most), avoid vomiting (risk of perforation!). Call a physician immediately. Do not attempt to neutralize.

Never give anything by mouth to an unconscious person.

### **Most important symptoms and effects, both acute and delayed**

Irritation and corrosion, bronchitis, Cough, Shortness of breath, Drowsiness, Nausea  
Risk of blindness!

### **Indication of any immediate medical attention and special treatment needed**

No information available.

---

## **SECTION 5. Fire-fighting measures**

### **Extinguishing media**

#### *Suitable extinguishing media*

Water, Foam, Carbon dioxide (CO<sub>2</sub>), Dry powder

#### *Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### **Special hazards arising from the substance or mixture**

Combustible.

Vapors are heavier than air and may spread along floors.

Forms explosive mixtures with air on intense heating.

Development of hazardous combustion gases or vapors possible in the event of fire.

Fire may cause evolution of:

nitrous gases, nitrogen oxides

### **Advice for firefighters**

#### *Special protective equipment for fire-fighters*

Stay in danger area only with self-contained breathing apparatus. Prevent skin contact by keeping a safe distance or by wearing suitable protective clothing.

#### *Further information*

Remove container from danger zone and cool with water. Suppress (knock down) gases/vapors/mists with a water spray jet. Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

## **SECTION 6. Accidental release measures**

### **Personal precautions, protective equipment and emergency procedures**

Advice for non-emergency personnel: Do not breathe vapors, aerosols. Avoid substance contact.

Ensure adequate ventilation. Keep away from heat and sources of ignition. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders:

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 800849  
Product name Ethanolamine for synthesis

Version 1.8

Protective equipment see section 8.

## Environmental precautions

Do not let product enter drains.

## Methods and materials for containment and cleaning up

Cover drains. Collect, bind, and pump off spills. Observe possible material restrictions (see sections 7 and 10). Take up with liquid-absorbent and neutralizing material (e.g. Chemizorb® OH<sup>-</sup>, Art. No. 101596). Dispose of properly. Clean up affected area.

---

## SECTION 7. Handling and storage

### Precautions for safe handling

Observe label precautions.

Work under hood. Do not inhale substance/mixture. Avoid generation of vapors/aerosols.

### *Advice on protection against fire and explosion*

Keep away from open flames, hot surfaces and sources of ignition. Take precautionary measures against static discharge.

### Conditions for safe storage, including any incompatibilities

Tightly closed.

Store at +15°C to +25°C (+59°F to +77°F).

---

## SECTION 8. Exposure controls/personal protection

### Exposure limit(s)

#### *Ingredients*

Basis	Value	Threshold limits	Remarks
<i>ethanolamine 141-43-5</i>			
ACGIH	Time Weighted Average (TWA):	3 ppm	
	Short Term Exposure Limit (STEL):	6 ppm	
NIOSH/GUIDE	Recommended exposure limit (REL):	3 ppm	
		8 mg/m <sup>3</sup>	
	Short Term Exposure Limit (STEL):	6 ppm	
		15 mg/m <sup>3</sup>	
OSHA_TRANS	PEL:	3 ppm	
		6 mg/m <sup>3</sup>	
Z1A	Time Weighted Average (TWA):	3 ppm	
		8 mg/m <sup>3</sup>	
	Short Term Exposure Limit (STEL):	6 ppm	
		15 mg/m <sup>3</sup>	

### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

800849

Version 1.8

Product name

Ethanolamine for synthesis

---

## Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

### *Hygiene measures*

Immediately change contaminated clothing. Apply skin- protective barrier cream. Wash hands and face after working with substance.

### *Eye/face protection*

Tightly fitting safety goggles

### *Hand protection*

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

### *Other protective equipment:*

Flame retardant antistatic protective clothing.

### *Respiratory protection*

required when vapors/aerosols are generated.

---

## SECTION 9. Physical and chemical properties

Physical state	liquid
Color	colorless
Odor	ammoniacal
Odor Threshold	No information available.
pH	12.1 at 100 g/l 20 °C (20 °C)
Melting point	10.5 °C
Boiling point/boiling range	171 °C (171 °C) at 1,013 hPa
Flash point	92.5 °C (92.5 °C) Method: DIN 51758
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	3.4 %(V)

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

800849

Version 1.8

Product name

Ethanolamine for synthesis

---

Upper explosion limit	27 %(V)
Vapor pressure	0.5 hPa at 20 °C (20 °C)
Relative vapor density	2.1
Density	1.02 g/cm <sup>3</sup> at 20 °C (20 °C)
Relative density	No information available.
Water solubility	at 20 °C (20 °C) soluble
Partition coefficient: n-octanol/water	log Pow: -1.91 (25 °C) OECD Test Guideline 107 Bioaccumulation is not expected.
Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	No information available.
Explosive properties	Not classified as explosive.
Oxidizing properties	none
Ignition temperature	410 °C (410 °C) Method: DIN 51794
Viscosity, kinematic	20 mm <sup>2</sup> /s at 23 °C (23 °C)

---

## SECTION 10. Stability and reactivity

### Reactivity

Forms explosive mixtures with air on intense heating.

A range from approx. 15 Kelvin below the flash point is to be rated as critical.

### Chemical stability

Sensitive to air.

hygroscopic

### Possibility of hazardous reactions

Exothermic reaction with:

Acrolein, Nitriles, chlorosulfonic acid, Hydrogen chloride gas, acetic acid, Acetic anhydride, fuming sulfuric acid, Nitric acid, sulfuric acid, mineral acids, vinyl acetate, Oxidizing agents

Risk of ignition or formation of inflammable gases or vapors with:

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

800849

Version 1.8

Product name

Ethanolamine for synthesis

---

sulfur, iron(III) compounds

Caution! In contact with nitrites, nitrates, nitrous acid possible liberation of nitrosamines!

## Conditions to avoid

Strong heating.

## Incompatible materials

rubber, Copper, Copper alloys

## Hazardous decomposition products

in the event of fire: See section 5.

---

## SECTION 11. Toxicological information

### Information on toxicological effects

#### *Likely route of exposure*

Inhalation, Eye contact, Skin contact

#### *Target Organs*

Eyes

Skin

Respiratory system

Central nervous system

#### *Acute oral toxicity*

LD50 Rat: 1,515 mg/kg (External MSDS)

absorption

Symptoms: If ingested, severe burns of the mouth and throat, as well as a danger of perforation of the esophagus and the stomach.

#### *Acute inhalation toxicity*

Acute toxicity estimate: 11.1 mg/l; vapor

Expert judgment

Corrosive to respiratory system.

Symptoms: mucosal irritations, Shortness of breath, Cough, Possible damages: bronchitis, damage of respiratory tract

absorption

#### *Acute dermal toxicity*

LD50 Rabbit: 1,025 mg/kg

(IUCLID)

absorption

#### *Skin irritation*

Rabbit

Result: Causes burns.

(IUCLID)

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

800849

Version 1.8

Product name

Ethanolamine for synthesis

---

Causes burns.

*Eye irritation*

Rabbit

Result: Causes burns.

(IUCLID)

Causes serious eye damage.

Risk of blindness!

*Genotoxicity in vivo*

In vivo micronucleus test

Mouse

Result: negative

Method: OECD Test Guideline 474

*Genotoxicity in vitro*

Ames test

Salmonella typhimurium

Result: negative

(IUCLID)

Mutagenicity (mammal cell test):

Human lymphocytes

Result: negative

(IUCLID)

*Specific target organ systemic toxicity - single exposure*

The substance or mixture is not classified as specific target organ toxicant, single exposure.

*Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

*Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

## Carcinogenicity

IARC

No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

OSHA

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

NTP

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

ACGIH

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

## Further information



# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 800849  
Product name Ethanolamine for synthesis

Version 1.8

After absorption:

Nausea, Drowsiness

Damage to:

Kidney, Liver

Under given conditions, contact with nitrites or nitric acid can lead to the formation of nitrosamines, which have shown themselves to be carcinogenic in animal experiments.

Handle in accordance with good industrial hygiene and safety practice.

---

## SECTION 12. Ecological information

### Ecotoxicity

#### *Toxicity to fish*

LC50 *Oncorhynchus mykiss* (rainbow trout): 150 mg/l; 96 h neutral (IUCLID)

#### *Toxicity to daphnia and other aquatic invertebrates*

EC5 *E.sulcatum*: 45 mg/l; 72 h neutral (IUCLID)

EC50 *Daphnia magna* (Water flea): 65 mg/l; 48 h (IUCLID)

#### *Toxicity to algae*

IC50 *Desmodesmus subspicatus* (green algae): 22 mg/l; 72 h (IUCLID)

IC5 *Scenedesmus quadricauda* (Green algae): 0.75 mg/l; 8 d neutral (IUCLID)

#### *Toxicity to bacteria*

EC50 activated sludge: > 1,000 mg/l; 3 h

OECD Test Guideline 209

### Persistence and degradability

#### *Biodegradability*

90 - 100 %; 28 d

OECD Test Guideline 301F

Readily biodegradable.

#### *Biochemical Oxygen Demand (BOD)*

800 mg/g (5 d)

(IUCLID)

#### *Theoretical oxygen demand (ThOD)*

1,310 mg/g

(IUCLID)

### Bioaccumulative potential

#### *Partition coefficient: n-octanol/water*

log Pow: -1.91 (25 °C)

OECD Test Guideline 107

Bioaccumulation is not expected.

### Mobility in soil

No information available.

#### *Additional ecological information*

Biological effects:

Harmful effect due to pH shift.

When discharged properly, no impairments in the function of adapted biological wastewater treatment plants are to be expected.

Discharge into the environment must be avoided.

SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 800849  
Product name Ethanolamine for synthesis

Version 1.8

---

**SECTION 13. Disposal considerations**

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

---

**SECTION 14. Transport information**

**Land transport (DOT)**

UN number UN 2491  
Proper shipping name ETHANOLAMINE  
Class 8  
Packing group III  
Environmentally hazardous --

**Air transport (IATA)**

UN number UN 2491  
Proper shipping name ETHANOLAMINE  
Class 8  
Packing group III  
Environmentally hazardous --  
Special precautions for user no

**Sea transport (IMDG)**

UN number UN 2491  
Proper shipping name ETHANOLAMINE  
Class 8  
Packing group III  
Environmentally hazardous --  
Special precautions for user yes  
EmS F-A S-B

---

**SECTION 15. Regulatory information**

**United States of America**

**SARA 313**

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

**SARA 302**

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SAFETY DATA SHEET  
according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 800849  
Product name Ethanolamine for synthesis

---

Version 1.8

**Clean Water Act**

This product does not contain any Hazardous Substances listed under the U.S. CleanWater Act, Section 311, Table 116.4A.

This product does not contain any Hazardous Chemicals listed under the U.S. CleanWater Act, Section 311, Table 117.3.

**DEA List I**

Not listed

**DEA List II**

Not listed

**US State Regulations**

**Massachusetts Right To Know**

*Ingredients*  
ethanolamine

**Pennsylvania Right To Know**

*Ingredients*  
ethanolamine

**New Jersey Right To Know**

*Ingredients*  
ethanolamine

**California Prop 65 Components**

This product does not contain any chemicals known to the State of California to cause cancer, birth, or any other reproductive defects.

**Notification status**

TSCA: All components of the product are listed in the TSCA-inventory.

DSL: All components of this product are on the Canadian DSL

---

**SECTION 16. Other information**

**Training advice**

Provide adequate information, instruction and training for operators.

**Labeling**

*Hazard pictograms*



*Signal Word*  
Danger

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

800849

Version 1.8

Product name

Ethanolamine for synthesis

---

## *Hazard Statements*

H227 Combustible liquid.

H302 + H312 + H332 Harmful if swallowed, in contact with skin or if inhaled.

H314 Causes severe skin burns and eye damage.

## *Precautionary Statements*

### Prevention

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

### Response

P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308 + P310 IF exposed or concerned: immediately call a POISON CENTER or doctor/ physician.

## **Full text of H-Statements referred to under sections 2 and 3.**

H227

Combustible liquid.

H302

Harmful if swallowed.

H312

Harmful in contact with skin.

H314

Causes severe skin burns and eye damage.

H318

Causes serious eye damage.

H332

Harmful if inhaled.

## **Key or legend to abbreviations and acronyms used in the safety data sheet**

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

Revision Date 12/17/2015

---

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*