



SAFETY DATA SHEET

VWR ANALYTICAL

Issuing Date 23-Jul-2012

Revision Date 12-Feb-2015

Revision Number 2

1. PRODUCT AND COMPANY IDENTIFICATION

Product Identifier

Product Name SODIUM FLUORIDE

Other means of identification

Product Code(s) BDH9290-500G

UN-No 1690

Synonyms No information available.

Recommended use of the chemical and restrictions on use

Recommended Use For Laboratory Use Only. Not for Drug, Food, or Household use.

Uses advised against Not for Human or Animal Drug Use

Details of the supplier of the safety data sheet

Company Address

VWR International, LLC
Radnor Corporate Center
100 Matsonford Road
Radnor, PA 19087-8660

Company Phone Number 610-386-1700

E-mail Address www.vwr.com

Emergency Telephone Number

Emergency Telephone Number CHEMTREC 800.424.9300
CANUTEC 613.996.6666

2. HAZARDS IDENTIFICATION

Classification

OSHA Regulatory Status

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Acute Toxicity - Oral	Category 3
Acute Toxicity - Dermal	Category 2
Skin Corrosion/Irritation	Category 2
Serious Eye Damage/Eye Irritation	Category 2

Label elements

Emergency Overview

Danger

Hazard Statements

H301 - Toxic if swallowed

H310 - Fatal in contact with skin

H315 - Causes skin irritation

H319 - Causes serious eye irritation

**Appearance** White to off-white**Physical State** Crystalline Powder**Odor** Odorless**Precautionary Statements**

P262 - Do not get in eyes, on skin, or on clothing

P264 - Wash face, hands and any exposed skin thoroughly after handling

P270 - Do not eat, drink or smoke when using this product

P280 - Wear protective gloves/protective clothing/eye protection/face protection

P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician

P302 + P350 - IF ON SKIN: Gently wash with plenty of soap and water

P302 + P352 - IF ON SKIN: Wash with plenty of soap and water

P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

P310 - Immediately call a POISON CENTER or doctor/physician

P321 - Specific treatment (see supplemental first aid instructions on this label)

P322 - Specific measures (see supplemental first aid instructions on this label)

P330 - Rinse mouth

P332 + P313 - If skin irritation occurs: Get medical advice/attention

P337 + P313 - If eye irritation persists: Get medical advice/attention

P361 - Remove/Take off immediately all contaminated clothing

P362 - Take off contaminated clothing and wash before reuse

P405 - Store locked up

P501 - Dispose of contents/ container to an approved waste disposal plant

Hazards not otherwise classified (HNOC)

Not applicable

Other Information**3. COMPOSITION/INFORMATION ON INGREDIENTS**

Chemical Name	CAS-No	EC No.	Weight %	Trade Secret
Sodium fluoride	7681-49-4	231-667-8	98-100	Not applicable

4. FIRST AID MEASURES**First Aid Measures****Eye Contact**

Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids. Consult a physician.

Skin Contact

Wash off immediately with soap and plenty of water removing all contaminated clothes and shoes. Consult a physician.

Inhalation

Move to fresh air. If not breathing, give artificial respiration.

Ingestion Clean mouth with water and afterwards drink plenty of water. Consult a physician.

Most important symptoms and effects, both acute and delayed

Most Important Symptoms/Effects No information available.

Indication of any immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

5. FIRE-FIGHTING MEASURES

Suitable Extinguishing Media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media No information available.

Specific Hazards Arising from the Chemical

No information available.

Explosion Data

Sensitivity to Mechanical Impact None.

Sensitivity to Static Discharge None.

Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal Precautions Ensure adequate ventilation, especially in confined areas.

Environmental Precautions

Environmental Precautions See Section 12 for additional Ecological information.

Methods and material for containment and cleaning up

Methods for Containment Prevent further leakage or spillage if safe to do so.

Methods for Cleaning Up Avoid dust formation. Pick up and transfer to properly labeled containers. Ventilate area and wash spill site after material pickup is complete.

7. HANDLING AND STORAGE

Precautions for Safe Handling

Handling Handle in accordance with good industrial hygiene and safety practice. Hygroscopic.

Conditions for safe storage, including any incompatibilities

Storage Keep containers tightly closed in a dry, cool and well-ventilated place. Hygroscopic.

Incompatible Products Acids. Oxidizing agents. Metals. Alkalis. Reacts with acids to form hydrogen fluoride.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

Exposure Guidelines ..

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
Sodium fluoride 7681-49-4	TWA: 2.5 mg/m ³	TWA: 2.5 mg/m ³	IDLH: 250 mg/m ³ TWA: 2.5 mg/m ³

Appropriate engineering controls

Engineering Measures Showers
 Eyewash stations
 Ventilation systems.

Individual protection measures, such as personal protective equipment

Eye/Face Protection Avoid contact with eyes. Safety glasses with side-shields.

Skin and Body Protection Wear protective gloves/clothing.

Respiratory Protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.

Hygiene Measures Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

Physical State	Crystalline Powder	Odor	Odorless
Appearance	White to off-white	Odor Threshold	No information available
Color	No information available		

<u>Property</u>	<u>Values</u>	<u>Remarks • Method</u>
pH	No information available.	
Melting point/freezing point	993 °C	
Boiling Point/Range	1700 °C	
Flash Point (High in °C)	No information available.	
Evaporation Rate	No information available	
Flammability (solid, gas)	No information available	
Flammability Limit in Air		
Upper flammability limit:	No information available.	
Lower flammability limit:	No information available.	
Vapor pressure	No information available.	
Vapor Density	No information available.	
Specific Gravity	No information available	
Water Solubility	No information available	
Solubility in other solvents	No information available.	
Partition coefficient	No information available	
Autoignition Temperature		
Decomposition Temperature	No information available	
Kinematic viscosity	No information available	
Dynamic viscosity	No information available.	
Explosive Properties	No information available	
Oxidizing Properties	No information available	

Other Information

Softening Point	No information available.
Molecular Weight	No information available
VOC Content	No information available.
Density	No information available
Bulk Density	No information available

10. STABILITY AND REACTIVITY

Reactivity

No data available

Chemical Stability

Stable under recommended storage conditions.

Possibility of Hazardous Reactions

None under normal processing.

Conditions to Avoid

Protect from heat. Hygroscopic. Keep away from glass and aluminum.

Incompatible Materials

Acids. Oxidizing agents. Metals. Alkalis. Reacts with acids to form hydrogen fluoride.

Hazardous Decomposition Products

Sodium oxides. Fluorides.

11. TOXICOLOGICAL INFORMATION**Information on likely routes of exposure**

Product Information	May be fatal if in contact with skin. Toxic if swallowed Irritating to eye, skin and respiratory system May be harmful if inhaled
Inhalation	There is no data available for this product.
Eye Contact	There is no data available for this product.
Skin Contact	There is no data available for this product.
Ingestion	There is no data available for this product.

Chemical Name	LD50 Oral	LD50 Dermal	LC50 Inhalation
Sodium fluoride 7681-49-4	= 52 mg/kg (Rat)	= 175 mg/kg (Rat)	-

Information on toxicological effects

Symptoms No information available.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Sensitization No information available.
Mutagenic Effects No information available.
Carcinogenicity The table below indicates whether each agency has listed any ingredient as a carcinogen.

Chemical Name	ACGIH	IARC	NTP	OSHA
Sodium fluoride 7681-49-4	-	Group 3	-	-

Reproductive Toxicity No information available.
STOT - single exposure No information available.
STOT - repeated exposure No information available.
Target Organ Effects Kidneys, Bone, Heart, Nerves, Gastrointestinal tract (GI), Teeth, Lungs.
Aspiration hazard No information available.

Numerical measures of toxicity - Product Information

The following values are calculated based on chapter 3.1 of the GHS document .

12. ECOLOGICAL INFORMATION**Ecotoxicity**

2% of the mixture consists of component(s) of unknown hazards to the aquatic environment

Chemical Name	Toxicity to Algae	Toxicity to Fish	Daphnia Magna (Water Flea)
Sodium fluoride 7681-49-4	272: 96 h Pseudokirchneriella subcapitata mg/L EC50 850: 72 h Desmodesmus subspicatus mg/L EC50 static	530: 96 h Lepomis macrochirus mg/L LC50 830: 96 h Lepomis macrochirus mg/L LC50 semi-static 38 - 68: 96 h Oncorhynchus mykiss mg/L LC50 static 180: 96 h Pimephales promelas mg/L LC50 semi-static	338: 48 h Daphnia magna mg/L EC50 98: 48 h Daphnia magna mg/L EC50 Static

Persistence and Degradability

No information available.

Bioaccumulation/Accumulation

No information available.

Chemical Name	Log Pow
Sodium fluoride 7681-49-4	-

Other Adverse Effects

No information available

13. DISPOSAL CONSIDERATIONS

Waste treatment methods**Waste Disposal Method**

Dispose of material in accordance with all federal, state, and local regulations.

Contaminated Packaging

Do not re-use empty containers.

Chemical Name	RCRA - Halogenated Organic Compounds	RCRA - P Series Wastes	RCRA - F Series Wastes	RCRA - K Series Wastes
Sodium fluoride 7681-49-4	-	-	-	-

Chemical Name	California Hazardous Waste Status
Sodium fluoride 7681-49-4	Toxic

14. TRANSPORT INFORMATION

DOT

UN-No 1690
 Proper Shipping Name Sodium fluoride, solid
 Hazard Class 6.1
 Packing Group III

IATA

UN-No 1690
 Proper Shipping Name Sodium fluoride, solid
 Hazard Class 6.1
 Packing Group III

15. REGULATORY INFORMATION

International Inventories

TSCA Complies
 DSL/NDSL Complies
 EINECS/ELINCS Complies

ENCS	Complies
IECSC	Complies
KECL	Complies
PICCS	Complies
AICS	Complies

Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

U.S. Federal Regulations**SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

Component	SARA 313 - Threshold Values %
Sodium fluoride 7681-49-4 (98-100)	-

SARA 311/312 Hazard Categories

Acute Health Hazard	Yes
Chronic Health Hazard	No
Fire Hazard	No
Sudden Release of Pressure Hazard	No
Reactive Hazard	No

Clean Water Act

This product contains the following substances which are regulated pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Chemical Name	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants	CWA - Hazardous Substances
Sodium fluoride 7681-49-4	1000 lb	-	-	X

CERCLA

This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302)

Chemical Name	Hazardous Substances RQs	Extremely Hazardous Substances RQs	RQ
Sodium fluoride 7681-49-4	1000 lb	-	RQ 1000 lb final RQ RQ 454 kg final RQ

U.S. State Regulations**California Proposition 65**

This product does not contain any Proposition 65 chemicals.

Chemical Name	California Prop. 65
Sodium fluoride - 7681-49-4	-

U.S. State Right-to-Know Regulations

Chemical Name	New Jersey	Massachusetts	Pennsylvania
Sodium fluoride 7681-49-4	X	X	X

U.S. EPA Label Information

EPA Pesticide Registration Number Not applicable

Canada

WHMIS Hazard Class

D2B Toxic materials D1B Toxic materials



16. OTHER INFORMATION

Issuing Date 23-Jul-2012

Revision Date 12-Feb-2015

Revision Note

No information available

Disclaimer

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. VWR International and its Affiliates shall not be held liable for any damage resulting from handling.

End of Safety Data Sheet