



# Material Safety Data Sheet

Catalog Number: 102147  
Revision date: 26-Apr-2006

## 1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND COMPANY INFORMATION

Catalog Number: 102147

Product name: LECITHIN SOY REFINED

**Supplier:**

MP Biomedicals, LLC  
29525 Fountain Parkway  
Solon, OH 44139  
tel: 440-337-1200

Emergency telephone number: CHEMTREC: 1-800-424-9300 (1-703-527-3887)

## 2. COMPOSITION/INFORMATION ON INGREDIENTS

Components	CAS Number	Weight %	ACGIH Exposure Limits:	OSHA Exposure Limits:
LECITHIN SOY REFINED	8002-43-5	90 - 100%	None	None

## 3. HAZARDS IDENTIFICATION

**EMERGENCY OVERVIEW:** May cause skin irritation and/or dermatitis

**Principle routes of exposure:** Skin

**Inhalation:** May cause irritation of respiratory tract

**Ingestion:** May be harmful if swallowed.

**Skin contact:** May cause allergic skin reaction

**Eye contact:** Avoid contact with eyes

**Statements of hazard** MAY CAUSE ALLERGIC SKIN REACTION.

**Statement of Spill or Leak - ANSI Label** Eliminate all ignition sources. Absorb and/or contain spill with inert materials (e.g., sand, vermiculite). Then place in appropriate container. For large spills, use water spray to disperse vapors, flush spill area. Prevent runoff from entering waterways or sewers.

## 4. FIRST AID MEASURES

**General advice:** In the case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

**Inhalation:** Move to fresh air. Call a physician immediately.

**Skin contact:** Rinse immediately with plenty of water and seek medical advice

**Ingestion:** Do not induce vomiting without medical advice.

**Eye contact:** In the case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

**Protection of first-aiders:** No information available

**Medical conditions aggravated by exposure:** None known

## 5. FIRE FIGHTING MEASURES

<b>Suitable extinguishing media:</b>	Use dry chemical, CO <sub>2</sub> , water spray or "alcohol" foam
<b>Specific hazards:</b>	Burning produces irritant fumes.
<b>Unusual hazards:</b>	None known
<b>Special protective equipment for firefighters:</b>	As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear
<b>Specific methods:</b>	Water mist may be used to cool closed containers.
<b>Flash point:</b>	Not determined
<b>Autoignition temperature:</b>	Not determined
<b>NFPA rating:</b>	
NFPA Health:	1
NFPA Flammability:	1
NFPA Reactivity:	0

## 6. ACCIDENTAL RELEASE MEASURES

<b>Personal precautions:</b>	Use personal protective equipment.
<b>Environmental precautions:</b>	Prevent product from entering drains.
<b>Methods for cleaning up:</b>	Sweep up and shovel into suitable containers for disposal.

## 7. HANDLING AND STORAGE

**Storage:**

0°C

**Handling:**

Use only in area provided with appropriate exhaust ventilation.

**Safe handling advice:**

Wear personal protective equipment.

**Incompatible products:**

Oxidising and spontaneously flammable products

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

**Engineering measures:** Ensure adequate ventilation.**PERSONAL PROTECTIVE EQUIPMENT****Respiratory protection:** Breathing apparatus only if aerosol or dust is formed.**Hand protection:** Pvc or other plastic material gloves**Skin and body protection:** Usual safety precautions while handling the product will provide adequate protection against this potential effect.**Eye protection:** Safety glasses with side-shields**Hygiene measures:** Handle in accordance with good industrial hygiene and safety practice.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

<b>Physical state:</b>	Powder
<b>Formula:</b>	C <sub>11</sub> H <sub>22</sub> NO <sub>8</sub> P
<b>Molecular weight:</b>	327.27
<b>Melting point/range:</b>	No data available at this time.
<b>Boiling point/range:</b>	No Data available at this time.
<b>Density:</b>	No data available
<b>Vapor pressure:</b>	No data available
<b>Evaporation rate:</b>	No data available

Vapor density:	No data available
Solubility (in water):	No data available
Flash point:	Not determined
Autoignition temperature:	Not determined

## 10. STABILITY AND REACTIVITY

<b>Stability:</b>	Stable under recommended storage conditions.
<b>Polymerization:</b>	None under normal processing.
<b>Hazardous decomposition products:</b>	Nitrogen oxides (NOx)/ammonia/CN- , Potassium oxides (K <sub>2</sub> O)
<b>Materials to avoid:</b>	Strong oxidising agents
<b>Conditions to avoid:</b>	Exposure to air or moisture over prolonged periods.

## 11. TOXICOLOGICAL INFORMATION

### Product Information

#### Acute toxicity

##### Components

LECITHIN SOY REFINED

#### RTECS Number:

Not Available

#### Selected LD50s and LC50s

Not Determined

<b>Chronic toxicity:</b>	Chronic exposure may cause nausea and vomiting, higher exposure causes unconsciousness.
<b>Local effects:</b>	Symptoms of overexposure may be headache, dizziness, tiredness, nausea and vomiting.
<b>Specific effects:</b>	May include moderate to severe erythema (redness) and moderate edema (raised skin), nausea, vomiting, headache.
<b>Primary irritation:</b>	No data is available on the product itself.
<b>Carcinogenic effects:</b>	No data is available on the product itself.
<b>Mutagenic effects:</b>	No data is available on the product itself.
<b>Reproductive toxicity:</b>	No data is available on the product itself.

## 12. ECOLOGICAL INFORMATION

<b>Mobility:</b>	No data available
<b>Bioaccumulation:</b>	No data available
<b>Ecotoxicity effects:</b>	No data available
<b>Aquatic toxicity:</b>	May cause long-term adverse effects in the aquatic environment.

<b>Components</b>	<b>U.S. DOT - Appendix B - Marine Pollutan</b>	<b>U.S. DOT - Appendix B - Severe Marine Pollutants</b>	<b>United Kingdom - The Red List:</b>
LECITHIN SOY REFINED	Not Listed	Not Listed	Not Listed
<b>Components</b>	<b>Germany VCI (WGK)</b>	<b>World Health Organization (WHO) - Drinking Water</b>	<b>Ecotoxicity - Fish Species Data</b>
LECITHIN SOY REFINED	Not Listed	Not Listed	Not Listed
<b>Components</b>	<b>Ecotoxicity - Freshwater Algae Data</b>	<b>Ecotoxicity - Microtox Data</b>	<b>Ecotoxicity - Water Flea Data</b>
LECITHIN SOY REFINED	Not Listed	Not Listed	Not Listed
<b>Components</b>	<b>EPA - ATSDR Priority List</b>	<b>EPA - HPV Challenge Program Chemical List</b>	<b>California - Priority Toxic Pollutants</b>
LECITHIN SOY REFINED	Not Listed	indicator 0; Not sponsored	Not Listed

Components  
LECITHIN SOY REFINED

California - Priority Toxic Pollutants  
Not Listed

California - Priority Toxic Pollutants  
Not Listed

### 13. DISPOSAL CONSIDERATIONS

**Waste from residues / unused products:**

Waste disposal must be in accordance with appropriate Federal, State, and local regulations. This product, if unaltered by use, may be disposed of by treatment at a permitted facility or as advised by your local hazardous waste regulatory authority. Residue from fires extinguished with this material may be hazardous.

**Contaminated packaging:**

Do not re-use empty containers

### 14. TRANSPORT INFORMATION

UN/Id No:

Not regulated

DOT:

Proper shipping name:

Not Regulated

Components  
LECITHIN SOY REFINED

U.S. DOT - Appendix A Table 1 - Reportable Quantities  
Not Listed

TDG (Canada):

WHMIS hazard class:

Non-controlled

IMDG/IMO

IMDG - Hazard Classifications

Not Applicable

Components  
LECITHIN SOY REFINED

U.S. DOT - Appendix B - Marine Pollutan  
Not Listed

U.S. DOT - Appendix B - Severe Marine  
Pollutants  
Not Listed

IMO-labels:

### 15. REGULATORY INFORMATION

**International Inventories**

Components

LECITHIN SOY REFINED

Inventory - United States TSCA - Sect. 8(b)

Present

Canada DSL Inventory List -

Present

Australia (AICS):

Present

Inventory - China:

Present

EU EINECS List -

232-307-2

Inventory - Japan:

8-493

Korean KECL:

KE-21956

Philippines PICCS:

Present

Catalog Number: 102147

Product name: LECITHIN SOY REFINED

Page 4 of 5

**U.S. regulations:****Components**

LECITHIN SOY REFINED

<b>California Proposition 65 -</b> Not Listed	<b>Massachusetts Right to Know List:</b> Not Listed	<b>New Jersey Right to Know List:</b> Not Listed	<b>Pennsylvania Right to Know List:</b> Not Listed
--	--	---	---

**Components**

LECITHIN SOY REFINED

<b>Florida substance List:</b> Not Listed	<b>Rhode Island Right to Know List:</b> Not Listed	<b>Illinois - Toxic Air Contaminants</b> Not Listed	<b>Connecticut - Hazardous Air Pollutants</b> Not Listed
--	---	--	---

**Components**

LECITHIN SOY REFINED

<b>SARA 313 Emission reporting/Toxic Release of Chemicals</b> Not Listed	<b>CERCLA/SARA - Section 302 Extremely Haz</b> Not Listed	<b>NTP:</b> None	<b>IARC:</b> None
---	--	---------------------	----------------------

**SARA 313 Notification:**

The above is your notification as to the SARA 313 listing for this product(s) pursuant to Section 313 of Title III of the Superfund Ammendments and Reauthorization Act of 1986 and 40 CFR Part 372.

If you are unsure if you are subject to the reporting requirements of Section 313, or need more information, please call the EPA Emergency Planning and Community Right-To-Know Information Hotline: (800) 535-0202 or (202) 479-2499 (in Washington, DC or Alaska).

**State Notification:**

The above information is your notice as to the Right-to-Know listings of the stated product(s). Individual states will list chemicals for a variety of reasons including, but not limited to, the compounds toxicity; carcinogenic, tumorigenic and/or reproductive hazards; and the compounds environmental impact if accidentally released.

**16. OTHER INFORMATION**

**Prepared by:** Health & Safety

**Disclaimer:** The information and recommendations contained herein are based upon tests believed to be reliable. However, MP Biomedicals does not guarantee the accuracy or completeness NOR SHALL ANY OF THIS INFORMATION CONSTITUTE A WARRANTY, WHETHER EXPRESSED OR IMPLIED, AS TO THE SAFETY OF THE GOODS, THE MERCHANTABILITY OF THE GOODS, OR THE FITNESS OF THE GOODS FOR A PARTICULAR PURPOSE. Adjustment to conform to actual conditions of usage maybe required. MP Biomedicals assumes no responsibility for results obtained or for incidental or consequential damages, including lost profits arising from the use of these data. No warranty against infringement of any patent, copyright or trademark is made or implied.

**End of Safety Data Sheet**