



NeoForm[®] X8

Neoprene Exam Gloves

PRODUCT DESCRIPTION

SPECIFICATIONS

PRODUCT:	NeoForm X8	
GRADE:	Powder-Free Exam	
MATERIAL:	Neoprene	
THICKNESS	Finger:	6.2mil 0.16mm
	Palm:	4.7mil 0.12mm
	Cuff:	3.8mil 0.10mm
COLOR	Outer:	Green
	Inner:	Green
SURFACE	Outer:	Textured
	Inner:	Smooth
LENGTH:	9.5"	242mm
WEIGHT*:	7.6g	
FIT TYPE:	Ambidextrous	
CUFF TYPE:	Beaded	
WATERTIGHT AQL:	1.5	
ORIGIN COUNTRY:	Indonesia	

PHYSICAL PROPERTIES

	PRE-AGING	POST-AGING
TENSILE STRENGTH:	18 MPa	18 MPa
ELONGATION:	750%	700%
FORCE AT BREAK:		

INTEGRATED TECHNOLOGIES

- Textured surface delivers excellent grip.
- Pure neoprene prevents latex allergic reactions.



*Specification based on size medium gloves

STANDARDS & CERTIFICATIONS

EN 420:2003
 EN 455:2000-2015
 ASTM D6977
 ASTM F1671-97
 ASTM F739-12
 REACH EU 1907/2006
 DIRECTIVE 2011/65/EU
 EU Reg 528/2012
 CAT III COMPLEX DESIGN

FDA 21 CFR 177.2600
 Component materials comply with FDA regulations for food contact.

Tested for use with Chemotherapy Drugs. Gloves used for protection against chemotherapy drug exposure should be selected specifically for the type of chemicals tested. For a complete list of the drugs and concentrations tested for, and testing method used, please call 510-429-8692.



Please visit swsafety.com/symbols for information on standards used on this document.

ORDERING

CAT. NO.	SIZE	PALM WIDTH		PACKAGING	
76051-220	S	7.5-8	3.0" 75mm	100 gloves per box	10 boxes per case
76051-222	M	8.5-9	3.4" 86mm	100 gloves per box	10 boxes per case
76051-224	L	9.5-10	3.8" 97mm	100 gloves per box	10 boxes per case
76051-156	XL	10.5-11	4.2" 107mm	100 gloves per box	10 boxes per case
76051-188	XXL	11.5-12	4.5" 114mm	100 gloves per box	10 boxes per case

APPLICATIONS

- Painting
- Janitorial Maintenance
- Chemical Handling
- Light Chemical Mixing

STORAGE

Keep out of sunlight. Store in a cool, dry place.
 Keep away from sources of ozone or ignition.

