



# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Revision Date 10/14/2015

Version 1.2

## SECTION 1. Identification

### Product identifier

Product number	102123
Product name	Calcium nitrate tetrahydrate 99.95 Suprapur®
CAS-No.	13477-34-4

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Reagent for analysis
-----------------	----------------------

### Details of the supplier of the safety data sheet

Company	EMD Millipore Corporation   290 Concord Road, Billerica, MA 01821, United States of America   General Inquiries: +1-978-715-4321   Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5)
---------	--

Emergency telephone	800-424-9300 CHEMTREC (USA) +1-703-527-3887 CHEMTREC (International) 24 Hours/day; 7 Days/week
---------------------	--

## SECTION 2. Hazards identification

### GHS Classification

Acute toxicity, Category 4, Oral, H302  
Serious eye damage, Category 1, H318

For the full text of the H-Statements mentioned in this Section, see Section 16.

### GHS-Labeling

#### Hazard pictograms



#### Signal Word

Danger

#### Hazard Statements

H302 Harmful if swallowed.  
H318 Causes serious eye damage.

#### Precautionary Statements

P264 Wash skin thoroughly after handling.

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

102123

Version 1.2

Product name

Calcium nitrate tetrahydrate 99.95 Suprapur®

---

P270 Do not eat, drink or smoke when using this product.

P280 Wear protective gloves/ eye protection/ face protection.

P301 + P312 IF SWALLOWED: Call a POISON CENTER or doctor/ physician if you feel unwell.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor/ physician.

P330 Rinse mouth.

P501 Dispose of contents/ container to an approved waste disposal plant.

## Other hazards

None known.

---

## SECTION 3. Composition/information on ingredients

Formula

$\text{Ca}(\text{NO}_3)_2 \cdot 4 \text{H}_2\text{O}$

$\text{CaN}_2\text{O}_6 \cdot 4 \text{H}_2\text{O}$  (Hill)

Molar mass

236.15 g/mol

## Hazardous ingredients

*Chemical Name (Concentration)*

CAS-No.

*calcium nitrate tetrahydrate (>= 90 % - <= 100 % )*

13477-34-4

Exact percentages are being withheld as a trade secret.

---

## SECTION 4. First aid measures

### Description of first-aid measures

*Inhalation*

After inhalation: fresh air.

*Skin contact*

After skin contact: wash off with plenty of water. Remove contaminated clothing.

*Eye contact*

After eye contact: rinse out with plenty of water. Call in ophthalmologist.

*Ingestion*

After swallowing: immediately make victim drink water (two glasses at most). Consult a physician.

Never give anything by mouth to an unconscious person.

### Most important symptoms and effects, both acute and delayed

Nausea, Vomiting, Irritation and corrosion

Risk of serious damage to eyes.

The following applies to nitrites/nitrates in general: methemoglobinemia after the uptake of large quantities.

### Indication of any immediate medical attention and special treatment needed

No information available.

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

102123

Version 1.2

Product name

Calcium nitrate tetrahydrate 99.95 Suprapur®

---

## SECTION 5. Fire-fighting measures

### Extinguishing media

#### *Suitable extinguishing media*

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

#### *Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### Special hazards arising from the substance or mixture

Not combustible.

Has a fire-promoting effect due to release of oxygen.

Ambient fire may liberate hazardous vapors.

Fire may cause evolution of:

nitrous gases, nitrogen oxides

### Advice for firefighters

#### *Special protective equipment for fire-fighters*

In the event of fire, wear self-contained breathing apparatus.

#### *Further information*

Suppress (knock down) gases/vapors/mists with a water spray jet. Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

## SECTION 6. Accidental release measures

### Personal precautions, protective equipment and emergency procedures

Advice for non-emergency personnel: Avoid inhalation of dusts. Avoid substance contact. Ensure adequate ventilation. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

### Environmental precautions

Do not empty into drains.

### Methods and materials for containment and cleaning up

Cover drains. Collect, bind, and pump off spills.

Observe possible material restrictions (see sections 7 and 10).

Take up dry. Dispose of properly. Clean up affected area. Avoid generation of dusts.

---

## SECTION 7. Handling and storage

### Precautions for safe handling

Observe label precautions.

### Conditions for safe storage, including any incompatibilities

Tightly closed. Dry. Do not store near combustible materials.

Store at +5°C to +30°C (+41°F to +86°F).

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

102123

Version 1.2

Product name

Calcium nitrate tetrahydrate 99.95 Suprapur®

---

## SECTION 8. Exposure controls/personal protection

### Exposure limit(s)

Contains no substances with occupational exposure limit values.

### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

### Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

### *Hygiene measures*

Change contaminated clothing. Application of skin- protective barrier cream recommended.

Wash hands after working with substance.

### *Eye/face protection*

Safety glasses

### *Hand protection*

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

### *Other protective equipment:*

protective clothing

### *Respiratory protection*

required when dusts are generated.

---

## SECTION 9. Physical and chemical properties

Physical state	solid
Color	colorless
Odor	odorless
Odor Threshold	Not applicable
pH	5.0 - 7.0 at 50 g/l 25 °C (25 °C)
Melting point	45 °C
Boiling point/boiling range	Not applicable
Flash point	does not flash

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

102123

Version 1.2

Product name

Calcium nitrate tetrahydrate 99.95 Suprapur®

---

Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	Not applicable
Upper explosion limit	Not applicable
Vapor pressure	No information available.
Relative vapor density	Not applicable
Density	1.82 g/cm <sup>3</sup> at 20 °C (20 °C)
Relative density	No information available.
Water solubility	1,470 g/l at 0 °C (0 °C)  2,710 g/l at 40 °C (40 °C)
Partition coefficient: n-octanol/water	No information available.
Autoignition temperature	No information available.
Decomposition temperature	> 130 °C (> 130 °C) Elimination of water of crystallization  561 °C (561 °C)
Viscosity, dynamic	No information available.
Explosive properties	Not classified as explosive.
Oxidizing properties	No information available.
Ignition temperature	Not applicable
Bulk density	ca.1,000 kg/m <sup>3</sup>

---

## SECTION 10. Stability and reactivity

### Reactivity

See below

### Chemical stability

releases water of crystallization when heated.

### Possibility of hazardous reactions

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

102123

Version 1.2

Product name

Calcium nitrate tetrahydrate 99.95 Suprapur®

---

Risk of explosion with:

organic combustible substances, Reducing agents, ammonium compounds

## Conditions to avoid

no information available

## Incompatible materials

no information available

## Hazardous decomposition products

in the event of fire: See section 5.

---

## SECTION 11. Toxicological information

### Information on toxicological effects

#### *Likely route of exposure*

Eye contact, Skin contact, Ingestion

#### *Acute oral toxicity*

LD50 Rat: > 300 - < 2,000 mg/kg

OECD Test Guideline 423

Symptoms: Nausea, Vomiting

absorption

#### *Acute inhalation toxicity*

Symptoms: Possible damages:, mucosal irritations

#### *Eye irritation*

Rabbit

Result: Irreversible effects on the eye

OECD Test Guideline 405

Causes serious eye damage.

#### *Genotoxicity in vitro*

Ames test

Result: negative

(anhydrous substance) (IUCLID)

#### *Specific target organ systemic toxicity - single exposure*

The substance or mixture is not classified as specific target organ toxicant, single exposure.

#### *Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

#### *Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

## Carcinogenicity

IARC

No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

102123

Version 1.2

Product name

Calcium nitrate tetrahydrate 99.95 Suprapur®

---

OSHA	human carcinogen by IARC. No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
NTP	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
ACGIH	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

## Further information

The following applies to nitrites/nitrates in general: methemoglobinemia after the uptake of large quantities.

Handle in accordance with good industrial hygiene and safety practice.

---

## SECTION 12. Ecological information

### Ecotoxicity

No information available.

### Persistence and degradability

No information available.

### Bioaccumulative potential

No information available.

### Mobility in soil

No information available.

### *Additional ecological information*

Depending on the concentration, phosphorus and/or nitrogen compounds may contribute to the eutrophication of drinking- water supplies.

Discharge into the environment must be avoided.

---

## SECTION 13. Disposal considerations

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

SAFETY DATA SHEET  
according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 102123  
Product name Calcium nitrate tetrahydrate 99.95 Suprapur®

Version 1.2

---

**SECTION 14. Transport information**

**Land transport (DOT)**

UN number UN 1454  
Proper shipping name CALCIUM NITRATE  
Class 5.1  
Packing group III  
Environmentally hazardous --

**Air transport (IATA)**

UN number UN 1454  
Proper shipping name CALCIUM NITRATE  
Class 5.1  
Packing group III  
Environmentally hazardous --  
Special precautions for user no

**Sea transport (IMDG)**

UN number UN 1454  
Proper shipping name CALCIUM NITRATE  
Class 5.1  
Packing group III  
Environmentally hazardous --  
Special precautions for user yes  
EmS F-A S-Q

---

**SECTION 15. Regulatory information**

**United States of America**

**SARA 313**

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

**SARA 302**

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

**Clean Water Act**

This product does not contain any Hazardous Substances listed under the U.S. CleanWater Act, Section 311, Table 116.4A.

This product does not contain any Hazardous Chemicals listed under the U.S. CleanWater Act, Section 311, Table 117.3.

This product does not contain any toxic pollutants listed under the U.S. Clean Water Act Section 307

**US State Regulations**



# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

102123

Version 1.2

Product name

Calcium nitrate tetrahydrate 99.95 Suprapur®

---

## Massachusetts Right To Know

### Remarks

No components are subject to the Massachusetts Right to Know Act.

## New Jersey Right To Know

### Ingredients

calcium nitrate tetrahydrate

## California Prop 65 Components

This product does not contain any chemicals known to the State of California to cause cancer, birth, or any other reproductive defects.

## Notification status

TSCA: All components of the product are listed in the TSCA-inventory.

DSL: All components of this product are on the Canadian DSL

---

## SECTION 16. Other information

### Training advice

Provide adequate information, instruction and training for operators.

### Labeling

#### Hazard pictograms



#### Signal Word

Danger

#### Hazard Statements

H302 Harmful if swallowed.

H318 Causes serious eye damage.

#### Precautionary Statements

##### Prevention

P280 Wear eye protection.

##### Response

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P313 Get medical advice/ attention.

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

102123

Version 1.2

Product name

Calcium nitrate tetrahydrate 99.95 Suprapur®

---

## Full text of H-Statements referred to under sections 2 and 3.

H302

Harmful if swallowed.

H318

Causes serious eye damage.

## Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

Revision Date 10/14/2015

---

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*