

**ACE GLASS**

# COATED GLASSWARE


Ace Glass offers many of our existing glassware vessels in an amber coated version. Flasks, Pressure Vessels, Bottles, and many other items listed in our catalog can be amber coated upon request. The coating is a proprietary process and yields excellent UV protection characteristics. If you need help choosing between clear and amber coated glass, contact us today for help.

## Why Amber Coat Your Glassware:

- Amber coating reduces the photo oxidation effects of sunlight and UV radiation.
- Amber coating reduces the transmission of energetic blue light wavelengths.
- Amber coating protects while affording transparency.

*Amber coating is applied to the outside surface only, preventing contamination of the vessel's contents.*

## AMBER



## Special request?

*Ace Glass can coat virtually any piece of glassware.*

**avantor**<sup>™</sup>  
delivered by **VWR**<sup>™</sup>

**VWR.COM**



**ACE GLASS**

# COATED GLASSWARE

## Epoxy Coating:

- Will greatly reduce the risk of breakage due to surface scratches, particularly in pressure application glassware such as pressure reactors, vessels and chromatography columns.
- More chemically resistant than plastic coating.
- Able to withstand a higher temperature (150°C) than plastic coating (120°C).

## Plastic Coating:

- When compared to standard glass vessels, plastic coating improves slip-resistance by lab personnel.
- May help to contain fluids and fragments in the event of accidental breakage.
- Ideal for flasks, bottles, rotary evaporator glassware, etc.

EPOXY  
or  
PLASTIC



## Special request?

*Ace Glass can coat virtually any piece of glassware.*

**avantor**<sup>™</sup>  
delivered by **VWR**<sup>™</sup>

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.