

VWR® CHEMICAL CAPTURE FILTRATION SYSTEM

Continuously Recirculates

Captures any
Potential Emissions

Eliminates Hazardous
Vapors






System Connects Directly to the Cabinet for Added Safety



Continuously draws emitted pollutants eliminating the possibility for combustion



VWR® Chemical Capture Filtration System

- Adaptable to any flammable or acid cabinet
- Connects directly to the rear top bung port
- Effectively removes noxious, odorous and dangerous emissions without the need for any HVAC connection
- Continuously draws 45 CFM of air, capturing all emitted pollutants and relieving any possible stagnation within the cabinet, eliminating the possibility for combustion.

Description	Electrical	Dimensions W x D x H, in.	Cat. No.	Each
Chemical Capture Filtration System	65 W	22½ x 15¾ x 10⅝	52428-821	1129.99

Accessories For Use With Captair 822 Smart and VWR Chemical Capture Filtration System

Description	Cat. No.	Each
AS: Carbon filter for use with volatile organic compounds	10812-238	241.90
BE: Carbon filter for use with Acids, Organics, In-organics and bases	10813-212	357.73
HP: Hepa filter provides 99.995% efficiency for trapping 0.1 µm particles or higher	10813-214	306.48



1.800.932.5000 | vwr.com

Prices and product details are current when published; subject to change without notice. Certain products may be limited by federal, state, provincial, or local regulations. VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.