

Product Specifications

VWR Part Number:	IR2438
ImmunoReagents, Inc. Part Number:	SP-132-VX10
Product:	NORMAL LLAMA SERUM 10ML
Concentration:	60.0 mg/ml (Bradford, IgG standards)
Amount:	10 ml
Conjugate:	Unconjugated
Form:	Lyophilized
Purification:	Lipid extracted and dialyzed against 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Purity:	N/A
Host:	Normal
Immunogen:	N/A
Buffer:	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Storage:	Store lyophilized material at 2-8 °C. For long term storage after reconstitution, prepare small aliquots and store at -20 °C. For storage at 2-8 °C, add a preservative to prevent growth of bacteria.
Shelf Life:	1 year from date of receipt. Prepare working dilution prior to use and then discard.
Specificity:	Non-immune
Cross Reactivity:	N/A
Country of Origin:	Normal Llama Serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.
Applications:	For use as a standard or control in most immunoassay formats.
Disclaimer:	For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.

Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.