

Certificate of Quality

VWR Catalogue Number	89031-278
Description	VWR® Critical Swab® Small Foam-Over-Cotton Tipped Swab, Wooden Handle, Handle Length 152 mm (5.984 inches), Handle width 2.5 mm (0.100 inches), Head Length 19.1 mm (0.750 inches), Head Width 6.4 mm (0.250 inches)
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2022-07-11

Quality System Compliance

Products are manufactured under the [ISO 13485:2016](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Wooden applicator with cotton & foam tipped. Foam tipped applicator is produced with reticulated 100ppi medical grade foam.
Sterilization	The product is non-sterile.
ATP Assay	Product is not tested for ATP Assay.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Product have been shown to be non-toxic.
Non-Pyrogenic Statement	Products were not routinely tested for pyrogens or endotoxins
Latex Statement	The product is Latex free.

BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP), DEHP (di (2- ethylhexyl) phthalate)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store in dry condition and away from oxidizing agents.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.