

Date Prepared: 03/01/2011

Reviewed On: 02/11/2011

1 Identification of the substance/mixture and of the company/undertaking



- **Product Identifier**
- **Product Name:** BBL Dry Slide PYR Kit
- **Catalog Number:** 231747
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
 BD Diagnostic Systems
 7 Loveton Circle
 Sparks, MD 21152
 Telephone.: (410) 771 - 0100 or (800) 638 – 8663
- **Information Department:** Technical Service
- **Emergency telephone number:**
 In case of a chemical emergency, spill, fire, exposure, or accident contact BD Diagnostic Systems (410) 771-0100 or (800)-638-8663, or ChemTrec at (800) 424-9300.

2 Composition/information on ingredients

- **Chemical characterization:** Mixture
- **Description:** Mixture consisting of the following components.

· Dangerous Components:		
109-86-4	2-methoxyethanol	5.0%
64-19-7	acetic acid	2.5%
75-12-7	formamide	2.0%

3 Hazards identification

- **Classification of the substance or mixture**
- **Classification according to Directive 67/548/EEC or Directive 1999/45/EC**
-  **T; Toxic**
 R60-61: May impair fertility. May cause harm to the unborn child.
-  **Xi; Irritant**
 R36/37/38: Irritating to eyes, respiratory system and skin.
- **Classification system:**
 The classification was made according to the latest editions of international substances lists, and expanded upon from company and literature data.
- **Label elements**
- **Labelling according to EU guidelines:**
 The product has been classified and marked in accordance with regulations on hazardous materials.
- **Code letter and hazard designation of product:**
 T Toxic
- **Hazard-determining components of labelling:**
 2-methoxyethanol
 formamide
 acetic acid
- **Risk phrases:**
 60 May impair fertility
 61 May cause harm to the unborn child

(Contd. on page 2)

Date Prepared: 03/01/2011

Reviewed On: 02/11/2011

Product Name: BBL Dry Slide PYR Kit

(Contd. of page 1)

36/37/38 Irritating to eyes, respiratory system and skin.

· Safety phrases:

- 7 Keep container tightly closed.
- 22 Do not breathe dust.
- 24/25 Avoid contact with skin and eyes.
- 36 Wear suitable protective clothing.
- 60 This material and its container must be disposed of as hazardous waste.

· NFPA ratings (scale 0-4)**· HMIS ratings (scale 0-4)**

HEALTH	1	Health = *1
FIRE	0	Fire = 0
REACTIVITY	0	Reactivity = 0

4 First aid measures

- **After inhalation** Seek medical treatment in case of complaints.
- **After skin contact** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact** Rinse opened eye for 15 minutes under running water. Then consult a doctor.
- **After swallowing** Immediately call a doctor.
- **Information for doctor** Show this product label or this MSDS.

5 Firefighting measures

- **Suitable extinguishing agents**
CO₂, ABC multipurpose dry chemical or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Protective equipment:** No special measures required.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Wipe up with damp sponge or mop.
- **Methods and material for containment and cleaning up:** No special measures required.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7 Handling and storage

- **Handling**
- **Precautions for safe handling** No special measures required.
- **Information about protection against explosions and fires:** No special measures required.

(Contd. on page 3)

Date Prepared: 03/01/2011

Reviewed On: 02/11/2011

Product Name: BBL Dry Slide PYR Kit

(Contd. of page 2)

- **Storage**
- **Requirements to be met by storerooms and receptacles:** 15 - 30 °C
- **Information about storage in one common storage facility:**
 Store away from oxidizing agents.
 Do not store together with acids.
 Do not store together with alkaline material.
- **Further information about storage conditions:** Store in cool, dry conditions in well sealed containers.

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see Section 7.

- **Components with limit values that require monitoring at the workplace:**

109-86-4 2-methoxyethanol

PEL	80 mg/m ³ , 25 ppm Skin
REL	0.3 mg/m ³ , 0.1 ppm Skin
TLV	16 mg/m ³ , 5 ppm Skin; BEI

64-19-7 acetic acid

PEL	25 mg/m ³ , 10 ppm
REL	Short-term value: 37 mg/m ³ , 15 ppm Long-term value: 25 mg/m ³ , 10 ppm
TLV	Short-term value: 37 mg/m ³ , 15 ppm Long-term value: 25 mg/m ³ , 10 ppm

75-12-7 formamide

REL	15 mg/m ³ , 10 ppm Skin
TLV	18 mg/m ³ , 10 ppm Skin

- **Additional information:** The lists that were valid during the creation were used as basis.
- **Personal Protective Equipment**
- **General protective and hygienic measures**
 Keep away from foodstuffs, beverages and feed.
 Immediately remove all soiled and contaminated clothing.
 Wash hands before breaks and at the end of work.
 Avoid contact with the eyes and skin.
- **Breathing equipment:**
 In case of brief exposure, use a chemical fume hood or a NIOSH/MSHA-approved respirator.
- **Protection of hands:**



Chemical resistant gloves (i.e. nitrile, or equivalent).

- **Eye protection:** Safety glasses
- **Body protection:** Protective work clothing (lab coat).

USA

(Contd. on page 4)

Date Prepared: 03/01/2011

Reviewed On: 02/11/2011

Product Name: BBL Dry Slide PYR Kit

(Contd. of page 3)

9 Physical and chemical properties

· **General Information**

· **Appearance:**

Form: Liquid
 Color: Yellow
 Odor: Characteristic

· **pH-value:** 3.6 +/-0.2

· **Change in condition:** Undetermined
 · **Melting point/Melting range:** Not determined
 · **Boiling point/Boiling range:** Not determined

· **Flash point:** Not applicable

· **Danger of explosion:** Product does not present an explosion hazard.

· **Density:** Not determined

· **Solubility in / Miscibility with**
 Water: Soluble

10 Stability and reactivity

- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Incompatible materials:** Incompatible materials: oxidizers, acids, bases, copper, brass, and bronze.
- **Hazardous decomposition products:**
 Nitrogen oxides (NO_x)
 Carbon dioxide (CO₂)
 Ammonia (NH₄)

11 Toxicological information

· **Acute toxicity:**

· **LD/LC50 values that are relevant for classification:**

109-86-4 2-methoxyethanol

Oral	LD50	3400 mg/kg (rat)
------	------	------------------

Dermal	LD50	2000 mg/kg (rbt)
--------	------	------------------

75-12-7 formamide

Oral	LD50	5570 mg/kg (rat)
------	------	------------------

· **Primary irritant effect:**

· **on the skin:** Irritant to skin and mucous membranes.

· **on the eye:** Irritating effect.

· **Sensitization:** No sensitizing effects known.

· **Subacute to chronic toxicity:**

Target organs: liver, kidney, bladder, testes, ureter, lungs, thorax,

blood, cardiovascular system, craniofacial abnormalities, gastrointestinal, central nervous system, musculoskeletal system.

· **Additional toxicological information:**

Product is suspected to cause birth defects.

The product shows the following dangers according to internally approved calculation methods for preparations:

(Contd. on page 5)

Date Prepared: 03/01/2011

Reviewed On: 02/11/2011

Product Name: BBL Dry Slide PYR Kit

Irritant

(Contd. of page 4)

12 Ecological information

- **Acquatic toxicity:** No further relevant information available.
- **Ecotoxicological effects:**
- **Other information:**
The ecological effects have not been thoroughly investigated, but currently none have been identified.
- **Additional ecological information:**
- **General notes:** Water hazard class 1 (Self-assessment): slightly hazardous for water.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**
Must not be disposed of with solid waste.
Dispose of material in accordance with federal (40 CFR 261.3), state and local requirements.
This product is not considered a RCRA hazardous waste.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to state and federal regulations.
- **Recommended cleansing agent:** Water, if necessary with cleansing agents.

14 Transport information

- **DOT regulations:**
- **Hazard class:** -
- **Land transport ADR/RID (cross-border)**
- **ADR/RID class:** -
- **Maritime transport IMDG:**
- **IMDG Class:** -
- **Marine pollutant:** No
- **Air transport ICAO-TI and IATA-DGR:**
- **ICAO/IATA Class:** -

- **Transport/Additional information:**
If a dashed line appears in the Hazard Class section for the type of transportation, this indicates the product is not regulated for transportation.

15 Regulatory information

- **SARA Section 355 (extremely hazardous substances)**
None of the ingredients is listed.
- **SARA Section 313 (specific toxic chemical listings)**
109-86-4 | 2-methoxyethanol

(Contd. on page 6)

Date Prepared: 03/01/2011

Reviewed On: 02/11/2011

Product Name: BBL Dry Slide PYR Kit

(Contd. of page 5)

· TSCA (Toxic Substances Control Act)	
109-86-4	2-methoxyethanol
64-19-7	acetic acid
75-12-7	formamide
· California Proposition 65 - Chemicals known to cause cancer	
None of the ingredients is listed.	
· California Proposition 65 - Chemicals known to cause reproductive toxicity for females:	
None of the ingredients is listed.	
· California Proposition 65 - Chemicals known to cause reproductive toxicity for males:	
109-86-4	2-methoxyethanol
· California Proposition 65 - Chemicals known to cause developmental toxicity:	
109-86-4	2-methoxyethanol
· Carcinogenic categories	
· NTP (National Toxicology Program)	
None of the ingredients is listed.	
· TLV (Threshold Limit Value established by ACGIH)	
None of the ingredients is listed.	

16 Other information

To the best of our knowledge, the information contained herein is accurate. However, neither Becton, Dickinson and Company or any of its subsidiaries assumes any liabilities whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we can not guarantee that these are the only hazards that exist.

- **Department issuing MSDS:**
Environmental, Health & Safety
Created by Michael J. Spinazzola
- **Contact:** Technical Service Representative