



Material	58934-540
Material description	Methyl ethyl ketone
Grade	ACS
CAS Number	78-93-3
Molecular formula	C ₄ H ₈ O
Molecular mass	72,11

Characteristics

Specifications

Assay (GC, corrected for water)	Min. 99,0 %
Titration acid	Max. 0,0005 meq/g
Color (Pt-Co)	Max. 15 APHA
Residue after evaporation	Max. 0,0025 %
Water	Max. 0,20 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada