



Material	84364-580
Material description	Sodium peroxide
Grade	ACS
CAS Number	1313-60-6
Molecular formula	Na2O2
Molecular mass	77.98

## Characteristics

## Specifications

Assay	Min. 93.0 %
Heavy metals (as Pb)	Max. 0.002 %
Cl (Chloride)	Max. 0.002 %
PO4 (Phosphate)	Max. 5 ppm
SO4 (Sulphate)	Max. 0.001 %
Fe (Iron)	Max. 0.005 %

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine  
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,  
QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.