



## VWR® Non-Sterile Safety Goggles

VWR Non-Sterile Safety Goggles provide optimum comfort for extended wear and fit over most prescription glasses. Available in two types of ventilation (indirect and direct) and two types of lens coating (uncoated and anti-fog). These economical goggles are designed to protect from hazards with a flexible ridge around a PVC frame and a fully adjustable headband. Goggles provide 99.9% U.V. protection.

Produced under stringent conditions and monitored closely to ensure high levels of quality for every pair. The goggles passed tests according to the American National Standard for Occupational and Educational Personal Eye and Face Protective Devices, ANSI Z87.1-2015.



|                                       | 76017-902                                | 75982-624 | 76017-944              | 75982-626 |
|---------------------------------------|--|-----------|------------------------|-----------|
|                                       | Indirect Ventilation                     |           | Direct Ventilation     |           |
|                                       | Quad vents to minimize lens fogging      |           | Vents maximize airflow |           |
| <b>Features</b>                       | Lightweight PVC frame                    |           |                        |           |
|                                       | Fully adjustable headband                |           |                        |           |
|                                       | Fits over most prescription glasses      |           |                        |           |
|                                       | Meets ANSI Z87.1-2015 standards          |           |                        |           |
| <b>Material:</b>                      |  |           |                        |           |
| <b>Lens Coating</b>                   | Uncoated                                 | Anti-Fog  | Uncoated               | Anti-Fog  |
| <b>Lens</b>                           | Polycarbonate (PC), Anti-UV 99.9%        |           |                        |           |
| <b>Frame</b>                          | PVC                                      |           |                        |           |
| <b>Elastic Band</b>                   | 60% Textured Polyester + 40% Rubber Band |           |                        |           |
| <b>Weight and Dimension:</b>          |  |           |                        |           |
| <b>Total Weight</b>                   | 71.5 g                                   |           | 77.8 g                 |           |
| <b>Lens Thickness</b>                 | 2.0 mm                                   |           | 2.0 mm                 |           |
| <b>Lens Base Curves</b>               | 0  |           | 0                      |           |
| <b>Bridge Width</b>                   | 30 mm                                    |           | 35 mm                  |           |
| <b>Bridge Height</b>                  | 15 mm                                    |           | 20 mm                  |           |
| <b>Overall Width (Hinge to Hinge)</b> | 160 mm                                   |           | 195 mm                 |           |
| <b>Overall Length</b>                 | 70 mm                                    |           | 70 mm                  |           |

## VWR® Non-Sterile Safety Goggles

| Ventilation | Lens Coating           | Cat. No.  | Unit                     |
|-------------|------------------------|-----------|--------------------------|
| Indirect    | Uncoated               | 76017-902 | Pack of 12 / Case of 144 |
| Indirect    | Anti-Fog, Anti-Scratch | 75982-624 | Pack of 12 / Case of 144 |
| Direct      | Uncoated               | 76017-944 | Pack of 12 / Case of 144 |
| Direct      | Anti-Fog, Anti-Scratch | 75982-626 | Pack of 12 / Case of 144 |

