



Item Number	U1018
Item	Urea, FCC
CAS Number	57-13-6
Molecular Formula	CH ₄ N ₂ O
Molecular Weight	60.06
	Assay (CH ₄ N ₂ O)..... 99.0-100.5%
MAXIMUM LIMITS	
	Alcohol-Insoluble Matter..... 0.04%
	Chloride..... 0.007%
	Lead..... 5 mg/kg
	Loss on Drying..... 1.0%
	Residue on Ignition..... 0.1%
	Sulfate..... 0.01%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters:	East Coast Facility:
14422 S. San Pedro St.	755 Jersey Ave.
Gardena, CA 90248	New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843	