

# Clear Media. Clean Packaging.

We know minimizing debris in your controlled environment is a priority, so we've made extensive improvements to our packaging.

15ml filled TSB and FTM needle-port septum vials are now packaged in fully recyclable thermoform plastic trays, and placed in a smooth SPS box with side opening for easy product access.



Description	UOM	VWR Cat. No.
Tryptic Soy Broth (TSB), USP, 15ml	50/bx	89408-060
Fluid Thioglycollate, USP, 15ml	50/bx	89407-582

Hardy Diagnostics' sterility media and rinse solutions are manufactured in an ISO 13485 environmentally controlled production facility. Each lot is certified for pH, sterility and growth promotion using ATCC® strains specified by the USP. Our approach to manufacturing ensures all sterility media and rinse fluids provide the highest level of clarity for improved reading accuracy, thus significantly reducing the risk of incorrect interpretation and false results.

