

# Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Sealing Mats for Square Well Plates
Country of Origin	Manufactured in United Kingdom
Date of Issue (yyyy-mm-dd)	2021-06-17

VWR Catalogue Number	Sterility	Colour
76196-420	Sterile	Clear
76196-422	Non-sterile	Clear

## Quality System Compliance

Products are manufactured under the [ISO 9001:2015](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

## Product Specifications

<b>Material</b>	Ethylene Vinyl Acetate (EVA).
<b>Sterilization</b>	Product labelled as sterile is gamma irradiated with minimum dosage range 5.00 kGy and maximum dosage 25.00 kGy.
<b>ATP Assay</b>	Not applicable.
<b>DNase &amp; RNase Free</b>	This product is free of any detectable DNase/RNase contamination.
<b>BSE/TSE</b>	Not applicable.
<b>Cytotoxicity</b>	Not applicable.
<b>Non-Pyrogenic Statement</b>	Not applicable.
<b>Latex Statement</b>	Not applicable.

<b>BPA Statement</b>	Not applicable.
<b>DEHP Statement</b>	Not applicable.
<b>RoHS</b>	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr <sup>6+</sup> ), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
<b>REACH Statement</b>	No substances as listed on >Substance of Very High Concern< candidate list published 01/19/2021 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
<b>Storage Conditions</b>	Avoid exposure to direct sunlight and store in a dry place.
<b>Shelf Life</b>	3 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.