



MATERIAL SAFETY DATA SHEET

according to the Global Harmonized System

Date of issue: 03/16/2013

Version 1.0

SECTION 1. Identification

Product identifier

Product number 818449
Product name 4,4'-Thiobis(2-tert-butyl-5-methylphenol) for synthesis

Relevant identified uses of the substance or mixture and uses advised against

Identified uses Chemical for synthesis

Details of the supplier of the safety data sheet

Company EMD Millipore Corporation | 290 Concord Road, Billerica, MA 01821,
United States of America | SDS Phone Support: +1-978-715-1335 |
General Inquiries: +1-978-751-4321 | Monday to Friday, 9:00 AM to
4:00 PM Eastern Time (GMT-5)

e-mail: mm_sds@merckgroup.com

Emergency telephone 613-996-6666 CANUTEC (Canada)
+1-703-527-3887 CHEMTREC (International)
24 Hours/day; 7 Days/week

SECTION 2. Hazards identification

GHS Classification

Skin irritation, Category 2, H315
Eye irritation, Category 2, H319

For the full text of the H-Statements mentioned in this Section, see Section 16.

GHS-Labeling

Hazard pictograms



Signal Word
Warning

Hazard Statements

H315 Causes skin irritation.
H319 Causes serious eye irritation.

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 818449
Product name 4,4'-Thiobis(2-tert-butyl-5-methylphenol) for synthesis

Version 1.0

Precautionary Statements

P260 Do not breathe dust.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Other hazards

None known.

SECTION 3. Composition/information on ingredients

Formula	(2-(CH ₃)-4-(OH)-5-((CH ₃) ₃ C)C ₆ H ₂) ₂ S	C ₂₂ H ₃₀ O ₂ S (Hill)
CAS-No.	96-69-5	
Molar mass	358.54 g/mol	

Hazardous ingredients

Chemical Name (Concentration)

CAS-No.

4,4'-thiobis[2-tert.-butyl-5-methylphenol] (>= 90 % - <= 100 %)

96-69-5

SECTION 4. First aid measures

Description of first-aid measures

Inhalation

After inhalation: fresh air.

Skin contact

After skin contact: wash off with plenty of water. Remove contaminated clothing.

Eye contact

After eye contact: rinse out with plenty of water. Call in ophthalmologist.

Ingestion

After swallowing: immediately make victim drink water (two glasses at most). Consult a physician.

Never give anything by mouth to an unconscious person.

Most important symptoms and effects, both acute and delayed

irritant effects, Cough, Shortness of breath

Indication of any immediate medical attention and special treatment needed

No information available.

SECTION 5. Fire-fighting measures

Extinguishing media

Suitable extinguishing media

Water, Carbon dioxide (CO₂), Foam, Dry powder

Unsuitable extinguishing media

For this substance/mixture no limitations of extinguishing agents are given.

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 818449
Product name 4,4'-Thiobis(2-tert-butyl-5-methylphenol) for synthesis

Version 1.0

Special hazards arising from the substance or mixture

Combustible material, Vapors are heavier than air and may spread along floors.
Forms explosive mixtures with air on intense heating.
Development of hazardous combustion gases or vapors possible in the event of fire.
Fire may cause evolution of:
Sulfur oxides

Advice for firefighters

Special protective equipment for fire-fighters

Stay in danger area only with self-contained breathing apparatus. Prevent skin contact by keeping a safe distance or by wearing suitable protective clothing.

Further information

Suppress (knock down) gases/vapors/mists with a water spray jet. Prevent fire extinguishing water from contaminating surface water or the ground water system.

SECTION 6. Accidental release measures

Personal precautions, protective equipment and emergency procedures

Advice for non-emergency personnel: Avoid substance contact. Avoid inhalation of dusts. Ensure adequate ventilation. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

Environmental precautions

Do not empty into drains.

Methods and materials for containment and cleaning up

Cover drains. Collect, bind, and pump off spills.
Observe possible material restrictions (see sections 7 and 10).
Take up dry. Dispose of properly. Clean up affected area. Avoid generation of dusts.

SECTION 7. Handling and storage

Precautions for safe handling

Observe label precautions.

Conditions for safe storage, including any incompatibilities

Tightly closed. Dry.

Store at +15°C to +25°C (+59°F to +77°F).

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 818449
Product name 4,4'-Thiobis(2-tert-butyl-5-methylphenol) for synthesis

Version 1.0

SECTION 8. Exposure controls/personal protection

Exposure limit(s)

Ingredients

Basis	Value	Threshold limits	Remarks
<i>4,4'-thiobis[2-tert.-butyl-5-methylphenol] 96-69-5</i>			
CAD AB OEL	Time Weighted Average (TWA):	10 mg/m ³	
CAD BC OEL	Skin designation:		Can be absorbed through the skin. Form of exposure: Inhalable
	Time Weighted Average (TWA):	1 mg/m ³	Form of exposure: Inhalable
CAD MB OEL	Time Weighted Average (TWA):	1 mg/m ³	Form of exposure: Inhalable fraction.
CAD ON OEL	Time Weighted Average (TWA):	10 mg/m ³	
OEL (QUE)	Time Weighted Average (TWA):	10 mg/m ³	

Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

Hygiene measures

Change contaminated clothing. Application of skin- protective barrier cream recommended.
Wash hands after working with substance.

Eye/face protection

Safety glasses

Hand protection

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

Other protective equipment:

protective clothing

Respiratory protection

required when dusts are generated.

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

SECTION 9. Physical and chemical properties

Physical state powder

Color white

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 818449 Version 1.0
Product name 4,4'-Thiobis(2-tert-butyl-5-methylphenol) for synthesis

Odor	weakly aromatic
Odor Threshold	No information available.
pH	No information available.
Melting point	158 - 161 °C
Boiling point	No information available.
Flash point	419 °F (215 °C)
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	No information available.
Upper explosion limit	No information available.
Vapor pressure	No information available.
Relative vapor density	No information available.
Relative density	No information available.
Water solubility	at 68 °F (20 °C) insoluble
Partition coefficient: n-octanol/water	log Pow: 8.24 (calculated) (Lit.) A remarkable bioaccumulation potential is expected (log Po/w >3).
Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	No information available.
Explosive properties	No information available.
Ignition temperature	756 °F (402 °C) Method: DIN 51794

SECTION 10. Stability and reactivity

Reactivity

Forms explosive mixtures with air on intense heating.

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 818449
Product name 4,4'-Thiobis(2-tert-butyl-5-methylphenol) for synthesis

Version 1.0

The following applies in general to flammable organic substances and mixtures: in correspondingly fine distribution, when whirled up a dust explosion potential may generally be assumed.

Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

Possibility of hazardous reactions

Violent reactions possible with:

Strong oxidizing agents

Conditions to avoid

Strong heating.

A range from approx. 15 Kelvin below the flash point is to be rated as critical.

Incompatible materials

no information available

Hazardous decomposition products

in the event of fire: See section 5.

SECTION 11. Toxicological information

Information on toxicological effects

Likely route of exposure

Inhalation, Eye contact, Skin contact, Ingestion

Acute oral toxicity

LD50 rat: 2,345 mg/kg (RTECS)

Symptoms: Irritations of mucous membranes in the mouth, pharynx, oesophagus and gastrointestinal tract.

Acute inhalation toxicity

Symptoms: mucosal irritations, Cough, Shortness of breath

Skin irritation

Causes skin irritation.

Eye irritation

Causes serious eye irritation.

Specific target organ systemic toxicity - single exposure

The substance or mixture is not classified as specific target organ toxicant, single exposure.

Specific target organ systemic toxicity - repeated exposure

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

Aspiration hazard

Regarding the available data the classification criteria are not fulfilled.

Carcinogenicity

IARC

No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 818449
Product name 4,4'-Thiobis(2-tert-butyl-5-methylphenol) for synthesis

Version 1.0

OSHA No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

NTP No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

ACGIH No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

Further information

Systemic effects:

Swallowing may result in damage to the following:

Liver, Kidney

Further data:

Other dangerous properties can not be excluded.

Handle in accordance with good industrial hygiene and safety practice.

SECTION 12. Ecological information

Ecotoxicity

No information available.

Persistence and degradability

No information available.

Bioaccumulative potential

Partition coefficient: n-octanol/water

log Pow: 8.24

(calculated)

(Lit.) A remarkable bioaccumulation potential is expected (log Po/w >3).

Mobility in soil

No information available.

Other adverse effects

Additional ecological information

We have no quantitative data concerning the ecological effects of this product.

Further information on ecology

Discharge into the environment must be avoided.

SECTION 13. Disposal considerations

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 818449
Product name 4,4'-Thiobis(2-tert-butyl-5-methylphenol) for synthesis

Version 1.0

SECTION 14. Transport information

Land transport (DOT)

Not classified as dangerous in the meaning of transport regulations.

Air transport (IATA)

Not classified as dangerous in the meaning of transport regulations.

Sea transport (IMDG)

Not classified as dangerous in the meaning of transport regulations.

SECTION 15. Regulatory information

Canada

WHMIS Classification

D2B Toxic Material Causing Other Toxic Effects
Skin irritant, Eye irritant

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all of the information required by the CPR.

Notification status

TSCA: On TSCA Inventory

DSL: All components of this product are on the Canadian DSL.

SECTION 16. Other information

Training advice

Provide adequate information, instruction and training for operators.

Full text of H-Statements referred to under sections 2 and 3.

H315 Causes skin irritation.
H319 Causes serious eye irritation.

Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at www.wikipedia.org.

Date of issue: 03/16/2013

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 818449
Product name 4,4'-Thiobis(2-tert-butyl-5-methylphenol) for synthesis

Version 1.0

All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.