

Product Specification

GMP Manufactured Product

Test	Specification
Appearance (Clear, colorless liquid)	Passes Test
Normality	1,98 – 2,02
Ammonium (NH ₄)	≤ 3 ppm
Trace Impurities – Heavy Metals (as Pb)	≤ 0,5 ppm
Trace Impurities – Iron (Fe)	≤ 0,1 ppm
Sulfate (SO ₄)	≤ 1 ppm

GMP Manufactured Product

CAUTION: For Manufacturing, processing or repackaging

Made from U.S.P. Purified Water and Hydrochloric Acid, N.F. (Multi-Compendial) which meets BP, Ph.Eur. and JP Chemical Specifications

No Class 1,2,3 or other solvents are used or produced in the manufacturing or purification of the product.

Filtered through a 0.2 micron filter.

Packaging Site: Paris Mfg Ctr & DC



Jamie Ethier
Vice President Global Quality