

# VWR LAB ORGANIZATION SOLUTIONS

Tube Racks

Pipette Storage

Histology & Microscopy  
Storage

Weighing Boats

Lab Organization



# TUBE RACKS

## VWR® Tube Rack

### Rotate With Ease to Accommodate Your Needs

This multi-purpose rack is designed with the flexibility to hold 20 - 5mL, 20 - 15mL or 17mm tubes or to hold 40 - 5mL, 40 - 15mL or 17mm tubes. The rack incorporates two moving hinged upper shelves that can be rotated to an up or down position to accommodate the size of the tube being used. In addition, the rack can also hold clinical tubes of a similar diameter (17mm) and height.

- Rack holds 5 or 15mL tubes (or tubes up to 17mm in diameter) or a combination of 20 - 5mL, 15mL
- Ergonomic handles for ease of use
- Imprinted grid on both sides of the upper shelves for quick identification and retrieval
- Shelves can be easily rotated to accommodate the different size tubes
- Autoclavable
- Exclusive design

Tubes	Wells	Wells diameter	L x W x H	Quantity	Color	Cat. No.
5/15 mL	40	17 mm	30.9 x 10.7 x 7.7 cm (12.16 x 4.21 x 3.03 ")	1/each	Purple, Natural	75876-392



NEW



Rotate one or both segments down to hold 5mL tubes

Unique tear drop shaped openings for better retention of tubes

## VWR® Customizable Chain Tube Rack


### Multiple Configurations for 5mL/15mL and 17mm Tubes

Work outside the traditional tube rack box shape, or don't. You decide. Now you can configure a tube rack to meet your requirements instead of the other way around. The Chain Rack has a unique linkable design that enables the possibility of multiple configurations to form different shaped tube racks. Each well has the flexibility to hold 5mL tubes on one side or rotate it 180° to hold 15mL or 17mm tubes (round and conical bottom). Mix and match the well sizes to create one rack that holds multiple sized tubes in a variety of fun and functional configurations.

- Chain Rack can be modified to be used in multiple applications including: in water or ice baths, in phlebotomy trays, in a Duraport® specimen transport case or anywhere a tube rack is needed
- Maximize your valuable bench space with the flexibility to hold multiple size tubes (5/15mL and 17mm tubes) in one configurable rack
- Open design allows exposure to water and ice baths and helps reduce identification errors
- Create different shapes and rack sizes with up to 16 wells depending on the need: straight line, a circle, 2 x 8, 4 x 4, etc.
- Link together multiple Chain Rack packs to accommodate larger tube rack configuration requirements
- Assorted colors for easy identification of tubes or testing (blue, green, yellow and natural)
- NOT Autoclavable
- Exclusive design

Tubes	Wells	Wells diameter	Tube Size	Quantity	Color	Cat. No.
5/15 mL	16	17.2 mm (0.7 ")	94 mm x 27 mm (3.7 x 1 ")	16/pk	Blue, Green, Yellow, Natural	75854-098

NEW



Rotate one or both segments up to hold 15mL or clinical tubes




Rotate one or both segments down to hold 5mL tubes

## VWR® Rota-Rack® Duo Tube Rack

### Dilution and Concentration Studies Made Easy

- Made from polypropylene
- Transfer solutions between 0.2mL PCR® tubes and 0.5 and 1.5mL micro centrifuge tubes
- Rotate and customize to suit your application needs. 4 Combinations possible
- Blue section holds 48 x 0.5mL and 48 x 1.5mL microfuge tubes on the reverse side
- Purple section holds 96 x 0.2mL tubes in a standard array pattern one and 48 x 1.5mL tubes on the reverse side
- Tapered wells securely suspend tubes and makes removal easy
- Alphanumeric well identification
- Autoclavable



Rotate and customize to suit your application needs

Tubes	Wells	Wells diameter	L x W x H	Quantity	Color	Cat. No.
0.2 mL	96	5.0 mm	29.6 x 11.8 x 3.4 cm (11.6 x 4.6 x 1.33 ")	1/each	Blue, Purple	89140-272
0.5 mL	48	8.0 mm				
1.5/2.0 mL	96	11 mm				

## VWR® Reversible Microtube Racks

### Removable Hinged Lid Fits Both Sides

- Side of the rack holds 8 or 12 x 0.2mL tube strips or individual tubes
- Reverse side holds 0.5mL tubes
- Both sides of the rack have 12 wells for 1.5mL tubes
- Autoclavable

Tubes	Wells	Wells diameter	L x W x H	Quantity	Color	Cat. No.
0.2 mL	168	5.5 mm	20.6 x 13.1 x 5.4 cm (8.11 x 5.16 x 2.13 ")	5/pk	Blue, Green, Pink, Yellow, Orange	82024-454
0.5 mL	40	7.6 mm	20.6 x 13.1 x 5.4 cm (8.11 x 5.16 x 2.13 ")	5/pk		
1.5/2.0 mL	12	11.4 mm	20.6 x 13.1 x 5.4 cm (8.11 x 5.16 x 2.13 ")	5/pk		



## VWR® Storage Racks for 80-Well Microtube Racks

### Strong, Clear Acrylic Construction

These racks help you keep large numbers of microtubes organized on your benchtop.

- Holds up to eight 80-well microtube racks filled with either 1.5 or 2.0mL microtubes
- Plenty of room to slide the racks in and out

Description	L x W x H	Cat. No.
Vertical rack	23.5 x 15.6 x 23.5 cm (9¼ x 6⅛ x 9¼ ")	82024-482
Horizontal rack	31.1 x 23.5 x 12.4 cm (12¼ x 9¼ x 4⅞ ")	82024-480




# TUBE RACKS

## VWR® Cube Racks

### Interlocking Multi-Tube Racks

- Made of polypropylene
- Allows you to work with different sized tubes all in the same plane
- Autoclavable



Tubes	Wells	Wells diameter	L x W x H	Quantity	Color	Cat. No.
50 mL	4	30 mm	10.8 x 11.3 x 11.3 cm (4.25 x 4.45 x 4.45 ")	5/pk	Blue, Green, Pink, Yellow, Orange	82010-746
15 mL	10	17 mm				
0.5 mL	12	11.5 mm				
1.5/2.0 mL	16	12.7 mm				

## VWR® 96-Place Microtube Racks





### Easy-to-Grasp Tabs for Sample Handling and Storage

Reverse 96 wells for 0.5mL tubes on one side and 96 wells for 1.5/2.0mL tubes on the reverse side.

- Made of polypropylene
- Wells are arranged in standard 8 x 12 array
- Tabs on each end of rack snap lid into place
- 0.5mL wells separated by molded-in grid lines, making them easier to see
- Raised rings around the 1.5 to 2.0mL wells make tubes easier to remove
- Autoclavable



Tubes	Wells	Wells diameter	L x W x H	Quantity	Color	Cat. No.
1.5/2.0 mL 0.5 mL	96	11 / 7.6	24.6 x 12.1 x 5 cm (9.6 x 4.7 x 2 ")	5/pk	Blue, Green, Pink, Yellow, Orange	82024-494
1.5/2.0 mL 0.5 mL	96	11 / 7.6	24.6 x 12.1 x 5 cm (9.6 x 4.7 x 2 ")	5/pk	Natural	82024-484
1.5/2.0 mL 0.5 mL	96	11 / 7.6	24.6 x 12.1 x 5 cm (9.6 x 4.7 x 2 ")	5/pk	Blue	82024-488
1.5/2.0 mL 0.5 mL	96	11 / 7.6	24.6 x 12.1 x 5 cm (9.6 x 4.7 x 2 ")	5/pk	Green	82024-490
1.5/2.0 mL 0.5 mL	96	11 / 7.6	24.6 x 12.1 x 5 cm (9.6 x 4.7 x 2 ")	5/pk	Pink	82024-492
1.5/2.0 mL 0.5 mL	96	11 / 7.6	24.6 x 12.1 x 5 cm (9.6 x 4.7 x 2 ")	5/pk	Yellow	82024-486





## VWR® Four-Way Microtube Racks

Link Together to Cater for Your Application Needs

Create as large a rack as you require for your bench top applications.

- Made from polypropylene
- Autoclavable

Tubes	Wells	Wells diameter	L x W x H	Quantity	Color	Cat. No.
50 mL	4	30 mm	17.4 x 9.5 x 5.2 cm (6.85 x 3.74 x 2.04 ")	5/pk	Clear	82010-752
15 mL	12	17.5 mm				
0.5 mL	32	8.2 mm				
1.5/2.0 mL	32	11.6 mm				
50 mL	4	30 mm	17.4 x 9.5 x 5.2 cm (6.85 x 3.74 x 2.04 ")	5/pk	Green	82024-460
15 mL	12	17.5 mm				
0.5 mL	32	8.2 mm				
1.5/2.0 mL	32	11.6 mm				
50 mL	4	30 mm	17.4 x 9.5 x 5.2 cm (6.85 x 3.74 x 2.04 ")	5/pk	Pink	82024-462
15 mL	12	17.5 mm				
0.5 mL	32	8.2 mm				
1.5/2.0 mL	32	11.6 mm				
50 mL	4	30 mm	17.4 x 9.5 x 5.2 cm (6.85 x 3.74 x 2.04 ")	5/pk	Yellow	82024-464
15 mL	12	17.5 mm				
0.5 mL	32	8.2 mm				
1.5/2.0 mL	32	11.6 mm				
50 mL	4	30 mm	17.4 x 9.5 x 5.2 cm (6.85 x 3.74 x 2.04 ")	5/pk	Orange	82024-458
15 mL	12	17.5 mm				
0.5 mL	32	8.2 mm				
1.5/2.0 mL	32	11.6 mm				
50 mL	4	30 mm	17.4 x 9.5 x 5.2 cm (6.85 x 3.74 x 2.04 ")	5/pk	Blue, Green, Pink, Yellow, Orange	82010-750
15 mL	12	17.5 mm				
0.5 mL	32	8.2 mm				
1.5/2.0 mL	32	11.6 mm				




## VWR® 80-Well Microtube Racks

Popular Fraction-Collector-Plate Style Rack

- Made of polypropylene
- Easily fill and transfer 80 tubes in compact arrangement compact arrangement (5 x 16 array)
- Molded-in coordinates to aid organization
- Autoclavable

Tubes	Wells	Wells diameter	L x W x H	Quantity	Color	Cat. No.
1.5/2.0 mL	80	11 mm	22.5 x 6.7 x 2.8 cm (8.9 x 2.6 x 1 ")	5/pk	Clear	82024-467
1.5/2.0 mL	80	11 mm	22.5 x 6.7 x 2.8 cm (8.9 x 2.6 x 1 ")	5/pk	Blue	82024-473
1.5/2.0 mL	80	11 mm	22.5 x 6.7 x 2.8 cm (8.9 x 2.6 x 1 ")	5/pk	Green	82024-475
1.5/2.0 mL	80	11 mm	22.5 x 6.7 x 2.8 cm (8.9 x 2.6 x 1 ")	5/pk	Pink	82024-477
1.5/2.0 mL	80	11 mm	22.5 x 6.7 x 2.8 cm (8.9 x 2.6 x 1 ")	5/pk	Yellow	82024-471
1.5/2.0 mL	80	11 mm	22.5 x 6.7 x 2.8 cm (8.9 x 2.6 x 1 ")	5/pk	Orange	82024-469
1.5/2.0 mL	80	11 mm	22.5 x 6.7 x 2.8 cm (8.9 x 2.6 x 1 ")	5/pk	Blue, Green, Pink, Yellow, Orange	82024-479



# TUBE RACKS

## VWR® 100-Well EPS Tube Racks

### High Storage, Minimal Footprint

- Two-piece racks are constructed of expanded polystyrene (eps) foam
- Lightweight microtube racks are ideal for long-term studies or tube storage
- Extra spacing between wells allows for easy gripping with fingertips
- Racks stack securely with nesting features on lid and base

Tubes	Wells	Wells diameter	L x W x H	Quantity	Color	Cat. No.
1.5/2.0 mL	50	11.5 mm	21.0 x 11.0 x 7.1 cm (8.26 x 4.33 x 2.7 ")	10/pk	White	89194-766
1.5/2.0 mL	100	11.5 mm	33.6 x 9.5 x 7.3 cm (13.22 x 3.74 x 2.87 ")	10/pk	White	82017-624



## VWR® Floating Microtube Racks

- Made from polypropylene
- Use on the bench, with water baths or ice containers
- Keeps microtubes afloat in a water bath or can be pushed into crushed ice - perfect for controlling sample temperature.
- Not recommended for autoclaving

Detachable legs provide two options:

- Attach legs so rack stands upright on bench top
- Detach legs so samples can be quickly popped out by pushing down onto any surface

Tubes	Wells	Wells diameter	Dimensions	Quantity	Color	Cat. No.
1.5/2.0 mL	16	10.9 mm	10.2 cm (square 4 ")	4/pk	White	82024-504
1.5/2.0 mL	16	10.9 mm	10.2 cm (square 4 ")	4/pk	Black	82024-506
1.5/2.0 mL	20	10.9 mm	9.8 cm dia. (round 3.8 ")	4/pk	White	82024-510




Detachable legs

## VWR® Floating Foam Tube Racks

- Made from high density polyethylene
- Reusable and flexible for a variety of applications
- Brightly colored, high density polyethylene foam racks are ideal for floating tubes in water baths, tubs, or beakers
- Closed-cell foam resists water absorption so racks can be washed and reused
- Racks include detachable carrying handles—easily lift samples out of baths
- The blue round rack features unique "x" slits to hold a variety of tube sizes
- Not recommended for autoclaving

Tubes	Wells	Well diameter	Size	Quantity	Color	Cat. No.
50 mL	4	29 mm	16.5 x 13.7 x 6.5 cm (6.4 x 5.3 x 2.5 ")	5/pk	Green	82017-628
15 mL	8	16 mm	17.2 x 5.7 x 6.5 cm (6.7 x 2.2 x 2.5 ")	5/pk	Blue	82017-630
1.5-2.0 mL	24	9.3 mm	15 x 8.3 x 7 cm (6 x 3.2 x 2.7 ")	5/pk	Yellow	82017-632
0.2/0.5 1.5-2.0 mL	18	9.3 mm	7 cm dia. (2.7 ")	5/pk	Blue	82017-634



Detachable handle


## VWR® Peg Tube Racks

Strong, Tough, Versatile Racks for Bench or Tube Storage

Fiberglass reinforced polypropylene racks are ideal for bench work or tube storage.

- Made from polypropylene
- Spacious peg design allows easy viewing of tube contents
- Use as a drying rack for glass plates or upside-down tubes
- Pegs remain straight even after repeated autoclaving
- Embossed indexing, ribbed grips, and drainage holes
- Autoclavable

Tubes	Wells	Wells diameter	L x W x H	Quantity	Color	Cat. No.
10 mL	50	14–17 mm	18.8 x 10.4 x 6.9 cm (7.4 x 4.1 x 2.71 ")	2/pk	Blue	82024-496



## VWR® Transformer™ Cube

For Storage of Large Samples and 15 and 50mL Tubes

Fiberglass reinforced polypropylene racks are ideal for bench work or tube storage.

- Polypropylene construction complete with two snap-in grids
- Imprinted grid and first tube orientation and keyed lid
- Use without grids for larger sized product storage
- Fill with ice and use as a temporary low temperature workstation
- Clear lid, translucent sides
- Autoclavable

Tubes	Wells	Wells diameter	L x W x H	Quantity	Color	Cat. No.
15 mL	16	18 mm	12.5 x 12.5 x 12.9 cm (4.9 x 4.9 x 5 ")	5/pk	Purple	89140-266
50 mL	9	31 mm				




Holds either:  
15mL or 50mL samples

## VWR® Rota-Rack® Tube Racks

Get Organized with this Fun, Interlocked Rack

- Made from polypropylene
- Rotate each of the four modules independently to configure for the tubes that you want to use
- Suitable for a massive variety of culture, centrifuge, micro centrifuge, round, flat or conical bottomed tubes
- Two versions to accommodate small and large tubes
- Autoclavable

Tubes	Wells	Wells diameter	L x W x H	Quantity	Color	Cat. No.
50 mL	2	29 mm	22.4 x 11 x 11 cm (8.81 x 4.33 x 4.33 ")	1/pk	Purple, Blue, Green, Yellow	89031-458
10-15 mL	8	20 mm				
15-20 mL	6	16.0 mm				
10-13 mL	10	13.5 mm				
15 mL	6	16.4 mm	21 x 10 x 10 cm (8.2 x 3.9 x 3.9 ")	1/pk	Purple, Blue, Green, Yellow	89031-460
1.5/2.0 mL	9	11.2 mm				
.5 mL	12	8.0 mm				
.2 mL	32	6.5 mm				



Twist and turn  
to suit your application

# TUBE RACKS

## VWR® Universal Tube Rack

Ideal for Holding a Wide Range of Tubes

- Made from polypropylene
- Holds up to 12 similar size tubes
- Imprinted references for easy identification
- Side handles for easy transport
- Not recommended for autoclaving



Adjustable slot adapts to different diameter tubes

Tubes	Wells	Wells diameter	L x W x H	Quantity	Color	Cat. No.
15-50 mL	12	16-30 mm	19.4 x 14 x 5.1 cm (7.6 x 5.5 x 2 ")	1/each	Purple/White	89234-744

## VWR® Pop-Up Tube Racks

Collapse for Compact Storage and "Pop-Up" for Quick Use

Fiberglass reinforced polypropylene racks are ideal for bench work or tube storage.

- Made from polypropylene
- Easy to open - easy to close
- Designed to help organize 15 and 50mL tubes
- Imprinted grid aids quick retrieval
- Stack racks in open or closed position
- Saves space when racks are not in use
- Not recommended for cleaning in autoclave



Tubes	Wells	Wells diameter	L x W x H (expanded)	L x W x H (collapsed)	Quantity	Color	Cat. No.
15 & 50 mL	21 / 12	17 / 30 mm	25.5 x 13.7 x 7.2 cm (10 x 5.39 x 2.83 ")	25.5 x 13.7 x 2.1 cm (10 x 5.39 x 0.83 ")	2/pk	Blue	89140-268
15 & 50 mL	21 / 12	17 / 30 mm	25.5 x 13.7 x 7.2 cm (10 x 5.39 x 2.83 ")	25.5 x 13.7 x 2.1 cm (10 x 5.39 x 0.83 ")	2/pk	Pink	89140-270



Easy to close



Easy to open



## VWR® Puzzle Tube Rack

Fun and Functional Racks Slide to Connect

- Made from polypropylene
- Two-sided puzzle-shaped pieces hold an assortment of 1.5/2.0, 15 and 50mL test tubes
- Pieces can be used separately for small sample applications or connect in a square or row configurations
- Aid with easy retrieval of tubes
- Autoclavable



Racks link together

Tubes	Wells	Wells diameter	L x W x H (expanded)	Quantity	Color	Cat. No.
1.5/ 2.0 mL	8/4	11.7 mm	15.1 x 13.5 x 6.6 cm (5.9 x 5.31 x 2.59 ")	4/pk	Yellow, Blue, Green, Orange	89234-740
15 mL	9	16.9 mm				
50 mL	4	29.5 mm				

## VWR® Drying Rack

Saves Glassware from Chips and Cracks

Drain and dry test tubes, beakers, flasks, and other labware.

- Place self-standing rack on a table or near a sink, or use hooks to hang rack on the wall
- Holds up to 90 pieces of labware on 18 x 127mm pins and 72 x 76mm pins
- Constructed of spot-welded, vinyl-coated steel wire for durability
- Easily accessible and removable bottom drainage pan

L x W x H	Quantity	Cat. No.
462 x 182 x 525 mm	1/each	82024-450



Removable pan helps remove excess water





## VWR® Universal Pipette Carousel Stand

Pipettes Never Looked so Good On Your Bench Top

NEW

Hold and store pipettes in style with the new universal pipette stand that holds up to 8 single channel, 4 multi channel, or a combination of both pipettes. A sleek and modern design enables this stand to securely hold virtually every known brand of pipette in the marketplace making it truly 'universal' while adding a pop of color to a bench top. Inserts on the upper holders can be added or removed to accommodate the fit of different of pipette widths. One handed removal or replacement of pipettes is made easy with a rotating head and stable weighted base. Moving the stand while keeping pipettes in place is simple with a quick grip handle allowing a user to grab and go. The pipette stands come with additional upper leaves to enable multiple color coding configuration options as well as a change in color to your work space, for a bit more fun.

- Ideal for labs that do not wish to be tied into a specific brand of pipette, or already have multiple brands in use
- Weighted base allows for easy removal or replacement of pipettes, with one hand
- Pipettes are held in a vertical position to help prevent contamination
- Conveniently rotate stand head to select or store pipettes
- Quick grip handle to move pipette stand (loaded or unloaded)
- Assorted color options allow for color coding pipette location
- The unit is not autoclavable but can be cleaned with soap/water or common laboratory disinfectants
- Made of ABS (Acrylonitrile Butadiene Styrene)
- Exclusive design



Inserts on the upper holders can be added or removed to accommodate the fit of different pipette widths

H x W	Quantity	Color	Lower holder	Upper holder	Cat. No.
Base 15.8 W cm; Top 22.9 W cm; 40.8 H cm (Base 6.2 W "; Top 9 "; Height 16.1 ")	1/each	White/Purple	Purple	Purple/Clear (8 of each)	75870-966

## VWR® Pipette Holders

Create an Instant Pipette Holding Point Within Your Lab

Secures to tabletops, benches and shelves with a twist.

- Chemical-resistant polypropylene and fits most major pipette brands
- Available in a single or triple format
- Autoclavable

L x W x H	Quantity	Color	Cat. No.
157 x 120 x 95 mm (6¼ x 4 ")	1/each	Blue	82007-226



Secures to tabletops, benches and shelves



## VWR® Acrylic Pipette Stands

Convenient and Safe Storage

3, 4 and 6-place pipette stands are designed to hold most major pipette brands.

Material	Holds	Slot upper	Slot lower	Color	Cat. No.
Acrylic	Three	25 mm	20 mm	Clear	82024-540
Acrylic	Four	25 mm	20 mm	Clear	82024-542
Acrylic	Six	25 mm	20 mm	Clear	82024-544



3-Place Stand

## VWR® 4-Place Pipette Racks

Protects and Organizes Glass and Plastic Pipettes

Angled four shelved compartments keep your 0.1 to 50mL pipettes or cans within easy reach while preventing damage and clutter.

- Heavy duty acrylic version provides crystal-clear viewing from any angle
- ABS version has cut-out viewing window and come with two removable attachments that hold pipette bulbs, fillers, and pumps
- Magnetized ABS version also has strong magnets on both sides and mounts securely to most metal surfaces

Material	L x W x H	Quantity	Color	Cat. No.
Acrylic	291 x 88 x 401 mm (3½ x 11½ x 16 ")	1/each	Clear	82024-536
ABS	291 x 88 x 401 mm (3½ x 11½ x 16 ")	1/each	Blue	82007-220



ABS version includes 2 removable pipette attachments




## VWR® Cargo Manual Pipet Rack

Customize Your Storage of Pipettes


Organize the cargo pipette rack to meet your laboratory or personal needs. Click together horizontally or vertically to build the desired configuration, or use as individual holders, utilizing the three lid sections provided.

- Magnets on both sides of units connect securely and safely to most metal surfaces
- Angled shelves store 0.1 to 50mL pipettes
- Protective lid encloses the unit to prevent dust build up and provides additional shelf space

Material	L x W x H	Quantity	Color	Cat. No.
Acrylic	301 x 100 x 127.5 mm (3¾ x 11⅞ x 5 ")	3/pk	Blue	89140-276



Magnets connect securely to most metallic surfaces



Units lock securely together horizontally or vertically

## PIPETTE STORAGE, BULBS

### VWR® Pipet Tip Disposal Box

#### Quick Flip of Port Cover Makes Tip Ejection Easy

Convenient for one-handed disposal of used pipette tips. Eliminates any need to handle pipette tips used with hazardous fluids.

- The 38" (10mm) thick acrylic provides protection from beta radiation
- The large, hinged lid allows easy disposal of contents
- Suggest using in conjunction with small biohazard disposal bags (page 81), if hazardous fluids are not acrylic friendly



L x W x H	Quantity	Cat. No.
154 x 171 x 206 mm (6 x 6 x 6 ")	1/pk	82010-754

### VWR® Small Latex Bulbs

- Designed for use with pasteur pipettes
- Small pipettes, and medicine droppers

Capacity	Length	OD	Hole diameter	Quantity	Color	Cat. No.
1 mL	43 mm	13 mm	5 mm	72/pk	Natural	82024-550
2 mL	52 mm	13 mm	6 mm	72/pk	Natural	82024-554



### VWR® Large-Volume, Heavy-Walled Bulbs

- Contoured shape and heavy-wall neck increase grip for easier use with pipettes and syringes

Capacity	Length	OD	Hole diameter	Quantity	Color	Cat. No.
12 mL	51 mm	34 mm	4 mm	1/each	Red	82024-564
15 mL	56 mm	34 mm	4 mm	1/each	Red	82024-566
30 mL	66 mm	46 mm	6 mm	1/each	Red	82024-568
60 mL	80 mm	54 mm	6 mm	1/each	Red	82024-570
150 mL	103 mm	73 mm	6 mm	1/each	Red	82024-572



### VWR® Vacuum Bulb, Synthetic Rubber

- Ideal for gas sampling, or for pumping gas through apparatus
- Two rubber valves with pvc connections


Capacity	Length	OD	Hole diameter	Quantity	Color	Cat. No.
85 mL	137 mm	53 mm	7.4 mm	1/each	Black	82024-578



### VWR® Safety Bulb Pipet Filler, Natural Rubber

- Features a white ABS attachment containing a silicone conical adapter that seals around pipettes without twisting
- Elastic band for hanging on reagent bottle or wrist
- Filler fits standard size pipettes

Capacity	Length	OD	Hole diameter	Quantity	Color	Cat. No.
15 mL	106 mm	41 mm	4.5-6.5 mm	3/each	Red	82024-574





## VWR® Premium Microscope Slide Boxes

Slide boxes are perfect for storage or transportation of standard size slides (25 x 75 mm and 1 x 3 ")

- Molded, rigid boxes securely hold slides in grooved slots
- Stackable design increases transit and storage convenience
- Both the lining and the slots are numbered to correspond to the slide inventory sheet



## VWR® 25-Place Microscope Slide Boxes

Capacity	Lining	Hinge	L x W x H	Quantity	Color	Cat. No.
25 slides	Cork	Friction fit	141 x 88 x 35 mm	72/case	Blue	82003-418
25 slides	Cork	Friction fit	141 x 88 x 35 mm	72/case	Green	82003-420
25 slides	Cork	Friction fit	141 x 88 x 35 mm	72/case	Red	82003-422
25 slides	Cork	Friction fit	141 x 88 x 35 mm	72/case	White	82003-426
25 slides	Cork	Friction fit	141 x 88 x 35 mm	72/case	Yellow	82003-424



Easy access slides with friction-fit cover

## VWR® 50-Place Microscope Slide Boxes

Capacity	Lining	Hinge	L x W x H	Quantity	Color	Cat. No.
50 slides	Cork	Thumb-latch	209 x 86 x 35 mm	72/case	Blue	82024-580
50 slides	Cork	Thumb-latch	209 x 86 x 35 mm	72/case	Red	82024-582
50 slides	Cork	Thumb-latch	209 x 86 x 35 mm	72/case	White	82024-584
50 slides	Cork	Thumb-latch	209 x 86 x 35 mm	72/case	Green	82024-588



Thumb-latch lock mechanism

## VWR® 100-Place Microscope Slide Boxes

Capacity	Lining	Hinge	L x W x H	Quantity	Color	Cat. No.
100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	72/case	Blue	82003-406
100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	72/case	Green	82003-408
100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	72/case	Red	82003-410
100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	72/case	Yellow	82003-412
100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	72/case	White	82003-414
100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	72/case	Gray	82003-416



Closes securely with a rust-resistant, nickel-plated clasp and pin hinges

## VWR® Tall 100-Place Microscope Slide Box

Capacity	Lining	Hinge	L x W x H	Quantity	Color	Cat. No.
100 slides	Cork	Nickle clasp	208 x 175 x 60 mm	72/case	Blue	89140-280



Stackable, holds 100-place slides with additional space for larger 3 x 2 in.

## VWR® Economy Microscope Slide Boxes

Constructed of Durable Polypropylene

Slide boxes hold either 25 or 100 standard slides (25 x 75mm and 1 x 3").

- Hinged lids are held securely closed with a latch lock
- Stacking feet maximize stability
- Slide identification card

Capacity	Hinge	L x W x H	Quantity	Color	Cat. No.
25 slides	Latch lock	141 x 92 x 30 mm	72/case 1/each 24/pk	Blue	82024-614
25 slides	Latch lock	141 x 92 x 30 mm	72/case 1/each 24/pk	Red	82024-616
25 slides	Latch lock	141 x 92 x 30 mm	72/case 1/each 24/pk	White	82024-612
100 slides	Latch lock	210 x 169 x 37 mm	72/case 1/each 12/pk	Blue	82024-608
100 slides	Latch lock	210 x 169 x 37 mm	72/case 1/each 12/pk	Red	82024-610
100 slides	Latch lock	210 x 169 x 37 mm	72/case 1/each 12/pk	White	82024-606



## VWR® Slide Rack for Slide Boxes

Create a File System for 50 and 100 Place Slide Boxes

Storage rack makes organization of samples easy.

- Clear acrylic construction allows the labels to be seen on the slide boxes
- Place boxes in rows or columns by color to create an easy coding system

L x W x H	Bins	Quantity	Cat. No.
711 x 171 x 228 mm (6¾ x 30 x 9 ")	15 bins	1/each	82024-620



## VWR® 20-Place Slide Tray Racks

Search, Identify, and Easily Locate the Specimen Slide You Need


Store multiple slide trays for easy retrieval - ideal for microbiology, histology and pathology fields.

- ABS rack holds five 20-place slide trays (not included see above)
- Eliminates the worry of breakage during retrieval - easy access to each of the trays
- Matching clasps provide area for identification labels and doubles as a handle
- Portable and refrigerator safe for cold sample storage
- Stack multiple units to create your own library

Capacity	L x W x H	Quantity	Color	Cat. No.
5 slide trays	226 x 340 x 125 mm (8¾ x 13¾ x 5 ")	1/each	Blue	89140-274



Includes 5 colored clasps



Stackable

## VWR® Micro Slide Trays

Cardboard Slide Trays Hold Standard Slides  
(25 x 75mm and 1 x 3")

- Horizontal compartments are recessed to protect the specimens when the cover is closed
- Features thumb cut for easy slide removal
- Available in 2 sizes with 5 color options for the 20-place trays
- Disposable

Capacity	L x W x H	Quantity	Color	Cat. No.
20 slides	208 x 341 x 9 mm (8 <sup>3</sup> / <sub>16</sub> x 13 <sup>3</sup> / <sub>16</sub> x 3 <sup>1</sup> / <sub>2</sub> ")	72/case	Black	82020-913
30 slides	281 x 328 x 10 mm (11 <sup>1</sup> / <sub>16</sub> x 12 <sup>1</sup> / <sub>16</sub> x 3 <sup>1</sup> / <sub>16</sub> ")	40/case	Black	82020-915




20 and 30-place

## VWR® Polyethylene Slide Holders

Hold Slides During Evaluation, Cleaning and Drying

- Made from chemical-resistant polyethylene
- Solid, strong construction with rubber feet for bench top stability

Capacity	L x W x H	Quantity	Cat. No.
15 slides	153 x 102 x 21 mm (6 x 4 x 5/8 ")	1/each	82024-524
40 slides	298 x 105 x 24 mm (11 <sup>3</sup> / <sub>4</sub> x 4 <sup>1</sup> / <sub>8</sub> x 1 ")	1/each	82024-526
78 slides	299 x 178 x 21 mm (11 <sup>3</sup> / <sub>4</sub> x 7 x 3/4 ")	1/each	82024-528



## VWR® Color-Coded Two-Place Slide Mailers

- Polypropylene protectors hold two standard (25 x 75mm and 1 x 3") slides securely
- Large, easy-to-secure snaps, provide safe storage for valuable slides
- Variety of colors for coding
- Disposable

Capacity	L x W x H	Quantity	Color	Cat. No.
2 place	93 x 68 x 5 mm	25/pk	Red, Yellow, Blue, Purple, Orange	82007-222
2 place	93 x 68 x 5 mm	25/pk	Natural	82007-224



Ideal for secure shipping and storing of slides

## VWR® 5-Place Slide Mailers

- Polypropylene containers hold five standard slides (25 x 75mm and 1 x 3")
- Inside wall features slotted channels to keep slides separated and to make removal easier
- Lid closes securely to keep slides safe
- Disposable

Capacity	L x W x H	Quantity	Color	Cat. No.
5 place	81 x 19 x 31 mm	25/pk	Natural	82024-598
5 place	43 x 24 x 88 mm	25/pk	Red, Yellow, Blue, Purple, Orange	82003-431



Side opening

End opening

Ideal for secure shipping, storing, or staining slides

# HISTOLOGY & MICROSCOPY STORAGE

## VWR® Freezer Racks

Corrosion-resistant, stainless steel horizontal racks feature handles for easy transport.

- Stainless steel fabrication
- Multiple shelf counts

Description	Config.	Box size	L x W x H	Quantity	Cat. No.
Four-Shelf Rack, 2" Boxes	4 x 4	5.3 cm (2")	55.6 x 14 x 22.5 cm (21.8 x 5.5 x 8.9")	1/each	82024-520
Three-Shelf Rack, 3" Boxes	3 x 4	7.6 cm (3")	55.6 x 1.40 x 2.45 cm (22 x 5.5 x 8.9")	1/each	82024-518



## VWR® Wire Petri Dish Holder

- Holds up to 6 x 100mm petri dishes
- Epoxy-coated steel wire
- Autoclavable

Description	L x W x H	Color	Quantity	Cat. No.
Wire Petri Dish Holder	127 x 86 x 107 mm	White	1/each	82008-018



## VWR® Drawer Organizer

Keeps Supplies Organized and Easily Accessible

Polystyrene organizer is perfect for storing and organizing small laboratory supplies.

- Protect small and fragile items
- Fit in most standard drawers or store on benches and shelves

Compartments	Overall Dimensions	Quantity	Cat. No.
3-Compartment Tray:	46 x 30.5 x 6.5 cm (18 x 12 x 2½")	10/case	82024-548
2 compartments measure 430 x 80 mm			
1 compartment measures 430 x 110 mm			



## VWR® Disposable Slide Mailers

Disposable, Safe, Secure Transport of Microscope Slides

Slide boxes hold either 25 or 100 standard slides (25 x 75mm and 1 x 3").

- Thumb grooves easy slide removal
- 1- and 2-place slide mailer versions
- Each will hold 25 x 75mm (1 x 3") and or 26 x 76mm sized slides
- Cardboard construction

### Cardboard slide mailers

Capacity	L x W x H	Quantity	Color	Cat. No.
1 place	95 x 42 x 5 mm (1½ x 3/16 x 3¼ x 3¹¹/₁₆")	25/pk	Natural	82020-908
2 place	103 x 79 x 5 mm (3½ x 3/16 x 4¼")	36/pk	Natural	82020-911



1-place

2-place

Heavy cardboard mailers for safe transport



## VWR® Disposable Square Weighing Boats

Flexible Polystyrene Boats Feature Smooth Surfaces to Provide Accurate Pour-Outs with Minimal Sample Loss

- Select white, black, or blue for best color contrast to ensure sample is completely removed
- Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases
- Temperature range -10 to 70°C (14 to 158°F)
- Can be used with microwave moisture analyzers
- All are FDA/food-grade compliant
- Boats are shallow with a wide, flat bottom to resist tipping
- Rounded corners simplify transfer



Smooth surfaces provide accurate pour-outs



Manufactured and packaged in a clean room

Capacity	L x W x H	Color	Quantity	Standard Cat. No.	Antistatic Cat. No.
7 mL	46 x 46 x 8 mm (1.8 x 1.8 x 0.3 ")	White	500/pk	10803-136	10770-442
7 mL	46 x 46 x 8 mm (1.8 x 1.8 x 0.3 ")	Black	500/pk	10803-138	10803-176
7 mL	85 x 85 x 24 mm (3.3 x 3.3 x 1 ")	Blue	500/pk	10803-116	10803-162
100 mL	85 x 85 x 24 mm (3.3 x 3.3 x 1 ")	White	500/pk	10803-148	10770-448
100 mL	85 x 85 x 24 mm (3.3 x 3.3 x 1 ")	Black	500/pk	10803-146	10803-164
100 mL	85 x 85 x 24 mm (3.3 x 3.3 x 1 ")	Blue	500/pk	10803-150	10770-406
250 mL	140 x 140 x 22 mm (5.5 x 5.5 x 0.9 ")	White	500/pk	10803-170	10770-450
250 mL	140 x 140 x 22 mm (5.5 x 5.5 x 0.9 ")	Black	500/pk	10803-172	10803-174
250 mL	140 x 140 x 22 mm (5.5 x 5.5 x 0.9 ")	Blue	500/pk	10803-168	10803-152



## VWR® Disposable Hexagonal Weighing Boats

- Bend flexible boats at points to form pour spouts
- Flat bottoms and sloping sides provide accurate transfer

Size	Base I.D.	Top I.D.	Color	Quantity	Cat. No.
Small	25 mm (1 ")	38 mm (1 ½ ")	White	500/pk	10770-452
Medium	51 mm (2 ")	64 mm (2 ½ ")	White	500/pk	10770-454
Large	76 mm (3 ")	102 mm (4 ")	White	500/pk	10770-456
X-Large	89 mm (3 ½ ")	127 mm (5 ")	White	500/pk	10770-458



## VWR® Pour-Boat Weighing Dishes

- Pour-spout design easily dispenses samples without flexing the dish
- Flat-bottom dishes are made from 12-mil thick polystyrene for better strength, stability, and handling

Size	W x D x H	Color	Quantity	Standard Cat. No.	Antistatic Cat. No.
Small	43 x 58 x 13 mm (1.6 x 2.2 x 0.5 ")	White	500/pk	10803-140	10770-408
Medium	83 x 132 x 26 mm (3.2 x 5.1 x 1 ")	White	500/pk	10803-160	10770-438
Large	108 x 183 x 26 mm (4.2 x 7.2 x 1 ")	White	500/pk	10803-166	10770-440



## VWR® Disposable Weighing Boats

- Diamond shape that easily flexes for accurate pouring
- Wide, flat bottom to resist tipping
- Flexible boats bend to form pour spout

Capacity	W x D x H	Color	Quantity	Standard Cat. No.	Antistatic Cat. No.
5 mL	55 x 35 x 6 mm (2.1 x 1.3 x 0.2 ")	White	500/pk	10803-114	10803-134
5 mL	55 x 35 x 6 mm (2.1 x 1.3 x 0.2 ")	Black	500/pk	10803-110	10803-112
25 mL	80 x 50 x 14 mm (3.1 x 1.9 x 0.5 ")	White	500/pk	10803-124	10803-142
25 mL	80 x 50 x 14 mm (3.1 x 1.9 x 0.5 ")	Black	500/pk	10803-128	10803-120
30 mL	80 x 60 x 14 mm (3.1 x 2.3 x 0.5 ")	White	500/pk	10803-132	10803-144
30 mL	80 x 60 x 14 mm (3.1 x 2.3 x 0.5 ")	Black	500/pk	10803-130	10803-122
100 mL	125 x 100 x 19 mm (4.9 x 3.9 x 0.7 ")	White	500/pk	10803-156	10803-126
100 mL	125 x 100 x 19 mm (4.9 x 3.9 x 0.7 ")	Black	500/pk	10803-158	10803-118



# LAB ORGANIZATION

## VWR® Glove Box Holders

Keep Gloves Organized and Ready for Immediate Use

Epoxy-coated steel wire racks are an economical solution.

- Free up bench top space by mounting holders on wall for easy dispensing of gloves
- All holders fit most major brands of boxes
- Mounting screws included

Description	W x H x D	Material	Quantity	Cat. No.
1-box holder	140 x 209 x 108 mm (4 x 5½ x 8¾ ")	Wire	1/each	89011-726
3-box holder	250 x 454 x 108 mm (9¼ x 3¾ x 17 ")	Wire	1/each	89011-728



## VWR® Acrylic Glove Box Holders

Includes a polyester foam insert to securely hold a range of glove box sizes—snug fit ensures easy removal of gloves.

Description	W x H x D	Material	Quantity	Cat. No.
1-box holder	141 x 260 x 110 mm (5½ x 10¼ x 4½ ")	Acrylic	1/each	89011-720
2-box holder	295 x 260 x 110 mm (11¾ x 10¼ x 4½ ")	Acrylic	1/each	89011-722
3-box holder	295 x 385 x 110 mm (11¾ x 15¾ x 4½ ")	Acrylic	1/each	89011-724



## VWR® Dispensing Bins

Work Great with Gloved Hands, Bright Colors

Ergonomically designed with large, angled openings and oversized lip.

- Lid remains open when tilted back
- Nonskid rubber feet and pre-drilled holes for wall mounting
- Green bin includes a removable divider
- Mounting screws included

Compartments	L x W x H	Color	Quantity	Cat. No.
Single	155 x 155 x 170 mm (4½ x 6 x 6½ ")	Blue	1/case	82024-530
Dual	180 x 155 x 165 mm (7 x 6 x 6½ ")	Green	1/case	82024-532



## VWR® Acrylic Dispenser for Parafilm Sealing Film

Heavy duty Construction for Increased Durability

Store, dispense, and cut Parafilm®, tape, labels, and Tough-Tags®.

- Strong, crystal-clear acrylic dispenser makes your film supply visible and features a serrated cutting blade
- Dispenser holds one or two 50mm (2") rolls, or one 100mm (4") roll



Overall Dimensions	Quantity	Cat. No.
119 x 168 x 175 mm	1/each	82024-547

## VWR® Mini Bin™ Foil Dispenser

### Safely Cut Foil to Any Size

The Mini Bin foil dispenser is a great alternative to wasteful giant-sized foil rolls and expensive precut foil squares.

- Compact to fit in limited work space areas and easy to transport
- Safe to use - load the foil roll, press the cover to activate the safety guard
- Foam strip prevents foil retracting back into unit
- For use with 102mm or 153mm
- (4 or 6") length foil rolls (sold separately)



W x H x D	Quantity	Cat. No.
190 x 100 x 112 mm	1/each	89234-754

## VWR® Extreme Temperature Mitt

### Protect Your Hands From Hot or Cold Labware

The silicone rubber mitt is designed to comfortably fit your hand and they have internal ribs in the pockets for added dexterity.

- Oversized surfaces with nubs for secure grip
- Mitts withstand temperatures from -70 to 500°F (-57 to 260°C)
- Autoclavable



Description	L x W x H	Color	Quantity	Cat. No.
Silicone full-hand mitt	177 x 100 x 32 mm	White	1/each	89026-680

## VWR® Laboratory Notebooks

For recording research data, experimental procedures, and design ideas, use these notebooks to keep information safe and well organized

- Numbered, acid-free, pages include spaces for title, project number, book number, date, and witness signature
- Lined notebooks feature college-ruled pages for easy report writing
- Covers for all versions are waterproof and chemical-resistant
- All notebooks contain periodic and conversion tables in the back



Pages	Page format	L x W	Quantity	Color	Cat. No.
100	Grid	218 x 281 mm (8½ x 11 ")	12/case	Black	82007-219
200	Grid	218 x 281 mm (8½ x 11 ")	12/case	Blue	82007-233
200	Lined	218 x 281 mm (8½ x 11 ")	12/case	Red	82026-407



1.800.932.5000 | [vwr.com](http://vwr.com)

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.

8422 BB LM