

Advance Protection for Critical Environments.



Because critical environments call for the highest standard of stainless steel.

316 Stainless Steel Tabletop and undercarriage, 16 ga.

304 Stainless Steel posts and undershelves, 14 ga.

316 Tabletops can withstand corrosive attacks by sodium, calcium brines, hypochlorite solutions, phosphoric acids, sulfate and sulfuric acid. 316 stainless steel specifically resists chloride corrosion and can stand up to harsh chemicals which prevents specific types of corrosion and the typical pitting that can be caused by industrial solvents and cleaners.



Metro 316
Stainless Steel
Worktables

316 protects 24/7/365

Easy to assemble and are shipped knock down!

Shelf

WVR Cat. No.	Description
76422-052	24x36 Mobile Worktable Shelf
76422-062	24x48 Mobile Worktable Shelf
76422-070	24x60 Mobile Worktable Shelf
76422-106	30x48 Mobile Worktable Shelf
76422-084	30x60 Mobile Worktable Shelf
76422-108	30x72 Mobile Worktable Shelf
76422-094	30x96 Mobile Worktable Shelf
76422-064	24x36 Worktable Shelf
76422-068	24x48 Worktable Shelf
76422-058	24x60 Worktable Shelf
76422-082	30x48 Worktable Shelf
76422-098	30x60 Worktable Shelf
76422-078	30x72 Worktable Shelf
76422-074	30x96 Worktable Shelf



3-Sided Frame

WVR Cat. No.	Description
76422-072	24x36 Mobile Worktable 3-Sided Frame
76422-060	24x48 Mobile Worktable 3-Sided Frame
76422-050	24x60 Mobile Worktable 3-Sided Frame
76422-110	30x48 Mobile Worktable 3-Sided Frame
76422-076	30x60 Mobile Worktable 3-Sided Frame
76422-096	30x72 Mobile Worktable 3-Sided Frame
76422-102	30x96 Mobile Worktable 3-Sided Frame
76422-056	24x36 Worktable 3-Sided Frame
76422-066	24x48 Worktable 3-Sided Frame
76422-054	24x60 Worktable 3-Sided Frame
76422-086	30x48 Worktable 3-Sided Frame
76422-088	30x60 Worktable 3-Sided Frame
76422-092	30x72 Worktable 3-Sided Frame
76422-112	30x96 Worktable 3-Sided Frame



H Frame

WVR Cat. No.	Description
76422-080	30x60 Mobile Worktable H Frame
76422-114	30x72 Mobile Worktable H Frame
76422-100	30x96 Mobile Worktable H Frame
76422-104	30x60 Worktable H Frame
76422-090	30x72 Worktable H Frame
75422-116	30x96 Worktable H Frame

U.S. part numbers listed.



We put space to work.

LO9-028
© 2022 InterMetro Industries Corporation, Wilkes-Barre, PA 18705
0822 Lit. No. 130256W

avantor[™]
delivered by **VWR**[™]

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.