



Product Specification

Catalog Number: BDH9230

Catalog Description: BDH DEXTROSE ANHYDRS ACS

Grade: ACS

Date Current Specification Set: NOV. 2012

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics	Specifications
-----------------	----------------

Chloride	<= 0.01 %
Heavy Metals (as Pb)	<= 0.0005 %
Insolubles	<= 0.005 %
Iron	<= 0.0005 %
Loss on Drying	<= 0.2 %
Residue on Ignition	<= 0.02 %
Specific Rotation	+52.5 - +53.0 degrees
Starch	PASS
Sulfate & Sulfite	<= 0.005 %
Titrate Acid	<= 0.002 meq/g

Signed on behalf of VWR:

Robin Crocker
Quality Control Manager