



NOMINAL CAPACITY	50 mL	1.7 oz.
OVERFLOW CAPACITY	62 mL	2.1 oz.
MATERIALS	CAP	VIRGIN PFA
	BOTTLE	VIRGIN PFA
THREAD	M33SP400	
NOMINAL WEIGHT	CAP	12 g
	BOTTLE	23 g
	TOTAL	35 g
OPENING I.D.	1.01"	2.57 cm
NOMINAL WALL	.040"	.10 cm
MANUFACTURING ENVIRONMENT	CLASS 10,000 (ISO 7)	
BOTTLE MARKINGS	SAVILLEX 50 mL OTHER PFA	
CAP MARKINGS	SAVILLEX PFA	
PACKAGING	INDIVIDUALLY CAPPED AND LDPE BAGGED	
LABEL ON PACKAGING	PART NUMBER: 150-01-0050 QTY: 1 MFG. LOT: _____	
CAP TORQUE	HAND TIGHT	
USP CLASS VI CERTIFIED		YES
ADCF	CERTIFIED FREE OF ANIMAL DERIVED COMPONENTS	
AUTOCCLAVABLE: 121°C		YES

- NOTES:
1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND CENTIMETERS UNLESS OTHERWISE NOTED.
 2. THIS IS NOT A PRODUCTION DRAWING BUT IS FOR INFORMATIONAL PURPOSES ONLY AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

300682-A.SLDPRT



Purilllex

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF SAVILLEX CORPORATION. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF SAVILLEX CORPORATION IS PROHIBITED.

UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: FRACTIONAL ± 1/32" ANGULAR ± .5 DEG. TWO PLACE DECIMAL ± .020" THREE PLACE DECIMAL ± .005"		DRAWN	NAME	DATE
MATERIAL		ENG. APPR.		
FINISH		TOOL APPR.		
DO NOT SCALE DRAWING		MFG APPR.		
		Q.A.		
		COMMENTS	PART NUMBER: 150-01-0050	
		TITLE: 50mL PFA BOTTLE SPECIFICATIONS		
SCALE: 1:1	WEIGHT: LBS.	10321 West 70th St E5544-3446 USA (952)935-4100		
		SIZE DWG. NO. 300682-A	REV. A	SHEET 2 OF 3