



Item Number ZI120
Item Zinc Sulfate, Heptahydrate, Granular, USP
CAS Number 7446-20-0
Molecular Formula $ZnSO_4 \cdot 7H_2O$
Molecular Weight 287.56

Assay.....99.0-108.7%
Acidity..... To pass test

MAXIMUM LIMITS

Alkalies and Alkali Earths.....0.9%
Arsenic (As).....14 ppm
Lead (Pb).....0.002%
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843