



Puritan
Quality since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

An affiliate of Hardwood Products Company LLC







Tel 800-321-2313 (US and Canada)
207-876-3311

Fax 800-323-4153 (US and Canada)
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

DNA-FREE, RNASE-FREE & DNA CONTROLLED APPLICATIONS Forensic professionals rely on Puritan swabs for clean and accurate sampling. Now you can feel even more confident. Our collection of DNA-free and RNase-free sampling swabs have been tested by an independent lab. These tests are done to the most stringent standards (less than 25 pg of human DNA). The result? Certified to be DNA-free at a threshold lower than any other swabs in the market today. Our collection of DNA-controlled swabs are manufactured in a DNA controlled and monitored clean room environment.

<i>DNA / RNase Free</i>	Item No.	Description	Pack	Handle Dimensions	Tip Dimensions
	25-806 1WC FDNA	Cotton tip, wood handle, <i>sterile</i>	10/100/1	<u>5.937 length</u> 0.1 diameter	<u>0.625 length</u> 0.188 diameter
	25-806 2WC TT FDNA	Cotton tip, wood handle, transport tube, <i>sterile</i>	10/50/2	<u>5.937 length</u> 0.1 diameter	<u>0.625 length</u> 0.188 diameter
	25-1805 1PF RND FDNA	Large foam paddle tip, polypropylene handle, <i>sterile</i>	10/50/1	<u>5.125 length</u> 0.25 diameter	<u>1.125 length</u> 0.906 diameter
	25-3306-U TT FDNA	PurFlock Ultra flocked tip, plastic handle, transport tube, <i>sterile</i>	10/50/1	<u>6.04 length</u> 0.1 diameter	<u>0.69 length</u> 0.2 diameter
<i>DNA Controlled</i>	Item No.	Description	Pack	Handle Dimensions	Tip Dimensions
	25-805 1PC TT FABUSA	Cotton tip, polystyrene handle, special patented tube design for product breathability, <i>sterile</i>	10/50/1	<u>4.970 length</u> 0.1 diameter	<u>0.687 length</u> 0.203 diameter
	25-825 1WCS TT FABUSA	Pointed small cotton tip, wood handle, special patented tube design for product breathability, <i>sterile</i>	10/50/1	<u>5.032 length</u> 0.083 diameter	<u>0.5 length</u> 0.12 diameter at widest point

Products shown are not actual size.

