



## Buy 2 Cases of WHEATON Vial Files Receive 1 FREE Case

### Sample Vials in Vial File®

- Sample Vials in Vial File provides a compact, easily accessible means of sample storage
- Polystyrene case with partitions, foam inserts and alphanumerical indexing card organize and protect vials
- Clear vials are manufactured from low extractable borosilicate glass that conforms to USP Type I and ASTM E 438 Type I, Class A requirements
- Amber vials made from borosilicate glass that conforms to USP Type I requirements for light transmission to protect light-sensitive products
- Choice of solid or open top black phenolic screw cap
- Open top screw cap provides for sample retrieval with a syringe
- Caps attached to vials
- Vials and caps are autoclavable

Cat. No.	Color	Size (mL)	Liner Material	Cap Size	Dia. x H (mm)	Vials / File
Black Phenolic Solid Top Cap						
66011-653	Clear	2	14B Rubber	8-425	12 x 38	60
80076-608	Clear	2	PTFE / 14B Rubber	8-425	12 x 38	60
80076-600	Amber	2	14B Rubber	8-425	12 x 38	60
80076-612	Amber	2	PTFE / 14B Rubber	8-425	12 x 38	60
66011-675	Clear	4	14B Rubber	13-425	15 x 48	40
80076-610	Clear	4	PTFE / 14B Rubber	13-425	15 x 48	40
80076-602	Amber	4	14B Rubber	13-425	15 x 48	40
Black Phenolic Open Top Cap						
66014-506	Clear	1.5	PTFE / Silicone	8-425	12 x 35	60
80076-620	Amber	1.5	PTFE / Silicone	8-425	12 x 35	60
66011-755	Clear	4	PTFE / Silicone	13-425	15 x 48	40

Purchase your Wheaton Vial Files by visiting [vwr.com](http://vwr.com) between May 1 and December 31, 2017. To redeem your FREE product: Have your VWR invoice or packing slip ready, visit [vwr.com/promotions](http://vwr.com/promotions), enter the promo code **4671** in the search box and fill in the online form before January 31, 2018. Please allow up to 8 weeks for your free product.

Offer valid in North America. Limited to two redemptions per customer. Free case must be of equal or lesser value. Not for resale. Offer valid while supplies last. All purchases must be on one purchase order.



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.

0517 Lit. No. 230012WREV Promo Code 4671

VPRO\_009 05/2017