



Item Number P1457
Item Propylparaben Sodium, NF
CAS Number 35285-69-9
Molecular Formula $C_{10}H_{11}NaO_3$
Molecular Weight 202.00

Assay (Anhydrous Basis).....98.5-101.5%
Completeness of Solution..... To pass test
pH of a 1 in 1,000 Solution..... 9.5-10.5

MAXIMUM LIMITS

Water..... 5.0%
Chloride.....0.035%
Sulfate.....0.12%
Residual solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843