



order on
VWR.COM

M-Zone Modular™ Cleaning System



Introducing Micronova's new M-Zone Modular™ Cleaning System

In an effort to make the disinfection process more cost efficient, ergonomic – and simple – Micronova offers the **MegaWringer**, and **BucketBinder** Multi Bucket kit. This unique core cleaning solution replaces the need for bulky carts and frames. A multi bucket system supporting a super simple downward press wringer. Efficiently wrings flat head mops with minimum effort.

Allows for easier usage, cleaning and storage – in both small and large controlled environments and sterile core.

M-Zone™ Flat Mops

A lightweight, easy to use cleanroom mop which can be used - wet or dry - to clean and disinfect in ISO 5 to 9 areas.



- Mops glide easily onto adapters and can be removed by pulling the looped tab
- A variety of fabrics available for added absorbency and cleaning strength. Generous foam interior allows optimum dwell time of disinfectant.
- Fits on QPSL plastic or QDSL stainless steel frame both with 360° swivel. Frame also locks at 180° for precise, controlled cleaning.
- Compatible with strong disinfectants and steam sterilization. Available gamma-irradiated for aseptic areas

M-Zone Flat Mop: MZSM11-14S
VWR Cat. No. 89498-078

M-Zone Flat Mop – Irradiated: MZSM11-14SIR
VWR Cat. No. 89498-080



MegaWringer™

Lightweight and portable downward press wringer for consistent wringing of flat head mops. The wringer boasts an ergonomic, one hand action. The mop handle slides into the wringer grip, providing leverage to create optimum wringing pressure.



- Lightweight and durable – weighs less than 5 lbs.
- Molded Polyketone and electropolished Stainless Steel.
- Fits on single, double or triple bucket systems – no cart required.

MegaWringer: W-94 VWR Cat. No. 76299-864

BucketBinder™

Polyketone 'hooks' firmly secure two or three rectangular buckets together to instantly form a multi bucket system, which can be scaled up or down as required.

- Less than 2lbs per pair, BucketBinders hold up to strong disinfectant, can be autoclaved and irradiated.
- Fit the B-7 AC series rectangular buckets

Double: BB-DBL-DBYE
VWR Cat. No. 76289-876

Triple: BB-TR-DBREYE
VWR Cat. No. 76289-878




Bucket Liners

- An impermeable barrier/liner for buckets and storage and transport containers
- LDPE/Nylon laminate is excellent for disinfectants and cleaning solutions
- Ideal for parts and raw material transfer
- Available sterile.

Regular: BKL7-3119 VWR Cat. No. 89173-176
Irradiated: BKL7-3119IR VWR Cat. No. 68400-780

Start Your New Cleanroom with the Essentials

Micronova[®]
MANUFACTURING INC.



Hood Cleaning

Effective and efficient.

Isolator Tool, Plastic: QPIT-09
VWR Cat. No. 89140-698

Isolator Tool Cover, Polyester: ITC6-gNS
VWR Cat. No. 76018-368

Isolator Tool Cover, Sterile: ITC6-gNSIR
VWR Cat. No. 76018-370

Plastic Handle: STAH-12WH
VWR Cat. No. 89428-626



Hand Care Made Simple

Touch Free MicroDispenser: ICDISP-4
VWR Cat. No. 89128-812

Aquahol, Sterile Alcohol: IC-220
VWR Cat. No. 68400-440

NovaClenz, Ethanol Sanitizer: IC-210
VWR Cat. No. 68400-031



Simple Collection and Removal of Waste

BinLiner, 40 gallon, 2 mil: BL2-40TRBB
VWR Cat. No. 75838-314

BinLiner, Sterile Version: BL2-40TRIRB
VWR Cat. No. 89092-152

Other sizes and thicknesses available.



Lab Markers – Suite Ready

IPA Resistant Lab Marker: PEN-40
VWR Cat. No. 89205-944

Irradiated Lab Marker: PEN-40IR
VWR Cat. No. 89205-946

VWR[™]
part of avantor

800 932 5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.