



**Item Number** B1102  
**Item** Benzoic Acid, FCC  
**CAS Number** 65-85-0  
**Molecular Formula**  $C_7H_6O_2$   
**Molecular Weight** 122.12

Assay ( $C_7H_6O_2$ ; Anhydrous Basis)..... 99.5-100.5%  
Solidification Point..... 121°-123°C

MAXIMUM LIMITS

Lead..... 2.0 mg/kg  
Readily Carbonizable Substances..... To pass test  
Readily Oxidizable Substances..... To pass test  
Residue on Ignition..... 0.05%  
Water..... 0.7%

***Spectrum Chemicals & Laboratory Products***

---

**Corporate Headquarters:**    **East Coast Facility:**  
14422 S. San Pedro St.    755 Jersey Ave.  
Gardena, CA 90248    New Brunswick, NJ 08901  
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843