

Product Name:	HCUCFS-0404
Description:	<p>VWR Catalog No. 10819-650  VWR Undercounter Freestanding Refrigerator  5 Cu. Ft. Capacity  One solid swing door, right hinged  Mechanical thermostat  Adjustable operating temperature range: 1°C – 10°C  Three shelves (2 adj/1 fixed)  Forced draft air circulation  Keyed door lock  Magnetic door gasket for positive seal  Recessed door handle  Leveling legs  Interior Light  Probe access port (3/8")  White cabinet and door  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Two year parts and labor warranty, plus an additional three year compressor parts only warranty  Exterior dimensions: 20 1/4"W x 22" D x 31 1/2"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Manual defrost  115V, 60 Hz, .90 Amps, 1/7 HP  UL/C-UL listed  Shipping weight: 90 lbs.</p>