

# Dry Wipes

(Cleanroom & Industrial)



## Cleanroom Wipes

Wipe name	ISO Class	Material	Key Attribute	VWR Item
Ultimate 60™	3-4	Polyester	Sealed-edge, Standard weight	10127-658
Ultimate 70™	3-4	Polyester	Sealed-edge, Heavy weight	10127-662
Ultimate 90™	3-4	Polyester	Sealed-edge, 2-Ply Ultrasonically Quilted	10127-664
Nova-Tech™	5	Poly/Cellulose blend	Lint-free (our most popular class 1000 wipe!)	43500-100
Nova-Tech 56™	5	Poly/Cellulose blend	Light Weight, Cost-Effective Wipe	10135-760
Vision 10™	5	Polyester	Standard weight	10141-364
Vision 20™	5	Polyester	Heavy weight	10141-366
Vision 50™	5	Polyester	ESD / Anti-static	75801-480
Nova-Cot™	6	Cotton	100% Woven Cotton Twill	24000-012
Synergy	6	Polypropylene & Wood pulp	Non-abrasive and solvent resistant	24000-004
Select-Sat™	6	Poly/Cellulose blend	Packaged dry in containers, just add solvent of your choice	10141-996

## Sterile Wipes

Wipe name	ISO Class	Material	Key Attribute	VWR Item
Cobalt Blue Dry Wipes	5	Poly/Cellulose blend	Validated Sterile, Cobalt-60 Gamma Irradiated	10210-956

## Lab Wipes

Wipe name	Material	Key Attribute	VWR Item
Friendly Green "The Little One"	100% recycled cellulose	Eco-friendly	89134-772
Friendly Green "The Big One"	100% recycled cellulose	Eco-friendly	89134-774

## Industrial Wipes

Wipe name	Material	Key Attribute	VWR Item
Blue Guard - 12"x15"	Poly-cellulose	Absorbent, low-linting	76169-636

## Cleanroom Wipe Dispenser

Product Name	Description	VWR Item
GRAB-EEZ Dispenser GE-DISP1	ISO Class 5/Class 100 Cleanroom Compatible Wipe Dispenser, ESD Safe, Eliminates wasteful wipe usage and cross-contamination	89492-640
Grab-EEZ Wipe Refills Poly-cellulose cleanroom wipes, 7"x 12" GE-NT1-712	Non-woven poly-cellulose cleanroom wipes, 7" x 12". 200 wipes/bag, 12 bags/case (2,400 wipes/case)	10146-288

