

Protide Pharmaceuticals, Inc. is devoted to the discovery and development of technologies and products in cell therapy, regenerative medicine, stem cell biology, and cell engineering. Protide has over 35 years of experience in R&D and manufacturing.

TCM™

General serum replacement for many cell types from a variety of species. TCM™ is a defined, multi-purpose serum replacement useful in studying the effects of growth factors and the production of proteins. TCM™ contains highly purified, heat-treated bovine serum albumin, heat-treated bovine transferrin and recombinant insulin. Does not contain growth factors, steroid hormones, glucocorticoids, cell adhesion factors, detectable Ig and mitogens. cGMP



Manufactured. FOR CELL CULTURE USE. FOR FURTHER MANUFACTURING USE. NOT FOR HUMAN USE.

Product Highlights	
<ul style="list-style-type: none"> Concentration: 50X 	<ul style="list-style-type: none"> Total Protein: 664ug/ml at 1X Concentrate
<ul style="list-style-type: none"> Endotoxin: <0.5EU/ml at 1X 	<ul style="list-style-type: none"> Storage: 2-8°C
<ul style="list-style-type: none"> Protein Origin: Bovine* 	<ul style="list-style-type: none"> Stability: Minimum 6 months at 2-8°C

*Bovine serum albumin and transferrin are highly purified and heat treated. These components are sourced from documented cattle born, raised and bled in the USA and inspected by the USDA.

VWR Cat. No.	Supplier Cat. No.	Product Volume	Final Packaging Configuration
TCM™			
76443-252	1001	40ml	Media Bottle
76443-254	1002	100ml	Media Bottle
76443-256	1005	500ml	Media Bottle