



Item Number SU104
Item Sucrose Octaacetate, NF
CAS Number 126-14-7
Molecular Formula $C_{28}H_{38}O_{19}$
Molecular Weight 678.60

Assay (Anhydrous Basis).....98.0-100.5%
Melting Temperature..... Min. 78°C
Acidity..... To pass test

MAXIMUM LIMITS

Water..... 1.0%
Residue on Ignition..... 0.1%
Residual solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843