

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Dark Amber Polypropylene Carboys with Versatile Caps
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2022-05-02

VWR Catalogue Number	Capacity	Cap Size
89233-024	5 L (1.3 gal.)	80 mm
89233-032	20 L (5.3 gal.)	120 mm

Quality System Compliance

Products are manufactured under the [ISO 13485](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	UV Block Dark Amber Polypropylene (PP); USP Class VI O-Ring Material Specification Platinum cured silicone; White
Sterilization	Product is non-sterile.
ATP Assay	Product is not tested for ATP Assay.
DNase & RNase Free	Product is DNase & RNase Free.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

Cytotoxicity	All component materials have been tested and met the requirements for United States Pharmacopoeia (USP) class VI Biological Test for Plastics, current edition. The plastics meet the requirements of the United States Food and Drug Administration (FDA) for food and beverage contact in 21 CFR 177.
Non-Pyrogenic Statement	This product is free from pyrogens.
Latex Statement	The product is Latex free.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP), DEHP (di(2- ethylhexyl) phthalate)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store at under normal warehouse storage conditions.
Shelf Life	10 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.