



CORIO PRODUCT LINE

Working Temperature Range: -40°C...+150°C

PRECISE. CONSISTENT. CONTROL.

When it comes to temperature control companies, JULABO USA offers unparalleled quality and service. JULABO manufactures and provides the most innovative liquid temperature control solutions. Science, research, and industry depend on JULABO to deliver highly precise temperature control for internal and external applications and rapid changes in temperature.

CORIO refrigerated and heated circulators offer a great price-to-performance ratio. They are ideally suited for standard tasks and routine work in laboratories and industrial facilities. They come equipped with the core characteristics you need to succeed with your daily applications. Models are available for internal and external applications and cover a working temperature range of -40°C to +150°C.

VWR Cat. No.	JULABO Model	Voltage	Temp Range	Heating Capacity	Cooling Capacity (20 °C)	Cooling Capacity (0 °C)	Cooling Capacity (-20 °C)	Bath Volume	Opening (WxL/D)	Price USD
10015-724	CD	115V/60Hz	+20...+150 °C	1 kW	N/A	N/A	N/A	Up to 50 L	N/A	\$1,629.01
10015-732	CD-B5	115V/60Hz	+20...+150 °C	1 kW	N/A	N/A	N/A	3.5 - 5 L	5.9x5.9/5.9"	\$2,278.16
10015-734	CD-B13	115V/60Hz	+20...+150 °C	1 kW	N/A	N/A	N/A	9 - 13 L	11.8x5.9/7.1"	\$2,620.94
10015-736	CD-B17	115V/60Hz	+20...+150 °C	1 kW	N/A	N/A	N/A	13 - 17 L	11.8x7/7.8"	\$2,762.60
10015-738	CD-B19	115V/60Hz	+20...+150 °C	1 kW	N/A	N/A	N/A	14 - 19 L	11.8x13.8/5.9"	\$2,742.21
10015-740	CD-B27	115V/60Hz	+20...+150 °C	1 kW	N/A	N/A	N/A	17 - 27 L	11.8x13.8/7.9"	\$2,996.54
10015-742	CD-B33	115V/60Hz	+20...+150 °C	1 kW	N/A	N/A	N/A	26 - 39 L	25.9x12.5/5.9"	\$4,797.02
10015-744	CD-B39	115V/60Hz	+20...+150 °C	1 kW	N/A	N/A	N/A	35 - 41 L	13x11.8/11.8"	\$4,943.22
10015-746	CD-BC4	115V/60Hz	+20...+150 °C	1 kW	N/A	N/A	N/A	3 - 4.5 L	5.1x5.9/5.9"	\$3,332.54
10015-748	CD-BC6	115V/60Hz	+20...+150 °C	1 kW	N/A	N/A	N/A	4.5 - 6 L	5.1x5.9/7.9"	\$3,846.96
10015-814	CD-BC12	115V/60Hz	+20...+150 °C	1 kW	N/A	N/A	N/A	8.5 - 12 L	8.7x5.9/7.9"	\$4,274.71
10015-816	CD-BC26	115V/60Hz	+20...+150 °C	1 kW	N/A	N/A	N/A	19 - 26 L	10.2x13.8/7.9"	\$4,916.46
10015-818	CD-200F	115V/60Hz	-20...+150 °C	1 kW	220W	170W	60W	3 - 4 L	5.1x5.9/5.9"	\$5,688.34
10015-750	CD-201F	115V/60Hz	-20...+150 °C	1 kW	220W	160W	60W	3 - 4 L	5.1x5.9/5.9"	\$5,688.34
10015-752	CD-300F	115V/60Hz	-25...+150 °C	1 kW	310W	280W	110W	3 - 4 L	5.1x5.9/5.9"	\$6,578.04
10015-754	CD-600F	115V/60Hz	-35...+150 °C	1 kW	600W	530W	220W	5 - 7.5 L	8.7x5.9/5.9"	\$8,309.97
75839-420	CD-601F	115V/60Hz	-40...+150 °C	1 kW	600W	500W	200W	8 - 10 L	8.7x5.9/7.9"	\$8,424.85
75839-422	CD-601F	230V/60Hz	-40...+150 °C	2 kW	600W	500W	200W	8 - 10 L	8.7x5.9/7.9"	\$8,424.85
75839-424	CD-900F	115V/60Hz	-40...+150 °C	1 kW	900W	500W	200W	21 - 30 L	10.2x13.8/7.9"	\$8,262.45
75839-566	CD-900F	230V/60Hz	-40...+150 °C	2 kW	900W	500W	200W	21 - 30 L	10.2x13.8/7.9"	\$8,262.45
75839-446	CD-1000F	115V/60Hz	-40...+150 °C	1 kW	1000W	900W	500W	5 - 7.5 L	8.7x5.9/7.9"	\$10,415.72
75871-448	CD-1000F	230V/60Hz	-40...+150 °C	2 kW	1000W	900W	500W	5 - 7.5 L	8.7x5.9/7.9"	\$10,415.72
75871-450	CD-1001F	230V/60Hz	-38...+100	2 kW	1000W	900W	350W	42 - 56 L	13.8x16.1/11.8"	\$12,802.99

Units include M16x1 male pump connections with 8 & 12 mm hose barbs // adjustable internal/external flow // USB connection for PC communication or data logging // Pump flow rate: 15 L/min // Max pressure: 5.1 psi // Stability ±0.03 °C