

# Certificate of Quality



|                          |                                    |
|--------------------------|------------------------------------|
| <b>Part Number</b>       | <b>MFLX30613-55</b>                |
| <b>Lot Number</b>        | <b>1377555</b>                     |
| <b>Part Description</b>  | <b>TEE L/S 17/24 + I/P26 10/PK</b> |
| <b>Country of Origin</b> | <b>USA</b>                         |

## Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Product Specifications & Compliance

- Compound: Natural Polypropylene (Raw material)
- The supplier has indicated that this material Natural Polypropylene meets the FDA requirements outlined in the 21 CFR 177.1520(a)(3)(i) and (c)3.1a.
- The supplier also indicates that according to information all other ingredients used in the formulation meet their respective FDA regulations and 21 CFR 177.1520(b).  
Specifically, this product meets the FDA requirements for food contact, except for cooking, as listed under conditions of use C through H in 21 CFR 176.170(c), Table 2.
- This product has been certified to USP Class VI Plastic (USP <88>, biological Reactivity Tests, In Vivo) and USP <661> requirements.
- Material is gamma irradiation stable.
- Product is RoHS compliant EU Directives 2011/65/EU R.H.S. (RoHS 2) & 2015/863/EU R.H.S. (RoHS 3).
- Tallow/BSE/TSE - supplier does not use tallow or any other animal derived materials in the manufacture or formulation of this product. Product is Animal Derivative Free.
- The manufacturer has stated this product has a 10 year shelf life if stored at ambient temperature and out of direct sunlight.
- Note - End user must make their own determinations ensuring the use of the products are safe, lawful, and suitable for their intended applications.

## Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

## Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.