



Item Number BE129
Item Benzethonium Chloride, USP
CAS Number 121-54-0
Molecular Formula $C_{27}H_{42}ClNO_2$
Molecular Weight 448.08

Assay (Dried Basis)..... 97.0-103.0%
Melting Range..... 158°-163°C

MAXIMUM LIMITS

Loss on Drying..... 5.0%
Residue on Ignition..... 0.1%
Ammonium Compounds..... To pass test
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843