

# Cole-Parmer®



Promo code  
6088

## CG-500, CG-450, CG-400, Freezer/Mill® Cryogenic Grinders



**Receive a complete starter setup package free with every CG-500, CG-450, or CG-400 purchased.**

Qualifying Equipment		
Item #	VWR #	Equipment Description
6875D-115	76182-012	Freezer/Mill® Dual Chamber Cryogenic Grinder, Model CG-500
6875A-115	76182-008	Freezer/Mill® High Capacity Cryogenic Grinder, Model CG-450
6875-115	76182-004	Freezer/Mill® High Capacity Cryogenic Grinder, Model CG-400
Complete Starter Package		
Item #	VWR #	Accessory Description
6801	76182-236	Large Grinding Vial Set - Contains everything necessary to make one complete grinding vial plus three spare Polycarbonate Center Cylinders. Includes one Large Stainless-Steel Impactor, two Large Stainless-Steel End Plugs, and one pack of Large Polycarbonate Center Cylinders. Sample capacity 10 to 100 g.
6870L	76182-442	Large Grinding Vial Accessory Package - Package includes Large Extractor/Vial Opener, Vial Rack for Large Grinding Vials
Choice of Cryogenic Gloves (Select ONE)		
Item #	VWR #	Size
6900S	76182-040	Small
6900M	76182-038	Medium
6900L	76182-036	Large
6900XL	76182-042	Extra Large

## Typical Applications Include

DNA/RNA extraction, forensic/pharmaceutical testing, RoHS/WEEE, consumer product testing, food safety, materials research, and medical research.

## Sample Types Include



Plastics & Polymers



Animal Tissue



Pharmaceuticals



Cannabis




Yeast

## Specifications – CG-500 Freezer/Mill

High-capacity cryogenic grinder with dual grinding and cooling chambers. It can grind and cool up to 200 g of samples (0.5 to 100 g per chamber). Each grinding chamber holds one large grinding vial, one mid-size vial, four small grinding vials, or four microvial sets each. Comes equipped with liquid nitrogen auto-fill.


Voltage	115 VAC, 60 Hz or 230 VAC, 50 Hz
Dimensions	50.8 x 54.6 x 64.8 cm (20.0 x 21.5 x 25.5")
Net Weight	35.4 kg (78 lb)
Warranty	3 years



## Specifications – CG-450 and CG-400 Freezer/Mill

Large cryogenic grinder that accommodates up to 100 g of samples. It has one dual-purpose chamber that pre-cools and grinds samples. Its grinding chamber holds one large grinding vial, one mid-size vial, four small grinding vials, or four microvial sets. The CG-450 comes equipped with liquid nitrogen auto-fill. The CG-400 is manually filled with liquid nitrogen.

Voltage	115 VAC, 60 Hz or 230 VAC, 50 Hz
Dimensions	52.0 x 55.0 x 46.0 cm (20.5 x 21.5 x 18.0")
Net Weight	20 kg (44 lb)
Warranty	3 years



Purchase by calling 1.800.932.5000, visiting [vwr.com](http://vwr.com) or contacting your Avantor sales representative. Have your VWR purchase order, VWR invoice or VWR packing slip ready and fill out the online form at [vwr.com/promotions](http://vwr.com/promotions) and search for **Promo Code 6088** before October 31, 2024. Offer valid on orders placed between January 1, 2024 and September 30, 2024. Redeem before October 31, 2024. While supplies last. Allow up to 8 weeks to receive your free product. Void where prohibited.