



Material	00884-700
Material description	Acetone
Grade	ACS
CAS Number	67-64-1
Molecular formula	C <sub>3</sub> H <sub>6</sub> O
Molecular mass	58,08
Storage	Store between 15° C and 30° C.

## Characteristics

## Specifications

Assay (GC, corrected for water)	Min. 99,5 %
Solubility in water	Passes test
Substances reducing permanganate	Passes test
Color (Pt-Co)	Max. 10 APHA
Titration acid	Max. 0,0003 meq/g
Titration base	Max. 0,0006 meq/g
Aldehyde (as HCHO)	Max. 0,002 %
Isopropyl alcohol	Max. 0,05 %
Methanol	Max. 0,05 %
Residue after evaporation	Max. 0,001 %
Water	Max. 0,5 %

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine  
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,  
QC, Canada