

Certificate of Quality



| | |
|-------------------|--------------------------------------|
| Part Number | MFLX31809-42 |
| Lot Number | 28565502 |
| Part Description | ELBOW 90-DEG 1IN LADISH CLAMP |
| Expiration Date | 08OCT2061 |
| Country of Origin | USA |

Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Product Specifications & Compliance

- Compound: PVDF
- Material grade: Sani-Pro® K PVDF.
- Material meets 21 CFR 177.2510 requirements.
- Material meets 3-A requirements for Sanitary Standards for Multiple-Use Plastic Materials Used as Product Contact Surfaces for Dairy Equipment.
- Material meets USDA requirements or use in processing or storage areas in contact with meat or poultry food products.
- Material meets USP Class VI requirements.
- No animal derived ingredients are used in the manufacturing of Sani-Pro® K PVDF material.
- Shelf life: 40 years.
- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.