

Think. Inside. The. Box.

PURELAB® Classic UF*



DI pack

LC186

Lifetime expectancy: 12 months

What is purified from the water?

Ions

Impact of not changing:

Weakly-bound ions released back into the water



Ultra Filter

LC169

Lifetime expectancy: 24 months

What is purified from the water?

Bacteria and microbial by-products such as endotoxins

Impact of not changing:

Bacterial and endotoxin specification will not be maintained

Reduced flow rate



Point of Use Filter**

LC134

Lifetime expectancy: 6 months

What is purified from the water?

Particulates and bacteria

Impact of not changing:

Potential recontamination of purified water and reduced flow rates

*A combination of treatments and technologies produce the final water specification

** For DNase, RNase & endotoxin specification at point of use, please use a LC197 Biofilter