


VWR® REAGENT RACK SYSTEMS

Flexible island shelving

Excellent for new or
existing workspaces

Variety of sizes and
accessories available



VWR® REAGENT RACK SYSTEMS

VWR® Reagent Rack systems provide excellent storage of laboratory equipment above the lab station. These systems are an ideal addition to new or existing laboratory workspaces and casework, where additional shelving is needed. VWR's Reagent Rack uprights are available in three heights and can support a range of shelf and power options.

Contact your local VWR Customer Service Center today, for more information!

Reagent Rack Components

Slotted Upright Posts

- Includes bottom mounting plate and hardware
- Mounting plates are 2.76" x 3.94". See drawing on page 4.

Upright Height	Mounting	VWR Cat No
31.10"	Left or Right	10808-826
31.10"	Center	10808-828
39.37"	Left or Right	10808-830
39.37"	Center	10808-832
47.24"	Left or Right	10808-834
47.24"	Center	10808-836

Component	VWR Cat No
30"W Horizontal Tube	10486-500
36"W Horizontal Tube	10486-502
48"W Horizontal Tube	10486-504
60"W Horizontal Tube	10486-506
72"W Horizontal Tube	10486-508
Corner bracket	10486-510
Straight/"T" bracket	10486-512



Preconfigured Reagent Rack

Top Mounted Reagent Rack

- Other heights and widths available upon request.

	30"W	36"W	48"W	60"W	72"W
Shelf Bay 31.10"	10808-838	10808-850	10808-862	10808-874	10808-886
Shelf Extension 31.10"	10808-840	10808-852	10808-864	10808-876	10808-888
Shelf Bay 39.37"	10808-842	10808-854	10808-866	10808-878	10808-890
Shelf Extension 39.37"	10808-844	10808-856	10808-868	10808-880	10808-892
Shelf Bay 47.24"	10808-846	10808-858	10808-870	10808-882	10808-894
Shelf Extension 47.24"	10808-848	10808-860	10808-872	10808-884	10808-896

Shelf Bays Include: (2) Uprights w/ mounting hardware, (1) Top horizontal rail, (2) Corner brackets

Shelf Bay Extensions Include: (1) Upright w/ mounting hardware, (1) Top horizontal rail, (1) T bracket

Shelves for Reagent Rack Systems

Steel Shelf with Standard Brackets

- 176 max lbs per shelf.

Shelf Depth	Shelf Width	VWR Cat No
11.81"	30"	89524-486
11.81"	36"	89524-488
11.81"	48"	89524-490
11.81"	60"	89524-492
11.81"	72"	89524-494
15.74"	30"	89524-496
15.74"	36"	89524-498
15.74"	48"	89524-500
15.74"	60"	89524-502
15.74"	72"	89524-504




Steel Shelf w/ standard bracket (included)

Bookend Brackets

- Includes pair of 2 brackets

Description	VWR Cat No
11.81" Bookend Bracket Pair	89524-486
15.74" Bookend Bracket Pair	89524-488



Shelf Lips

- Lip height: 1.5" (.32" above shelf)

Description	VWR Cat No
30"W Shelf lip	89524-510
36"W Shelf lip	89524-512
48"W Shelf lip	89524-514
60"W Shelf lip	89524-516
72"W Shelf lip	89524-518



Accessories for Reagent Rack Systems

LED Lighting

- Includes 13' power cord
- Daisy chain up to 10 fixtures with one power source

Description	Size	Power	VWR Cat No
24" dual intensity LED light fixture, magnetic attachment	22.44" x 1.26" x 1.81"	1200Lm 5000K	89524-520
36" dual intensity LED light fixture, magnetic attachment	33.46" x 1.26" x 1.81"	1800Lm 5000K	89524-522
48" dual intensity LED light fixture, magnetic attachment	44.48" x 1.26" x 1.81"	2400Lm 5000K	89524-524
Extension cable (for daisy-chaining lights)*	12"L	-----	89524-526
Extension cable (for daisy-chaining lights)*	36"L	-----	89524-528



High/Low/Off switch allows for varying light intensity



Magnetic mounting attachment under metal shelves.

Power Rails

- Includes 10' power cord.
- Vertical power rails can also be mounted underneath the worksurface.

Description	VWR Cat No
30" Bin rail w/6 outlet power	10611-600
30" Bin rail w/6 outlet power, surge	10611-602
36" Bin rail w/6 outlet power	10611-604
36" Bin rail w/6 outlet power, surge	16011-606
48" Bin rail w/6 outlet power	89524-530
48" Bin rail w/6 outlet power, surge	89524-532
60" Bin rail w/6 outlet power	89524-534
60" Bin rail w/6 outlet power, surge	89524-536
72" Bin rail w/6 outlet power	89524-538
72" Bin rail w/6 outlet power, surge	89524-540
6 outlet vertical power bar with clamp	89524-542
6 outlet vertical power bar with clamp, surge	89524-544
12 outlet vertical power bar with clamp	89524-546
12 outlet vertical power bar with clamp, surge	89524-548



Vertical power rail



Horizontal Power Rail

Additional accessories available, visit vwr.com for more information.

Dimensions and Mounting Hole Locations

