

Material No.: 5635-02
Revision No.: 0

Product Specification

Test	Specification
Appearance (Clear, colorless liquid)	Passes Test
Normality (eq/L)	0.995 – 1.005
Chloride (Cl)	≤ 5 ppm
ACS – Heavy Metals (as Pb)	≤ 1 ppm
Trace Impurities – Iron (Fe)	≤ 0.5 ppm

For Laboratory, Research, or Manufacturing Use
Standardization at 25°C traceable to NIST Standard Reference Material. Filtered through a 0.2 micron filter.
Storage Condition: Protect from air to avoid absorption of carbon dioxide.

Packaging Site: Paris Mfg Ctr & DC



Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700

Page 1 of 1