



**Item Number** TR149  
**Item** Tromethamine, USP  
**CAS Number** 77-86-1  
**Molecular Formula**  $C_4H_{11}NO_3$   
**Molecular Weight** 121.14

Assay (Dried Basis)..... 99.0-101.0%  
Melting Range..... 168°-172°C  
pH of a 5% Solution @ 25°C..... 10.0-11.5

MAXIMUM LIMITS

Loss on Drying..... 1.0%  
Residue on Ignition..... 0.1%  
Heavy Metals.....0.001%  
Residual solvents..... To pass test

***Spectrum Chemicals & Laboratory Products***

---

**Corporate Headquarters:**      **East Coast Facility:**  
14422 S. San Pedro St.      755 Jersey Ave.  
Gardena, CA 90248      New Brunswick, NJ 08901  
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843