

**avantor™**  
delivered by **VWR™**

VWR.COM

**ELGA**  **VEOLIA**

# PURELAB® Dispenser

## Uninterrupted Discovery



**WATER TECHNOLOGIES**

# Proven reliability

Unrivalled purity in design, unmatched flexibility & effortless use. Together with PURELAB® Chorus & Quest, the PURELAB® Dispenser enables uninterrupted discovery.



## Accuracy:

Flexible dispensing with unrivalled built in purity monitoring sensor

- Delivering proof of water purity
- Repeated volume dispensing
- Improving lab efficiency
- Reproducible results

## Streamlined footprint:

- Maximized lab space
- Optimization of storage space
- Up to 4 dispensers per water purifier\*

\*Dependent on water purifier purchased.



# Ultimate flexibility



## Simple use and maintenance:

- Plug and Play installation with PURELAB® range
- Height adjustable, 180° rotating arm and flexible dispenser arm
- IPX7 rating dispenser handset with menu navigation
- Optional connection to point of use filtration

## Efficient

- From precise drop by drop up to 2 litres per minute allowing for the right speed including calibration work
- Variable flow rate for easier filling of different sizes of labware
- Volumetric dispensing freeing you to continue working while water is being dispensed
- Profile & locked dispensing allowing you to choose your desired amount of water and prevent accidental flooding

# Product Specifications

**Height**  
580mm

**Depth**  
375mm

**Width**  
264mm


Intuitive dispenser handset

Multiple dispensing options

Point of use filter

Height adjustable button

Flexible dispensing arm



**Dispatch weight:** 7.7kg

**Operational Weight:** 4.5kg

**Installation:** Free-standing on a worktop or secured to a worktop using the mounting bracket provided.  
(See quick start guide for details)

**Mains input**

100-240 VAC, 50-60 Hz

**Power required**

24 V DC

**Power consumption**

120 VA

**Noise Output**

dBA - < 40

**Inlet**

8mm OD Tube

**Return**

8mm OD Tube

**Flow bends**

8mm

(provided in the basic installation kit with the main product)

## Dedicated to Discovery

ELGA is the worldwide specialist in the innovation, engineering, service & support of water purification systems.

Unrivalled product designs have achieved international recognition and awards.

Worldwide technical service teams support science & healthcare globally. Water treatment experts provide in-depth understanding of analytical techniques.

Hubgrade, ELGA's digital platform, monitors equipment performance to ensure laboratory work continues uninterrupted.

A global supply chain supports clients from regional centres on every continent.

ELGA is part of Veolia, the world's largest environmental management organisation.

Supplier part number: LA824  
VWR Catalog number: 76514-420

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.