

# SAFETY DATA SHEET

## 1. Identification

**Product identifier:** Tetrahydrofuran

### Other means of identification

**Product No.:** 2858, 8498, 9432, 9439, 9440, 9441, 9446, 9447, 9450, V530, V558, 31800

### Recommended use and restriction on use

**Recommended use:** For Laboratory, Research or Manufacturing Use.

**Restrictions on use:** Not determined.

### Details of the supplier of the safety data sheet

	Avantor Performance Materials, LLC 100 Matsonford Rd, Suite 200 Radnor, PA 19087
Telephone:	
	Customer Service: 855-282-6867
Fax:	
Contact Person:	Product Information Compliance
E-mail:	info@avantormaterials.com

### Emergency telephone number:

CHEMTREC: 1-800-424-9300 within US and Canada (24 hrs/day, 7 days/week)

## 2. Hazard identification

### Hazard Classification

#### Physical Hazards

Flammable liquids	Category 2
Static-accumulating flammable liquid	Category 1

#### Health Hazards

Acute toxicity (Oral)	Category 4
Skin Corrosion/Irritation	Category 2
Serious Eye Damage/Eye Irritation	Category 2A
Carcinogenicity	Category 2
Specific Target Organ Toxicity - Single Exposure	Category 3 <sup>1</sup>

#### Target Organs

1. Respiratory tract irritation.

### Label Elements

**Hazard Symbol:**



<b>Signal Word:</b>	Danger
<b>Hazard Statement:</b>	Highly flammable liquid and vapor. Harmful if swallowed. Causes skin irritation. Causes serious eye irritation. Suspected of causing cancer. May cause respiratory irritation. May cause drowsiness or dizziness.
<b>Precautionary Statements</b>	
<b>Prevention:</b>	Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Use personal protective equipment as required. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. Ground and bond container and receiving equipment. Use explosion-proof [electrical/ventilating/lighting] equipment. Use non-sparking tools. Take action to prevent static discharges. Wear protective gloves/protective clothing/eye protection/face protection. Avoid breathing dust/mist/vapors/spray. Use only outdoors or in a well-ventilated area. Do not eat, drink or smoke when using this product. Wash thoroughly after handling.
<b>Response:</b>	In case of fire: Use water spray, foam, dry powder or carbon dioxide for extinction. IF exposed or concerned: Get medical advice/attention. IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell. Rinse mouth. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/doctor if you feel unwell. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower]. If skin irritation occurs: Get medical advice/attention. Take off contaminated clothing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. Collect spillage.
<b>Storage:</b>	Keep container tightly closed. Store in a well-ventilated place. Keep cool. Store locked up.
<b>Disposal:</b>	Dispose of contents/container to an appropriate treatment and disposal facility in accordance with applicable laws and regulations, and product characteristics at time of disposal.
<b>Other hazards which do not result in GHS classification:</b>	Static accumulating flammable liquid can become electrostatically charged even in bonded and grounded equipment.

<b>3. Composition/information on ingredients</b>
--

**Substances**

Chemical name	Common name and synonyms	CAS number	Content in percent (%)*
Tetrahydrofuran		109-99-9	99 - 100%

\* All concentrations are percent by weight unless ingredient is a gas. Gas concentrations are in percent by volume.

**4. First-aid measures**

<b>General information:</b>	Get medical advice/attention if you feel unwell. Show this safety data sheet to the doctor in attendance.
<b>Ingestion:</b>	Call a physician or poison control center immediately. Do not induce vomiting without advice from poison control center. Never give liquid to an unconscious person.
<b>Inhalation:</b>	Move to fresh air. Get medical attention if symptoms persist. If breathing is difficult, give oxygen. If breathing stops, provide artificial respiration.
<b>Skin Contact:</b>	Immediately flush with plenty of water for at least 15 minutes while removing contaminated clothing and shoes. Get medical attention. Wash contaminated clothing before reuse. Destroy or thoroughly clean contaminated shoes.
<b>Eye contact:</b>	Immediately flush with plenty of water for at least 15 minutes. If easy to do, remove contact lenses. Get medical attention.

**Most important symptoms/effects, acute and delayed**

<b>Symptoms:</b>	Irritating to eyes, respiratory system and skin.
<b>Hazards:</b>	Suspected of causing cancer.

**Indication of immediate medical attention and special treatment needed**

<b>Treatment:</b>	Treat symptomatically. Symptoms may be delayed.
-------------------	---

**5. Fire-fighting measures**

<b>General Fire Hazards:</b>	Vapors may cause a flash fire or ignite explosively. Vapors may travel considerable distance to a source of ignition and flash back.
------------------------------	--

**Suitable (and unsuitable) extinguishing media**

**Suitable extinguishing media:** Water spray, foam, dry powder or carbon dioxide.

**Unsuitable extinguishing media:** Avoid water in straight hose stream; will scatter and spread fire.

<b>Specific hazards arising from the chemical:</b>	Vapors may cause a flash fire or ignite explosively. Vapors may travel considerable distance to a source of ignition and flash back. Prevent buildup of vapors or gases to explosive concentrations. May form explosive peroxides.
--	--

**Special protective equipment and precautions for firefighters**

- Special fire fighting procedures:** Use water spray to keep fire-exposed containers cool. Water may be ineffective in fighting the fire. Fight fire from a protected location. Move containers from fire area if you can do so without risk. Cool containers exposed to flames with water until well after the fire is out. Vapors may travel considerable distance to a source of ignition and flash back. Vapors may cause a flash fire or ignite explosively.
- Special protective equipment for fire-fighters:** Firefighters must use standard protective equipment including flame retardant coat, helmet with face shield, gloves, rubber boots, and in enclosed spaces, SCBA.

**6. Accidental release measures**

- Personal precautions, protective equipment and emergency procedures:** Use personal protective equipment. Keep unauthorized personnel away. Keep upwind. Ventilate closed spaces before entering them. ELIMINATE all ignition sources (no smoking, flares, sparks or flames in immediate area). See Section 8 of the SDS for Personal Protective Equipment. Do not touch damaged containers or spilled material unless wearing appropriate protective clothing.
- Methods and material for containment and cleaning up:** In case of leakage, eliminate all ignition sources. Absorb spill with vermiculite or other inert material, then place in a container for chemical waste. Clean surface thoroughly to remove residual contamination. Dike far ahead of larger spill for later recovery and disposal.
- Notification Procedures:** Dike for later disposal. Prevent entry into waterways, sewer, basements or confined areas. Stop the flow of material, if this is without risk. Inform authorities if large amounts are involved.
- Environmental Precautions:** Do not contaminate water sources or sewer. Prevent further leakage or spillage if safe to do so. Avoid discharge into drains, water courses or onto the ground.

**7. Handling and storage**

- Precautions for safe handling:** DO NOT handle, store or open near an open flame, sources of heat or sources of ignition. Protect material from direct sunlight. Take precautionary measures against static discharges. Contact with air and light may form explosive peroxides. If peroxide formation is suspected, do not open or move container. Use personal protective equipment as required. Avoid breathing mists or vapors. Do not taste or swallow. Use only with adequate ventilation. Avoid contact with eyes, skin, and clothing. Wash hands thoroughly after handling. See Section 8 of the SDS for Personal Protective Equipment.
- Conditions for safe storage, including any incompatibilities:** Keep away from food, drink and animal feeding stuffs. Prolonged contact with air may cause formation of explosive peroxides. Nitrogen blanketing of containers is recommended. Keep container tightly closed in a cool, well-ventilated place. Ground container and transfer equipment to eliminate static electric sparks. Comply with all national, state, and local codes pertaining to the storage, handling, dispensing, and disposal of flammable liquids.

**8. Exposure controls/personal protection**

**Control Parameters**

**Occupational Exposure Limits**

Chemical Identity	Type	Exposure Limit Values	Source
-------------------	------	-----------------------	--------

Tetrahydrofuran	STEL	100 ppm 295 mg/m3	Canada. Alberta OELs (Occupational Health & Safety Code, Schedule 1, Table 2) (07 2009)
	TWA	50 ppm 147 mg/m3	Canada. Alberta OELs (Occupational Health & Safety Code, Schedule 1, Table 2) (07 2009)
Tetrahydrofuran	TWA	50 ppm	Canada. British Columbia OELs. (Occupational Exposure Limits for Chemical Substances, Occupational Health and Safety Regulation 296/97, as amended) (07 2007)
	STEL	100 ppm	Canada. British Columbia OELs. (Occupational Exposure Limits for Chemical Substances, Occupational Health and Safety Regulation 296/97, as amended) (07 2007)
Tetrahydrofuran	TWA	50 ppm	Canada. Manitoba OELs (Reg. 217/2006, The Workplace Safety And Health Act) (03 2011)
	STEL	100 ppm	Canada. Manitoba OELs (Reg. 217/2006, The Workplace Safety And Health Act) (03 2011)
Tetrahydrofuran	STEL	100 ppm	Canada. Ontario OELs. (Control of Exposure to Biological or Chemical Agents) (11 2010)
	TWA	50 ppm	Canada. Ontario OELs. (Control of Exposure to Biological or Chemical Agents) (11 2010)
Tetrahydrofuran	15 MIN ACL	100 ppm	Canada. Saskatchewan OELs (Occupational Health and Safety Regulations, 1996, Table 21) (05 2009)
	8 HR ACL	50 ppm	Canada. Saskatchewan OELs (Occupational Health and Safety Regulations, 1996, Table 21) (05 2009)
Tetrahydrofuran	TWA	100 ppm 300 mg/m3	Canada. Quebec OELs. (Ministry of Labor - Regulation Respecting the Quality of the Work Environment) (09 2017)
Tetrahydrofuran	TWA	50 ppm	US. ACGIH Threshold Limit Values (2011)
	STEL	100 ppm	US. ACGIH Threshold Limit Values (2011)

**Biological Limit Values**

Chemical Identity	Exposure Limit Values	Source
Tetrahydrofuran (tetrahydrofuran: Sampling time: End of shift.)	2 mg/l (Urine)	ACGIH BEI (03 2013)

**Appropriate Engineering Controls** No data available.

**Individual protection measures, such as personal protective equipment**

**General information:** Good general ventilation (typically 10 air changes per hour) should be used. Ventilation rates should be matched to conditions. If applicable, use process enclosures, local exhaust ventilation, or other engineering controls to maintain airborne levels below recommended exposure limits. If exposure limits have not been established, maintain airborne levels to an acceptable level. Use explosion-proof ventilation equipment.

**Eye/face protection:** Wear safety glasses with side shields (or goggles) and a face shield.

**Skin Protection**

**Hand Protection:** Chemical resistant gloves

**Other:** Wear suitable protective clothing.

**Respiratory Protection:** In case of inadequate ventilation use suitable respirator. Chemical respirator with organic vapor cartridge and full facepiece.

**Hygiene measures:** Provide eyewash station and safety shower. Observe good industrial hygiene practices. Do not eat, drink or smoke when using the product. Wash hands before breaks and immediately after handling the product. Wash contaminated clothing before reuse. Avoid contact with eyes, skin, and clothing.

## 9. Physical and chemical properties

### Appearance

Physical state:	Liquid
Form:	Liquid
Color:	Colorless
Odor:	Ether-like
Odor threshold:	No data available.
pH:	No data available.
Melting point/freezing point:	-108,3 °C
Initial boiling point and boiling range:	65 °C
Flash Point:	-14 °C (Closed Cup)
Evaporation rate:	8 (butyl acetate=1)
Flammability (solid, gas):	Class IB Flammable Liquid
<b>Upper/lower limit on flammability or explosive limits</b>	
Flammability limit - upper (%):	11,8 %(V)
Flammability limit - lower (%):	1,8 %(V)
Explosive limit - upper (%):	No data available.
Explosive limit - lower (%):	No data available.
Vapor pressure:	21,60 kPa (25 °C)
Vapor density:	2,56 (Air=1)
Density:	0,88 g/ml (25 °C)
Relative density:	0,88 (25 °C)
<b>Solubility(ies)</b>	
Solubility in water:	Miscible
Solubility (other):	No data available.
Partition coefficient (n-octanol/water):	0,46
Auto-ignition temperature:	321 °C
Decomposition temperature:	No data available.
Viscosity:	No data available.
<b>Other information</b>	
Minimum ignition energy:	0,54 mJ
Molecular weight:	72,11 g/mol (C <sub>4</sub> H <sub>8</sub> O)

## 10. Stability and reactivity

Reactivity:	No dangerous reaction known under conditions of normal use.
Chemical Stability:	Material is stable under normal conditions.
Possibility of hazardous reactions:	Hazardous polymerization does not occur.
Conditions to avoid:	Heat, sparks, flames. Protect against direct sunlight.
Incompatible Materials:	Strong oxidizing agents. Strong oxidizing agents. Acids. Bases, alkalies (organic). Air. May attack some plastics, rubber and coatings.
Hazardous Decomposition Products:	Thermal decomposition may release oxides of carbon.

**11. Toxicological information**

**Information on likely routes of exposure**

- Inhalation:** May cause respiratory irritation.
- Skin Contact:** Causes skin irritation.
- Eye contact:** Causes serious eye irritation.
- Ingestion:** Harmful if swallowed. May cause irritation of the gastrointestinal tract.

**Information on toxicological effects**

**Acute toxicity (list all possible routes of exposure)**

- Oral**
  - Product:** LD 50 (Rat): 1.650 mg/kg
- Dermal**
  - Product:** LD 50 (Rat): > 2.000 mg/kg
- Inhalation**
  - Product:** LC 50 (Rat): 80975 ppm  
LC 50 (Rat): 18000 - 22000 ppm  
NOAEL (Rat): 15,9 - 16,8 mg/l

**Repeated dose toxicity**

- Product:** None known.

**Skin Corrosion/Irritation**

- Product:** Causes skin irritation.

**Serious Eye Damage/Eye Irritation**

- Product:** Causes serious eye irritation.

**Respiratory or Skin Sensitization**

- Product:** Not a skin nor a respiratory sensitizer.

**Carcinogenicity**

- Product:** Suspected of causing cancer.

**IARC Monographs on the Evaluation of Carcinogenic Risks to Humans:**

No carcinogenic components identified

**US. National Toxicology Program (NTP) Report on Carcinogens:**

No carcinogenic components identified

**ACGIH Carcinogen List:**

No carcinogenic components identified

**Germ Cell Mutagenicity**

- In vitro**
  - Product:** No mutagenic components identified
- In vivo**
  - Product:** No mutagenic components identified

**Reproductive toxicity**

**Product:** No components toxic to reproduction

**Specific Target Organ Toxicity - Single Exposure**

**Product:** Narcotic effect. Respiratory tract irritation.

**Specific Target Organ Toxicity - Repeated Exposure**

**Product:** None known.

**Aspiration Hazard**

**Product:** Not classified

**Other effects:** None known.

**12. Ecological information**

**Ecotoxicity:**

**Acute hazards to the aquatic environment:**

**Fish**

**Product:** No data available.

**Specified substance(s):**

Tetrahydrofuran LC 50 (Fathead minnow (Pimephales promelas), 96 h): 1.970 - 2.360 mg/l  
EC 50 (Fathead minnow (Pimephales promelas), 96 h): 1.930 mg/l

**Aquatic Invertebrates**

**Product:** No data available.

**Specified substance(s):**

Tetrahydrofuran LC 50 (Water flea (Daphnia magna), 24 h): > 10.000 mg/l  
EC 50 (Daphnia magna, 24 h): 5.930 mg/l

**Chronic hazards to the aquatic environment:**

**Fish**

**Product:** No data available.

**Aquatic Invertebrates**

**Product:** No data available.

**Toxicity to Aquatic Plants**

**Product:** No data available.

**Persistence and Degradability**

**Biodegradation**

**Product:** The product is moderately biodegradable.

**BOD/COD Ratio**

**Product:** No data available.

**Bioaccumulative potential**

**Bioconcentration Factor (BCF)**

**Product:** No data available on bioaccumulation.

**Partition Coefficient n-octanol / water (log Kow)**

**Product:** Log Kow: 0,46



**Mobility in soil:** The product is water soluble and may spread in water systems.

**Other adverse effects:** The product components are not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment.

### 13. Disposal considerations

**Disposal instructions:** Discharge, treatment, or disposal may be subject to national, state, or local laws.

**Contaminated Packaging:** Since emptied containers retain product residue, follow label warnings even after container is emptied.

### 14. Transport information

#### TDG

UN Number:	UN 2056
UN Proper Shipping Name:	TETRAHYDROFURAN
Transport Hazard Class(es)	
Class:	3
Label(s):	3
Packing Group:	II
Marine Pollutant:	No
Special precautions for user:	Not determined.

#### IMDG

UN Number:	UN 2056
UN Proper Shipping Name:	TETRAHYDROFURAN
Transport Hazard Class(es)	
Class:	3
Label(s):	3
EmS No.:	F-E, S-D
Packing Group:	II
Marine Pollutant:	No
Special precautions for user:	Not determined.

#### IATA

UN Number:	UN 2056
UN Proper Shipping Name:	Tetrahydrofuran
Transport Hazard Class(es):	
Class:	3
Label(s):	3
Packing Group:	II
Marine Pollutant:	No
Special precautions for user:	Not determined.
Cargo aircraft only:	Allowed.

**Transport in bulk according to Annex II of MARPOL and the IBC Code:** Not applicable

### 15. Regulatory information

#### Canada Federal Regulations

##### List of Toxic Substances (CEPA, Schedule 1)

Not Regulated

##### Export Control List (CEPA 1999, Schedule 3)

Not Regulated

**National Pollutant Release Inventory (NPRI)**

**Canada. National Pollutant Release Inventory (NPRI) Substances, Part 5, VOCs with Additional Reporting Requirements**

NPRI PT5                      Tetrahydrofuran

**Canada. National Pollutant Release Inventory (NPRI) (Schedule 1, Parts 1-4)**

NPRI                              Not Regulated

**Greenhouse Gases**

Not Regulated

**Controlled Drugs and Substances Act**

CA CDSI                        Not Regulated

CA CDSII                      Not Regulated

CA CDSIII                     Not Regulated

CA CDSIV                     Not Regulated

CA CDSV                      Not Regulated

CA CDSVII                    Not Regulated

CA CDSVIII                 Not Regulated

**Precursor Control Regulations**

Not Regulated

**International regulations**

**Montreal protocol**

Not applicable

**Stockholm convention**

Not applicable

**Rotterdam convention**

Not applicable

**Kyoto protocol**

Not applicable

**Inventory Status:**

Australia AICS:	On or in compliance with the inventory
Canada DSL Inventory List:	On or in compliance with the inventory
China Inv. Existing Chemical Substances:	On or in compliance with the inventory
Japan (ENCS) List:	On or in compliance with the inventory
Japan ISHL Listing:	On or in compliance with the inventory
Korea Existing Chemicals Inv. (KECI):	On or in compliance with the inventory
Mexico INSQ:	On or in compliance with the inventory
New Zealand Inventory of Chemicals:	On or in compliance with the inventory
Philippines PICCS:	On or in compliance with the inventory
Taiwan Chemical Substance Inventory:	On or in compliance with the inventory
US TSCA Inventory:	On or in compliance with the inventory
EINECS, ELINCS or NLP:	On or in compliance with the inventory

<b>16. Other information</b>
------------------------------

**Revision Date:** 20.02.2020

**Version #:** 1.2

**Source of information:** Sources of information used in preparing this SDS included one or more of the following: results from in house or supplier toxicology studies, information from the Toxicology Data Network (TOXNET), European Chemical Agency (ECHA) substance dossiers, IARC Monographs, US National Toxicology Program data, the Agency for Toxic Substances and Disease Registry, other manufacturer's SDSs and other sources, as appropriate.

**Further Information:** No data available.

**Disclaimer:** The information provided in this Safety Data Sheet (SDS) was prepared based on data believed to be accurate as of the date of this SDS. TO THE GREATEST EXTENT PERMITTED BY LAW, AVANTOR PERFORMANCE MATERIALS ("AVANTOR") EXPRESSLY DISCLAIMS ANY AND ALL REPRESENTATIONS AND WARRANTIES REGARDING THE INFORMATION CONTAINED HEREIN INCLUDING, WITHOUT LIMITATION, AS TO ACCURACY, COMPLETENESS, FITNESS FOR PURPOSE OR USE, MERCHANTABILITY, NON-INFRINGEMENT, PERFORMANCE, SAFETY, SUITABILITY AND STABILITY. This SDS is intended as a guide to the appropriate use, handling, storage and disposal of the product to which it relates by properly trained personnel, and is not intended to be comprehensive. Users of Avantor's products are advised to perform their own tests and to exercise their own judgment to determine the safety, suitability and appropriate use, handling, storage and disposal of each product and product combination for their own purposes and uses. TO THE GREATEST EXTENT PERMITTED BY LAW, AVANTOR DISCLAIMS LIABILITY FOR, AND BY USING AVANTOR'S PRODUCTS PURCHASER AGREES THAT UNDER NO CIRCUMSTANCES SHALL AVANTOR BE LIABLE FOR, SPECIAL, INDIRECT, INCIDENTAL, PUNITIVE OR CONSEQUENTIAL DAMAGES OF ANY TYPE OR KIND, INCLUDING WITHOUT LIMITATION, FOR LOSS OF PROFITS, REPUTATIONAL DAMAGE, PRODUCT RECALL OR BUSINESS INTERRUPTION.