

Product Name:	PR041WWW/0
Description:	<p>Corepoint Scientific™ General Purpose Undercounter Refrigerator Freestanding 5 Cu. Ft. Capacity  One solid swing door, right hinged, not reversible  Mechanical thermostat  Temperature range : 0°C to 10°C  Three shelves (1 fixed/2 adjustable)  Keyed door lock  Magnetic door gasket for positive seal  Recessed door handle  Leveling legs  Interior light  Probe access port (3/8")  White cabinet and door  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging One year parts and labor warranty  Exterior dimensions: 20 1/4"W x 22" D x 31 1/2"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Cycle defrost  115V, 60 Hz, .90 Amps, 1/7 HP  UL/C-UL listed  Shipping weight: 86 lbs.</p>