

Cubis® II Laboratory Balance

The ultimate in precision weighing.

Now backed by a 10-year limited warranty.

Why the Cubis® II lab balance?

The Cubis® II platform is a laboratory balance with fully customizable hardware, software, and connectivity. It offers modern user interfaces, pharmaceutical and GxP compliance, including data handling, data integrity and connectivity, ergonomic sample handling, easy process integration, and unlimited communication at the highest level of accuracy and precision.

- 3rd Generation Monolithic Weighing System
- Designed for demanding workflows with a strong focus on, connectivity, data integrity, data handling and 21 CFR Part 11 compliant
- Motorized Auto-Leveling
- Gesture control opening/closing of the draft shield
- State of the Art Connectivity: easily integrated into existing IT architecture over standard IT Protocols: HTTPS, FTPS, File Share (SMB), SQL
- For all critical environments e.g. isolator, glove box, etc.



10-Year Limited Warranty
at No Additional Cost!

Get 10-year limited warranty at no additional cost

- Industry-leading 10-year limited warranty.
- Easy Warranty Registration www.sartorius.com/cubis-ii-warranty-registration
- 10-year parts and labor www.sartorius.com/en/general-terms-and-conditions

Order Information Cubis® II preconfigured models

Model Number	Weighing capacity (g)	Readability (mg)	VWR Part no.
Cubis® II Premium Balance			
MCA10.6SM-S00	10.1	0.001	CA76377-394
MCA1203SE-S00	1,200	1	CA76377-162
MCA1203SR-S00	1,200	1	CA76377-164
MCA1220IS0-S00	12,200	100	CA76376-998
MCA124SI-S00	120	0.1	CA76377-264
MCA124SU-S00	120	0.1	CA76377-082
MCA125SI-S00	120	0.01	CA76377-344
MCA2.7SF-S00	2.1	0.0001	CA76377-398
MCA2.7SM-S00	2.1	0.0001	CA76377-396
MCA2020IS0-S00	20,200	100	CA76333-800
MCA2202S0-S00	2,200	10	CA76333-802
MCA2203SE-S00	3,200	1	CA76333-804
MCA2203SR-S00	3,200	1	CA76333-806
MCA224SI-S00	220	0.1	CA76377-314
MCA224SU-S00	220	0.1	CA76333-808
MCA225SI-S00	221	0.01	CA76377-352
MCA225SU-S00	220	0.01	CA76377-336
MCA3.6PM-S00	1.1 2.1 3.1	0.001 0.002 0.005	CA76377-364
MCA3203SE-S00	3,200	1	CA76333-912
MCA3203SR-S00	3,200	1	CA76333-810
MCA32202P0-S00	4,200 32,200	10 100	CA76333-812
MCA324SI-S00	320	0.1	CA76333-922
MCA324SU-S00	320	0.001	CA76333-924
MCA3620IS0-S00	36,200	100	CA76377-198
MCA4202S0-S00	4,200	10	CA76333-926
MCA5203SE-S00	5,200	1	CA76377-288
MCA5203SR-S00	5,200	1	CA76333-928
MCA524SI-S00	520	0.1	CA76333-930
MCA524SU-S00	520	0.1	CA76333-932

Model Number	Weighing capacity (g)	Readability (mg)	VWR Part no.
MCA6.6SF-S00	6.1	0.001	CA76377-376
MCA6.6SM-S00	6.1	0.001	CA76377-374
MCA6202S0-S00	6,200	10	CA76377-034
MCA623SE-S00	620	1	CA76377-062
MCA623SR-S00	620	1	CA76334-142
MCA7020IS0-S00	70,200	100	CA76334-144
MCA820IS0-S00	8,200	100	CA76334-146
MCE1203SR-S00	1,200	1	CA76334-148
MCE1220IS0-S00	12,200	100	CA76334-150
MCE124SU-S00	120	0.1	CA76334-152
MCE125SU-S00	120	0.01	CA76333-940
MCE2020IS0-S00	20,200	100	CA76333-942
MCE2202S0-S00	2,200	10	CA76333-944
MCE2203SR-S00	3,200	1	CA76333-946
MCE224SU-S00	220	0.1	CA76333-948
MCE225SU-S00	221	0.01	CA76333-950
MCE3.6PM-S00	1.1 2.1 3.1	0.001 0.002 0.005	CA76334-154
MCE3203SR-S00	3,200	1	CA76334-156
MCE32202P0-S00	4,200 32,200	10 100	CA76334-158
MCE324SU-S00	320	0.1	CA76333-814
MCE3620S0-S00	36,200	1,000	CA76333-816
MCE3620IS0-S00	36,200	100	CA76333-818
MCE4202S0-S00	4,200	10	CA76333-820
MCE6202S0-S00	6,200	10	CA76333-822
MCE623SR-S00	620	1	CA76333-824
MCE820IS0-S00	8,200	100	CA76333-826

Customizable configurations of Cubis® II Balances are available for special order. Please contact your sales representative for details.

*Free 10-years limited warranty offer is valid for Cubis® II lab balances purchased and located in Canada. The 10-year limited warranty is subject to the terms and limitations set out in the Sartorius Terms and Conditions of Sale (Canada), relating to the sale of Goods and Services. The warranty consists of an extended Goods Warranty Period of 10-years, and is applicable only to Cubis®II balances sold and located in Canada.



VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in Canadian dollars unless otherwise noted. Offers valid in Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.